Aqua Screen®

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
Instruction manual

Primer/ Nucleotide Mix

Lyophilized primers and deoxynucleotide triphosphates
dATP, dCTP, dGTP and dUTP, aliquoted for 25 reactions

red cap

PCR Reaction Buffer

10x buffer, 500 µl

blue cap

Positive Control DNA

DNA-fragments of *Legionella pneumophila* genome, prepared by PCR,
non-infectious, lyophilized

green cap

Internal Control DNA

plasmid DNA, lyophilized, non-infectious

yellow cap

1.2 Stability and Storage

Kit components are stable during shipping. Upon receipt, store at +2 to +8°C. After rehydration of the controls, *Primer/Nucleotide Mix*, store below -18°C and avoid repeated freezing and thawing. For repeated testing of low sample numbers, controls and the *Primer/Nucleotide Mix* should be aliquoted after rehydration. By following these recommendations, the kit is stable until the expiration date stated on the Guarantee Certificate (online: <http://www.minerva-biolabs.com>).

1.3 Supplemental Requirements

- PCR thermal cycler
- mineral oil, if required for the particular thermal cycler used
- PCR reaction tubes
- agarose gel electrophoresis apparatus
- microcentrifuge, micro pipettes and filtered tips
- deionized, DNA-free water
- polymerase

The test provides excellent results with MB Taq DNA-Polymerase (Cat # 53-0050; 53-0100; 53-200; 53-0250).

Note: We can neither guarantee a high level of sensitivity nor compatibility with other polymerases. If you intend to test our MB Taq DNA-Polymerase in parallel with your in-house polymerase, please feel free to contact us and get a gratis MB Taq sample (10 units). However, if you want to use your own polymerase, it may be necessary to use the specific buffer provided with this polymerase.

2. Application and Test Principle

The Aqua Screen® PCR system is a test for the qualitative detection of *Legionella pneumophila* in water samples. The test is based on the polymerase chain reaction. The supplied primer set is specific for a segment of the *mip* region of the *Legionella pneumophila* genome. The amplified PCR product is 275 bp long and can be rendered visible directly in the agarose gel. By using the supplied internal control, false-negative results (e.g. due to inhibition of the reaction) can be excluded individually for each sample. The internal control amplicon is 191 bp in size. Detailed studies confirm the high specificity and sensitivity of Aqua Screen®. The selected template is highly preserved within the *Legionella pneumophila* species. Cross activity with other *Legionella* species or species typical for water contamination is not known. The detection limit of the PCR unit in combination with the Aqua Screen® filtration unit is 300 genome equivalents per filtered volume.

The heat-labile Uracil-DNA Glycosylase (UNG) is suitable for preventing carry-over contamination between PCRs. This technique relies on the incorporation of deoxyuridine triphosphate (dUTP) during all amplification reactions and the pretreatment of all successive PCR mixtures with the heat-labile UNG. The UNG cleaves DNA at any site where a deoxyuridylate residue has been incorporated. Subsequently, the resulting abasic sites are hydrolyzed due to the high temperature during the initial denaturation step, and cannot serve as PCR templates any longer. The heat-labile UNG is inactivated at the same time. Native DNA (e.g., the template DNA) does not contain uracil and is therefore not degraded by this procedure. Therefore dTTP is replaced by dUTP in this kit. UNG is not provided with this kit.

3. Test Protocol

3.1 Preparation of Sample Material

The preparation of the sample material is explained in the Aqua Screen® extraction kit manual. The extracts can be stored at a temperature of at least -18°C for a period of one year. Repeated freezing and thawing, or storage in the refrigerator for more than 12 hours should be avoided. The sample should not contain more than 100 µg/ml DNA.

3.2 Rehydration of the Reagents

1. centrifuge tubes with lyophilized components (5 sec at maximum speed)
2. add appropriate amount of deionized, DNA-free water:

<i>primer/nucleotide mix</i> (per portion of 25 reactions)	65 µl
<i>positive control</i>	300 µl
<i>internal control</i>	300 µl
3. incubate for 5 minutes at room temperature
4. vortex and centrifuge again

Keep reagents on ice and store below -18°C after rehydration.

3.3 Thermal Profile

The programming process of your cycler is explained in the manual of the instrument.

Program:

1 cycle 94°C for 2 min
35 cycles 94°C for 30 sec
60°C for 30 sec
72°C for 30 sec
cool down to 4 to 8 °C

The incubation time depends on the polymerase used. Hot start enzymes need to be activated at 94°C. Please see polymerase data sheet for duration.

3.4 The PCR Master Mix

Total volume per reaction is 25 µl. When setting up reactions, calculations should also include positive and negative controls. Pipet mastermix on ice into a 1.5 ml reaction tube and mix gently.

Pipetting Schemes:

	for 1 reaction	for 25 reactions*
water	12.3 µl	307.5 µl
10x reaction buffer (blue cap)	2.5 µl	62.5 µl
primer/nucleotide mix (red cap)	2.5 µl	62.5 µl
internal control (yellow cap)	2.5 µl	62.5 µl
polymerase (5 U/µl)	0.2 µl	5.0 µl

*the equivalent of the content of one red-capped vial.

For other polymerase concentrations the amount of water needs to be adjusted. The mastermix can also be prepared for 25 reactions, aliquoted as needed and stored below -18°C for up to 3 months.

Aliquot 20 µl of master mix into each PCR reaction tube.

Add 5 µl of sample (as described above) to PCR reaction tube per sample being tested. After pipetting the negative control, the tube must be sealed before proceeding with the samples. Also pipetting of the samples and sealing the tubes must be completed before proceeding with the positive control (green cap) in order to avoid cross contamination.

3.5 Agarose Gel Run

- 1.5% standard agarose gel with 5 mm-comb
- load 5 μ l of each PCR reaction, mixed with bromophenol blue loading buffer per lane

Only bromophenol blue in a low concentration should be used as a run marker.

- stop electrophoresis after 2 cm run distance (depending on the electrophoresis chamber used e.g. run for 20 minutes at 100 V)

3.6 Gel Evaluation

If internal control DNA was used, a distinct 191 bp band should appear in every lane indicating a successfully performed PCR. This band may fade out with increased amount of amplicons formed, caused by *Legionella pneumophila* DNA loads of $> 5 \times 10^6$ copies/ml. The initial concentration of positive control DNA exceeds 5×10^6 copies/ml in order to account for DNA loss resulting from repeated freezing and thawing.

Relevant amplicon sizes:

Internal control 191 bp

Legionella pneumophila 275 bp

Results of a successfully performed analysis:

negative control band at 191 bp

positive control band at 275 bp, possibly an additional band at 191 bp

No amplification of control DNA may be due to the following reasons:

- activity of polymerase is insufficient
- control DNA tubes have not been spun down before rehydration
- programming mistake
- pipetting mistake


Before rerun of a negative and a positive control please check thermocycler protocol and pipetting scheme. When using polymerases other than the MB TAQ DNA Polymerase, please note the comments made under chapter 1.3. The concentration can then be raised up to 2.5 U/reaction. Please note the complete change of the pipetting scheme.

Interpretation of possible band patterns:

band at 191 bp	negative sample
band at 275 bp and at 191 bp	<i>Legionella pneumophila</i> -positive sample with weak contamination
strong band at 275 bp	<i>Legionella pneumophila</i> -positive sample, highly contaminated
no band	PCR inhibition

With conventional PCR a quantification of the DNA is only conditionally possible. This test provides information regarding the absence or presence of a contamination. The detection limit of conventional PCR averages 6 copies per PCR. Hence it follows that with a positive sample (positive band is visible in the gel), at least 120 *Legionella* particles were in the sample water.

With Aqua Screen® designed for high sensitivity and therefore prone to nonspecific annealing, bands of various length that are less intensive can be produced, but not indicate positive results. Possible primer self-annealing produces another band of 80-90 bp in length, but also does not affect the precision or results of the test.

This page is intentionally blank.

Appendix

Limited Product Warranty

This warranty limits our liability for replacement of this product. No warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided. Minerva Biolabs shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

Notice to Purchaser

This product is optimized for use in the Polymerase Chain Reaction („PCR“) covered by patents owned by Hoffmann-La Roche, Inc., and F. Hoffmann-La Roche Ltd. („Roche“). No license under these patents to use the PCR process is conveyed expressly or by implication to the purchaser of this product. Minerva Biolabs does not encourage or support the unauthorized or unlicensed use of the PCR process. Use of this product is restricted to persons that either have a license to perform PCR or are not required to obtain a license.

Trademarks

Aqua Screen is a registered trademark of Minerva Biolabs.

Minerva Biolabs' International Aqua Screen® Distributors

Belgium, Luxemburg

Lucron Bioproducts bvba
Tel: +32 92 82 05 31
Fax: +32 92 82 05 32
Email: info@lucronbioproducts.com
Web: www.lucronbioproducts.com

Canada

Medicorp. Inc
Tel: +1 514-7331 900
Fax: +1 514-7331 212
Email: mktg@medicorp.com
Web: www.medicorp.com

China

Beijing Vian-Saga Ltd.
Tel: +86-10 8411 8493
Fax: +86-10 8411 8494
Email: Michael@vian_saga.com
Web: www.vian-saga.com

Croatia

Gorea plus d.o.o.
Tel: +385-1 3369 610
Fax: +385-1 3369 611
Email: management@gorea-plus.hr
Web: www.gorea-plus.hr

Czech Republic

Bio-Consult spol.sro.
Tel: +420-02 4447 1239
Fax: +420-02 4447 1239
Email: info@bioconsult.cz
Web: www.bioconsult.cz

Finland

Immuno Diagnostic Oy
Tel: +358 3 615 370
Fax: +358 3 682 2039
Email: info@immunodiagnostic.fi
Web: www.immunodiagnostic.fi

Great Britain

Cambio Ltd.
Tel: +44 1954 210 200
Fax: +44 1954 210 300
Email: support@cambio.co.uk
Web: www.cambio.co.uk

Greece

Varelas S.A.
Tel: +30 210-5281-901
Fax: +30 210-5220-926
Email: varelas@otenet.gr
Web: www.otenet.gr

Bioanalytica S.A.
Tel: +30 210 6400 318
Fax: +30 210 6462 748
Email: bioanalyt@hol.gr
Web: www.bioanalytica.gr

Ireland

Medical Supply Company
Tel: +353 1 8224 222
Fax: +323 1 8224 100
Email: support@cambio.co.uk
Web: www.cambio.co.uk

Israel

Origolab Ltd.
Tel: +972-2-566 9285
Fax: +972-2-561 2120
Email: Origolab@netmedia.net.il

Italy

Biospa
Tel: +39 2 891 391
Fax: +39 2 891 20996
Email: biospa@spaspa.it
Web: www.spaspa.it/biospa.htm

Japan

Funakoshi Co. Ltd.
Tel: +81 3 5684 1615
Fax: +81 3 5684 1775
Email: info@funakoshi.co.jp
Web: www.funakoshi.co.jp

Lithuania

Interlux
Tel: +370 5 27 86 850
Fax: +370 5 27 96 728
Email: spirit@interlux.lt
Web: www.interlux.lt

Malaysia

Helix Biotech Sdn Bhd
Tel: +60-3-9076 8010
Fax: +60-3-9076 8007
Email: helixbio@tm.net.my
Web: www.helixbiot.com

Netherlands

Lucron Bioproducts B.V.
Tel: +31 485 51 16 75
Fax: +31 485 51 20 52
Email: lucron@lucron.nl
Web: www.lucron.nl

Norway

E. Pedersen & Son
Tel: +47 22 95 59 59
Fax: +47 22 95 59 40
Email: eped@eped.com
Web: www.eped.com

Korea

Bio and Information Corp
Tel: +82 31 7139 439
Fax: +82 31 7139 438
Email: sales@bioninfo.com
Web: www.bioninfo.com

Poland

STI
Tel: +48-61-6417 759
Fax: +48-61-6417 758
Email: office@sti.biz.pl
Web: www.sti.biz.pl

Portugal

Quilaban Lda.
Tel: +21 923 63 50
Fax: +21 923 63 89
Email: quilaban@quilaban.pt
Web: www.quilaban.pt

Spain

LacClinics S.A.
Tel: +34 3 446 47 00
Fax: +34 3 348 10 39
Email: info@labclinics.com
Web: www.labclinics.com

Sweden

ANL-Produkter AB
Tel: +46 8 99 00 90
Fax: +46 8 99 20 40
Email: info@anl.se
Web: www.anl.se

Switzerland

Socochim SA
Tel: +41 21 721 04 50
Fax: +41 21 721 04 51
Email: info@socochim.ch
Web: www.socochim.ch

Thailand

Biomed Diagnostics Co. Ltd
Tel: +66-2-8796 026
Fax: +66-2-8796 065
Email: somboon@biomedthai.com

Turkey

Genomed Saglik Hizmetleri A.S.
Tel: +90 212 248 20 00
Fax: +90 212 220 15 64
Email: muratyazici@genomedtr.com
Web: www.genomed-biotech.com

United Arab Emirates

Mix Max Trading CCC
Tel: +971-4-3535 756
Fax: +971-4-3538 891
Email: manoj@mixmaxe.com
Web: www.mixmaxe.com

Related Products

Supplements

53-0050	MB TAQ DNA Polymerase	50 units
53-0100	MB TAQ DNA Polymerase	100 units
53-0200	MB TAQ DNA Polymerase	200 units
53-0250	MB TAQ DNA Polymerase	250 units

Filtration Unit

31-0300	Filtration Rack	1 unit
---------	-----------------	--------

Aqua Screen® Fast Extract Kit

32-1010	Extraction Kit	10 extractions
32-1050	Extraction Kit	50 extractions
32-1200	Extraction Kit	200 extractions

Diagnostic Kits for Conventional PCR

33-1025	<i>Legionella</i> spp. Detection Kit for conv. PCR	25 tests
33-1100	<i>Legionella</i> spp. Detection Kit for conv. PCR	100 tests
33-1250	<i>Legionella</i> spp. Detection Kit for conv. PCR	250 tests
34-1025	<i>Legionella pneumophila</i> Detection Kit for conv. PCR	25 tests
34-1100	<i>Legionella pneumophila</i> Detection Kit for conv. PCR	100 tests
34-1250	<i>Legionella pneumophila</i> Detection Kit for conv. PCR	250 tests

Diagnostic Kits for Real-Time PCR

33-2025	<i>Legionella</i> spp. Detection Kit for qPCR type 1	25 tests
33-2100	<i>Legionella</i> spp. Detection Kit for qPCR type 1	100 tests
33-2250	<i>Legionella</i> spp. Detection Kit for qPCR type 1	250 tests
33-4025	<i>Legionella</i> spp. Detection Kit for qPCR type 2	25 tests
33-4100	<i>Legionella</i> spp. Detection Kit for qPCR type 2	100 tests
33-4250	<i>Legionella</i> spp. Detection Kit for qPCR type 2	250 tests
34-2025	<i>L. pneumophila</i> Detection Kit for qPCR type 1	25 tests
34-2100	<i>L. pneumophila</i> Detection Kit for qPCR type 1	100 tests
34-2250	<i>L. pneumophila</i> Detection Kit for qPCR type 1	250 tests
34-4025	<i>L. pneumophila</i> Detection Kit for qPCR type 2	25 tests
34-4100	<i>L. pneumophila</i> Detection Kit for qPCR type 2	100 tests
34-4250	<i>L. pneumophila</i> Detection Kit for qPCR type 2	250 tests

Genomic DNA Extracts 100 µl each

51-0101	<i>Legionella pneumophila</i> , ATCC 33152	+/- 10 ng / 100 µl
51-1514	<i>Legionella pneumophila</i> subsp. <i>fraseri</i> , DSMZ 7514	+/- 10 ng / 100 µl
51-1515	<i>Legionella pneumophila</i> subsp. <i>pascullei</i> , DSMZ 7515	+/- 10 ng / 100 µl
51-1368	<i>Legionella bozemanii</i> , NC 011368	+/- 10 ng / 100 µl
51-1370	<i>Legionella dumoffii</i> , NC 011370	+/- 10 ng / 100 µl

Quantification Standards

52-0101	<i>Legionella pneumophila</i> DNA Standard	1x10 ⁶ genomes/µl
---------	--	------------------------------

matching funnels

31-3050	250ml, polycarbonate, for single use, steril	50 funnels
---------	--	------------