**Onar®Lp-QP**  
**Legionella pneumophila Diagnostic Kit for qPCR**  
**- Type 2 -**

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## 1. Reagents and Materials

### 1.1 Test Kit Components

*Primer/Probe/Nucleotide Mix (PPN Mix)*  
Primer set, probes and deoxynucleotide triphosphates dATP, dCTP, dGTP and dUTP, aliquoted for 25 reactions, lyophilized

red cap  

REAC	B
------	---

*Rehydration Buffer*  
1.8 ml

blue cap  

REAC	A
------	---

*Positive Control DNA*  
Fragment of *Legionella pneumophila* DNA prepared by PCR, non-infectious, lyophilized

green cap  

CONTROL	+
---------	---

*Internal Control DNA*  
Plasmid DNA, non-infectious, lyophilized

yellow cap  

CONTROL	i
---------	---

*PCR Grade Water*  
2 ml

white cap  

REAC	H
------	---

Please find the guarantee certificate and the conformity declaration on our website.

### 1.2 Stability and Storage

Kit components are stable during shipping. Upon receipt, store at +2°C to +8°C. After rehydration of the *Primer/Probe/Nucleotide Mix* store below -18°C and avoid repeated freezing and thawing. Protect the *Primer/Probe/Nucleotide Mix* from light. For repeated testing of low sample numbers, the *Primer/Probe/Nucleotide Mix* and the controls should be aliquoted. By following these recommendations, the kit is stable until the expiration date stated on the box label.

### 1.3 Supplemental Requirements

- qPCR machine
- corresponding PCR reaction tubes
- microcentrifuge, micropipettes and filtered tips (1 - 1000 µl)
- polymerase

**Please do not use the product after the date of expiry. Discard remaining reagents.**

**The test provides excellent results with MB Taq DNA Polymerase (Cat # 53-0050/0100/0200/0250). We can neither guarantee a high level of sensitivity nor compatibility with other polymerases. If you intend to test our MB Taq DNA Polymerase in parallel with your in-house polymerase, please feel free to contact us and get a gratis MB Taq sample (10 units).**

## 2. Application and Test Principle

*Legionella pneumophila* ubiquitously colonizes surface water and wet soil and is the cause of the legionnaires' disease. Transmission is most likely aerogenic. Typical symptoms are high fever, cough, thorax ache, diarrhea, discomposure, sometimes with heavy progression.

The Onar®Lp-QP test system is an in vitro test for the qualitative and the quantitative diagnosis of *Legionella pneumophila* in clinical samples. The test is based on the polymerase chain reaction. The supplied primer set is specific for a segment of the *mip* region of the Legionella genome. The target probe emits fluorescent light at 520 nm.

By using the supplied internal control, false-negative results (e.g. due to inhibition of the reaction by the sample matrix) can be excluded individually for each sample. The internal control is detected by another probe emitting fluorescent light at 610 nm.

Detailed studies confirm the high specificity and sensitivity of the Onar®Lp-QP test system. The selected template is highly preserved within *Legionella pneumophila*. Cross activity with other Legionella species or species typical for water contamination is not known.

The heat-labile Uracil-DNA Glycosylase (UNG) is suitable for preventing carry-over contamination between PCRs. This technique relies on the incorporation of deoxyuridine triphosphate (dUTP) during all amplification reactions and the pretreatment of all successive PCR mixtures with the heat-labile UNG. The UNG cleaves DNA at any site where a deoxyuridylate residue has been incorporated. Subsequently, the resulting abasic sites are hydrolyzed due to the high temperature during the initial denaturation step, and cannot serve as PCR templates any longer. The heat-labile UNG is inactivated at the same time. Native DNA (e.g., the template DNA) does not contain Uracil and is therefore not degraded by this procedure. UNG is not provided with this kit.

The detection limit of Onar®Lp-QP is 10 GC/PCR for *L. pneumophila*. Due to the broad linear detection range, reliable results can be obtained with Onar®Lp-QP without timeconsuming conditioning of the sample material. The samples are markedly more stable compared to the antigen test.

### 3. Test Protocol

#### 3.1 Preparation of Sample Material

Nasopharyngeal swabs, nasal and pharyngeal secretions, sputum, provoked sputum, bronchial lavage, tissue and infected cells or cultures can be used as test materials. The quality of sample-taking has a major influence on the reliability of the test results. Material from the lower respiratory tract is particularly suitable for detecting *Legionella pneumophila*. A fine catheter and a vacuum pump with secretion trap can be used for nasopharyngeal secretions. The oral cavity can be rinsed beforehand, and saliva and rinsing water discarded. Taking a sample from the oral cavity should not be preceded by eating or gargling. The swab is turned upwards behind the arch of the palate and wiped. When taking a sample through the nasal cavity, the thin swab is advanced into the nasopharynx and wiped repeatedly. Ensure that adequate material is obtained. The swab should be washed out in 1 ml transfer buffer (e.g. 50 mM Tris-HCL, pH 7.5, 0,1% sodium azide, 5% sodium dodecylsulfate). The sample material is stable for a few days and should be transported cooled.

DNA extraction with a commercially available DNA extraction kit (e.g. MB DNA Extraction Kit Cat # 56 -1100) is always advisable when preparing the samples in order to remove inhibitors of the PCR safely and to concentrate the *Legionella* DNA at greater sample volumes. The obtained DNA extract can be used directly for the Onar<sup>®</sup>Lp-QP test. The extracts can be stored at a temperature of at least -18°C for a period of one year. Repeated freezing and defrosting, or storage in the refrigerator for more than 12 hours should be avoided. The sample should not contain more than 100 µg/ml DNA.

#### 3.2 Rehydration of the Reagents

1. centrifuge tubes with lyophilized components (5 sec at maximum speed)
2. add 365 µl of *Rehydration Buffer* to *Primer/Probe/Nucleotide Mix*
3. add *appropriate amount of deionized, DNA-free water*  
*Positive Control DNA* 300 µl  
*Internal Control DNA* 300 µl
4. incubate for 10 minutes at room temperature
5. vortex and centrifuge again

**Keep reagents on ice and store below -18°C after rehydration.**

### 3.3 Experimental Protocols and Result Reading

for **ABI 7500**

#### Detector Settings:

*L. pneumophila* Target Probe: Reporter - FAM Quencher - none  
Internal Control Probe: Reporter - VIC Quencher - none

**The ROX Reference needs to be disabled. Activate both detectors for each well.**

**Measurement of fluorescence during extension.**

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#### Result Reading:

- Enter the following basic setting at the right task bar:  
Data: Delta RN vs. Cycle  
Detector: FAM and VIC  
Line Color: Well Color
- Open a new window with for the *Graph settings* by clicking the right mouse button  
Select the following setting and confirm with *ok*:  
Real Time Settings: Linear  
Y-Axis Post Run Settings: Linear and Auto Scale  
X-Axis Post Run Settings: Auto Scale  
Display Options: 2
- Initiate the calculation of the ct-values and the graph generation by clicking on *Analyze* within the report window.
- Pull the threshold line into the graph. Adapt the threshold line to the initial linear section of the positive control reaction.
- Samples showing no ct-value can be considered as negative.

### 3.4 PCR Master Mix Setup

Total volume per reaction is 25  $\mu$ l. When setting up reactions, calculations should also include positive (PC) and negative controls (NC). Pipet master mix on ice into a 1.5 ml reaction tube and mix gently.

Pipetting scheme:

	for 1 reaction	for 25 reactions
primer/probe/nucleotide mix	14.0 $\mu$ l	350.0 $\mu$ l
Internal Control DNA (yellow cap)	1.0 $\mu$ l	25.0 $\mu$ l
polymerase (5 U/ $\mu$ l)	0.2 $\mu$ l	5.0 $\mu$ l
+ template DNA, NC or PC	10.0 $\mu$ l	

For other polymerase concentrations the amount of enzyme and the amount of water added to the mix need to be adjusted.

Aliquot 15  $\mu$ l of master mix into each PCR reaction tube.

**The total duration from master mix preparation to PCR cycling must not exceed 60 minutes to avoid a decrease in the fluorescent signal.**

Add 10  $\mu$ l of sample (as described above) to PCR reaction tube per sample being tested. After pipetting the negative control (10  $\mu$ l of water or negative control of DNA extraction/reaction), the tube must be sealed before proceeding with the samples. Also pipetting of the samples and sealing the tubes must be completed before proceeding with the positive control (10  $\mu$ l /reaction) in order to avoid cross contamination.

### 3.5 Result Interpretation

The presence of *L. pneumophila* in the sample is indicated by an increasing fluorescence signal in the *L. pneumophila* channel during PCR. The concentration of the contaminant can be calculated by a software comparing the crossing cycle number of the sample with a standard curve created in the same run.

A successfully performed PCR without inhibition is indicated by an increasing fluorescence signal in the internal control channel, provided the Internal Control was added to the master mix. *L. pneumophila* DNA and Internal Control DNA are competitors in PCR. Because of the very low concentration of Internal Control in the PCR mix, the signal strength in this channel is reduced with increasing *L. pneumophila* DNA loads in the sample.

<i>L. pneumophila</i> PCR	Internal Control	Interpretation
positive	irrelevant	<i>L. pneumophila</i> positive
negative	negative	PCR inhibition
negative	positive	<i>L. pneumophila</i> negative

#### 4. Trouble shooting

Before repeating a negative and a positive control run please check the cycler protocol and pipetting scheme. When using polymerases other than the MB TAQ DNA Polymerase, please note the comments made under chapter 1. The concentration can then be raised up to 2.5 U/reaction. Please note the complete change of the pipetting scheme.

Insufficient (inhibition) or no amplification of the internal control may be due to the following reasons:

- activity of *Taq* polymerase is insufficient or enzyme not compatible with the kit buffer
- programming mistake
- pipetting mistake
- sample elution buffer and sample used to complete the negative control are different

#### 5. Instrument Compatibility

Instrument	type 1	type 2
LightCycler® 1.2	+ <sup>v</sup>	–
LightCycler® 1.5	+	–
LightCycler® 2.0	+	o
LightCycler® 480	+	o
Rotorgene 3000	+	o
Rotorgene 6000 (5-plex)	+ <sup>v</sup>	o
ABI Prism® 7000	–	+
ABI Prism® 7300	–	+
ABI Prism® 7500	–	+ <sup>v</sup>
ABI Prism® 7700	–	+
ABI Prism® 7900	–	+
iCycler iQ®	o	o
iQ™5	+	o
Opticon 2	–	o
Chromo 4	o	o
MX3000P	o	o
MX4000	o	o

+ = recommended Kit version

– = the internal control is not detectable

o = untested but presumed to be compatible

v = validated

## **Appendix**

### Limited Product Warranty

This warranty limits our liability for replacement of this product. No warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided. Minerva Biolabs shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

### Notice to Purchaser

This product is optimized for use in the Polymerase Chain Reaction („PCR”) covered by patents owned by Hoffmann-La Roche, Inc., and F. Hoffmann-La Roche Ltd. („Roche”). No license under these patents to use the PCR process is conveyed expressly or by implication to the purchaser of this product. Minerva Biolabs does not encourage or support the unauthorized or unlicensed use of the PCR process. Use of this product is restricted to persons that either have a license to perform PCR or are not required to obtain a license.

### Limited License

The use of this product for the detection of Legionella infections is covered by a Scorpions license from DxS Ltd. Further information on Scorpions licenses can be obtained from DxS Ltd, 48 Grafton St, Manchester, UK M13 9XX.

## Related Products

### MB Taq DNA Polymerase

53-0050/-0100/-0200/-0250 MB Taq DNA Polymerase 50/100/200/250 units

### Clinical Diagnostic Kits for qPCR

21-4025/-4100/-4250	Onar <sup>®</sup> Ls-QP <i>Legionella</i> species	25/100/250 tests
21-5025/-5100/-5250	Onar <sup>®</sup> Lp-QP <i>Legionella pneumophila</i>	25/100/250 tests
20-3025/-3100/-3250	Venor <sup>®</sup> Mp-QP <i>Mycoplasma pneumoniae</i>	25/100/250 tests
24-3025/-3100/-3250	Onar <sup>®</sup> Ct <i>Chlamydia trachomatis</i>	25/100/250 tests
23-2025/-2100	Onar <sup>®</sup> Tp <i>Treponema pallidum</i>	25/100 tests
25-3025/-3100/-3250	Onar <sup>®</sup> Pertussis <i>Bordetella pertussis/parapertussis</i>	25/100/250 tests

### Quantification Standards, 100 µl each, 1x10<sup>8</sup> genomes/µl

52-0101	<i>Legionella pneumophila</i> DNA Standard
52-0119	<i>Mycoplasma pneumoniae</i> DNA Standard

### Genomic DNA Extracts, 100 µl each, +/- 10 ng / 100 µl

51-0566	<i>Streptococcus pneumoniae</i> , DSMZ 20566
51-5571	<i>Bordetella pertussis</i> , DSMZ 5571
51-3415	<i>Bordetella parapertussis</i> , DSMZ 13415
51-0177	<i>Ureaplasma urealyticum</i> , NC 10177
51-0111	<i>Mycoplasma hominis</i> , NC 010111
51-0112	<i>Mycoplasma orale</i> , NC 010112
51-0117	<i>Mycoplasma fermentans</i> PG19, NC 010117
51-0119	<i>Mycoplasma pneumoniae</i> , NC 010119
51-0129	<i>Mycoplasma arginini</i> , NC 010129
51-0162	<i>Mycoplasma arthritidis</i> , NC 010162
51-0195	<i>Mycoplasma genitalium</i> , NC 010195
51-1746	<i>Mycoplasma penetrans</i> , NC 11746

### DNA Remover™

15-2025	DNA Decontamination Reagent, spray bottle	250 ml
15-2200	DNA Decontamination Reagent, refill bottles	4 x 500 ml