Art. Nr. / Order No. 24-2025/-2100/-2250

Onar[®]Ct

***Chlamydia trachomatis* Diagnostic Kit for qPCR**
version 2.0

Gebrauchsinformation / Instructions for Use

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Reagents for 25/100/250 reactions

Hersteller / Manufacturer:

Minerva Biolabs GmbH, Koepenicker Strasse 325, 12555 Berlin, Germany

ZUR IN VITRO DIAGNOSTIK

FOR IN VITRO DIAGNOSTIK USE

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ABI Prism® 7900

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INDICATION

Onar[®]Ct is used for quantitative and direct detection of *Chlamydia trachomatis* in urine, ejaculates and swabs from the eye as well as endocervical and urethral areas.

EXPLANATION OF THE TEST

This kit is a clinical in vitro test system based on the PCR method. The test can be performed within 2 to 3 hours and allows for the quantification of the pathogen for significant clinical valuation and infection process monitoring.

TEST PRINCIPLE

The primer set and probes are specific for a segment of the chromosomal *momp* gene (*major outer-membrane protein*). The genome of *Chlamydia trachomatis* contains only one copy of the *momp* gene locus. In contrast to other test methods amplifying the multicopy *cryptic plasmid*, this test design allows for a quantitative interpretation of the results.

False negative results due to PCR inhibitors or improper DNA extraction are detected by the internal amplification control. The Internal Control DNA can be added to the sample prior to DNA extraction and analysis for verification of the complete process (DNA extraction and PCR reaction). The Internal Control DNA can also be added directly to the PCR master mix to act as a PCR control only. The amplification of the control reaction is detected at 610 nm (ROX[™] channel) and the pathogen-specific sequence at 520 nm (FAM[™] channel).

The Trachomatis Mix contains dUTP instead of dTTP, so the option is available to degrade amplicons from previous analysis by use of uracil-DNA glycosylase (UNG). Thus the occurrence of false-positive result can be minimized. UNG is not part of Onar[®]Ct.

REAGENTS

Each kit contains reagents for 25, 100 or 250 reactions. The expiry date of the unopened package is marked on the package label. The kit components are stored until use at +2 to +8 ° C. The lot specific quality control certificate (Certificate of Analysis) can be downloaded from our website (www.minerva-biolabs.com).

Kit Component Label Information	Quantity			Cap Colour
	25 reactions Order No. 24-2025	100 reactions Order No. 24-2100	250 reactions Order No. 24-2250	
Trachomatis Mix	1 x lyophilized	4 x lyophilized	10 x lyophilized	red
Rehydration Buffer A	1 x 1.8 ml	2 x 1.8 ml	3 x 1.8 ml	blue
Rehydration Buffer B	1 x 1.8 ml	2 x 1.8 ml	3 x 1.8 ml	blue
Positive Control DNA	1 x lyophilized	1 x lyophilized	2 x lyophilized	green
Internal Control DNA	1 x lyophilized	1 x lyophilized	2 x lyophilized	yellow
PCR grade Water	1 x 2.0 ml	1 x 2.0 ml	2 x 2.0 ml	white

NEEDED, BUT NOT INCLUDED IN THE KIT

The Onar[®]Ct kit contains the reagents for the specific detection of *Chlamydia trachomatis*. General industrial supplies and reagents, usually available PCR laboratories are not included:

- qPCR device with filter sets for the detection of the fluorescence dyes FAM[™] and ROX[™]
- PCR reaction tubes for the specific qPCR device
- 1,5 ml reaction tubes, DNA- and RNA-free
- Micro centrifuge for 1.5 ml and PCR reaction tubes
- Pipettes with corresponding filter tips to prepare and dispense the reaction mix (10, 100 und 1000 µl)

SPECIMEN - Collection and Storage

Urine Specimens

1. The patient should be advised not to urinate for at least two hours prior to specimen collection
2. Collect 10 ml of first-catch urine in a clean container of polypropylene without preservatives.

- The urine specimens may be stored at room temperature (19 - 23 °C) if testing of samples can be completed within 24 hours. Samples have to be stored at 2 - 8 °C in the fridge, if the testing is done within 7 days. Store the samples at < - 20 °C for up to 30 days.

Ejaculate Specimens

1. The ejaculates must be collected the same day to be going to be shipped to the laboratory.
2. If a ejaculates cannot be processed directly after receipt in the laboratory, it must be stored at below - 20 ° C and tested within 30 days.

Swab Specimens

For the collection of eye, end cervical and urethral swabs, Darcon® or rayon swabs can be used:

1.
 - MASTASWAB, Mast Diagnostica, Cat. No. 800108 MD 508
 - Sterile Swab Applicators, mico Rheologics/Copan, Cat. No. 502CS01
 - Forensic swabs predetermined break zone, sterile wrapped, Copan, Cat. No. 3520CS01-4N6
 - Flocked swabs, microRheologics/Copan, Cat. No. 552C nylon
-

For the collection of *endocervical* specimens: The subject must refrain from any kind of sexual activity, douching, and inserting any intravaginal products for at least 48 hours prior to the collection of the specimens. If excess mucus or menses clot has accumulated near the cervical os, a large cotton-tipped cotton swab may be used to gently remove this material before inserting the swab. A swab is inserted and rotated 360 degrees in all four quadrants of the vaginal vault. The swab is then placed into a sterile tube, sealed and stored for further analysis.

2.

For the collection of *urethral* specimens: 1 hour prior to specimen collection should not be urinating. Insert the swab 2-4 cm into the urethral chanel turn for about 3 - 5 seconds under slight pressure,

Samples must be shipped refrigerated. If the swab specimens are not processed directly after receipt in the laboratory, they have to be stored at 2 - 8 °C and have to be processed within seven days. If they are not tested within seven days after collection, the samples have to be stored at < -20 °C and tested within up to 30 days.

3.

SPECIMEN - Preparation

It is necessary that templates for PCR analysis are prepared by DNA extraction using a suitable DNA extraction kit. Minerva Biolabs DNA Extraction Kit (Cat. No. 56-1100) provides protocols for urine and pool urine, ejaculates and swab samples. When using extraction kits from other manufacturers or applying own procedures compatibility of the elution buffer with the diagnostic kit must be ensured.

The extracts can be stored at a temperature of below -18 °C for a period of one year. Repeated freezing and thawing or storage in the refrigerator for more than 12 hours should be avoided. The sample should not contain more than 100 µg/ml of DNA.

The internal amplification control can be used to verify the DNA extraction. 10 µl of the Internal Control should be added directly to the sample volume used for extraction. The master mix needs to be prepared without Internal Control. An adequate volume of water equals for the loss of volume.

For urine as sample material, this alternative manual protocol has been proved of value:

1. Centrifuge 1 ml urine or 1 ml transport medium (swabs) at 18.000 x g for 20 min and discard the supernatant.
 2. Wash the sediment in 500 µl 1 x PBS pH 7.3 (phosphate-buffered saline), centrifuge (18.000 x g for 20 min) and resuspend the sediment in 50 µl 10 mM Tris-HCl, pH 7.2.
 3. Add 5 µl proteinase K (1 mg/ml !!) and incubate at 56 °C for 30 min.
 4. Incubate at 95 °C for 10 min to lyse the cells and inactivate the proteinase K .
 5. After a short centrifugation, the sample is ready for PCR.
 6. If not tested immediately, samples should be stored at -20 °C.
-

PRECAUTIONS

For in vitro diagnostic use. This kit should be used only by trained persons. All samples should be considered potentially infectious and handled at the local or national regulations. This kit does not contain hazardous substances and may be disposed of according to local regulations.

TEST PROCEDURE

The test should be carried out with negative and positive controls and patient samples in duplicate. For quantification, a dilution series of an appropriate standard, e. g. Minerva Biolabs *Chlamydia trachomatis* Quantification Standard (Order No. 52-0440) should be processed in parallel. All reagents and samples must be equilibrated to +2 to +8 ° C prior use.

1. Rehydration of the Reagents

After reconstitution, the reagents should be stored at below -18 ° C. Repeated freezing and thawing should be avoided and reconstituted controls (internal amplification control and the positive control) stored in aliquots.

1.	Trachomatis Mix Internal Control Positive Control	red cap yellow cap green cap	Spin all lyophilized components for 5 sec at maximum speed of the mini centrifuge
2.	Trachomatis Mix	red cap	Add 365 µl Rehydration Buffer A or B (blue cap). The buffer can be selected according to Annex 2 (page 26).
3.	Internal Control	yellow cap	Add 300 µl PCR grade Water (white cap)
4.	Positive Control	green cap	Add 300 µl PCR grade Water (white cap)
5.	Trachomatis Mix Internal Control Positive Control	red cap yellow cap green cap	Incubate 5 min at room temperature
6.	Trachomatis Mix Internal Control Positive Control	red cap yellow cap green cap	Vortex briefly and spin for 5 sec

2. Preparation of the Reaction Mix

Prepare the required amount of master mix at room temperature in a 1.5 ml reaction tube for all control and test reactions.

1.	for 1 reaction	for 25 reactions
Trachomatis Mix	14,0 µl	350,0 µl
Internal Control DNA	1,0 µl	25,0 µl

2. Homogenize the reaction mix by carefully snipping the tube

3. Add 15 µl to each PCR tube, discard remaining material.

3. Loading the test tubes

- Negative Control:
1. Add 10 μ l elution buffer from DNA extraktion (ref. Chapter sample material) or PCR grade Water (white cap) to each negative control reaction.
 2. Sample: Add 10 μ l of each sample.
 3. Positive Control: Add 10 μ l Positive Control (green cap) to each positive control reaction.
 4. Spin all PCR tubes briefly, load the qPCR cycler and start the program.
-

Preparation of the master mix and sample loading should not take more than 45 minutes to avoid a reduction in the fluorescent signal. The pipetting sequence should be respected and the tubes closed after each sample load.

4. Starting the reaction

1. Load the cycler, check each PCR tube and the cycler lid for tight fit.
 2. Program the qPCR cycler or check stored temperature profiles. See Appendix for temperature profiles of selected qPCR cyclers. Programs for additional cyclers might be available on request.
 3. Start the program and data reading.
-

5. Result reading

1. Save the data at the end of the run.
 2. Read the channels for the fluorescence dyes FAMTM and ROXTM and show the 2nd deviation of the data.
 3. Read the calculation of the Ct values for the negative controls, the positive controls and the samples.
 4. Evaluate in accordance with the laboratory and instrument-specific reference ranges.
-

Reference range for selected qPCR devices:

qPCR device	Manufacturer	Evaluating the threshold cycle		
		positive	borderline	negative
LightCycler [®] 1.2	Roche Diagnostics	< 38	38 - 40	> 40
ABI Prism [®] 7500	Applied Biosystems	< 38	38 - 40	> 40
Rotor-Gene TM 6000	Corbett	< 38	38 - 40	> 40
Mx3005P [®]	Agilent Technologies	< 38	38 - 40	> 40

NOTES ON THE TEST PROCEDURE

1. This leaflet must be widely understood for a successful use of Onar[®]Ct. The reagents supplied should not be mixed with reagents from different lots and used as an integral unit. The reagents of the kit should not be used beyond its shelf life.
2. Any deviation from the test method can affect the results.
3. Inhibition may be caused by the sample matrix, but also by sample elution buffer of DNA extraction kits not recommended or validated. Negative controls should always be completed with the use sample elution buffer. Extract should be eluted in water or 10 mM Tris-HCl buffer, pH 7.2.
4. For each test setup, at least one negative control may be included also including sample preparation. Positive controls facilitate the evaluation of the test. Typical Ct values for the internal control and positive control are shown on the Certificate of Analysis and can be used for orienting quality control.
5. The use of control samples is advised to secure the day-to-day validity of results. The controls should be carried out in the same manner as the samples. It is recommended to run laboratory specific control samples with a high, medial and low (e.g. 3x LOD₉₅) level or establish commercial controls, for example, Minerva Biolabs Quantification Standard *Chlamydia trachomatis* (Order No. 52-0440). Participation in external quality control programs, such as offered by the INSTAND e.V., the Society for the promotion of quality assurance in medical laboratories, is recommended.
6. The high sensitivity of the PCR method increases the risk of false-positive results by cross-contamination or sloppy work. All work should be carried out in compliance with good laboratory practise.

INTERPRETATION OF RESULTS


As with all diagnostic tests the final clinical assessment should not be based on the results of a single test, but should be done by doctors only after the evaluation of all clinical and laboratory findings. The values given here are only a guide. Each laboratory should establish its own reference range for the Ct values of controls and samples.

The presence of *C. trachomatis* in the sample is indicated by an increasing fluorescence signal in the target channel. When using the internal control, an error-free PCR is indicated by an increasing fluorescence signal in the channel for internal control. The Internal control is present at a very low concentration that it is no longer amplified at high concentrations of target DNA. A successful PCR is then detected by the signal in the target channel.

Detection of <i>C. trachomatis</i> FAM [™] channel	Internal control ROX [™] channel	Interpretation
positive	irrelevant	<i>C. trachomatis</i> positive
negative	negative	PCR inhibition
negative	positive	<i>C. trachomatis</i> negative

ANALYTICAL CHARACTERISTICS OF THE TEST

A typical standard curve for Onar®Ct is shown below. The linear regression was determined from four individual dilution series.

Performance Parameter	Performance Characteristics
Analytical sensitivity including purification	0.9 IFU/ μ l, corresponding to 2.4 GU/ μ l
Analytical sensitivity without purification	2.5 IFU/ μ l, corresponding to 8.1 GU/ μ l
Analytical specificity	100 % for serovars A - K and L1 - L3 no cross reactivity with relevant commensals
Diagnostic sensitivity urines: n = 61 ejaculates n = 21 swabs n = 61	urines: 93.6 % ejaculates: 100 % swabs: 100 %
Diagnostic specificity n = 143	100 %
Precision (intra assay)	4 % at 2000 GU/PCR
Precision (inter assay)	6 % at 2000 GU/PCR
Robustness	cross contamination rate: 0 % fatal error rate: 0 % matrix effects: 0 %
Additional Parameters	lot-to-lot variations: 0.9 % linearity: 200 to 2.0×10^7 GU/PCR pipetting time: < 60 minutes without negative effect inter laboratory comparison: Instand RV531 completed

APPENDIX 1: PROGRAMMING AND DATA RECORDING OF DIFFERENT QPCR DEVICES

1. LightCycler® 1.2, 1.5 and 2.0

Program 1: Pre-incubation

Cycles	1
Analysis Mode	None
Temperature Targets [°C]	Segment 1
Target Temperature [°C]	95
Incubation time [min]	2:00
Temperature Transition Rate [°C/s]	20.0
Secondary Target Temperature [°C]	0
Step Size [°C]	0.0
Step Delay [Cycles]	0
Acquisition Mode	None

LightCycler® 2.0:

Please check the correct settings for the „seek temperature“, this should be adjusted on at least 90 °C.

Program 2: Amplification

Cycles	45
Analysis Mode	Quantification

Temperature Targets [°C]	Segment 1	Segment 2	Segment 3	Segment 4
Target Temperature [°C]	95	55	60	72
Incubation time [s]	0	5	7	5
Temperature Transition Rate [°C/s]	20.0	20.0	20.0	20.0
Secondary Target Temperature [°C]	0	0	0	0
Step Size [°C]	0.0	0.0	0.0	0.0
Step Delay [Cycles]	0	0	0	0
Acquisition Mode	None	None	Single	None

Program 3: Cooling

Cycles	1
Analysis Mode	None
Temperature Targets [°C]	Segment 1
Target Temperature [°C]	40
Incubation time [s]	30
Temperature Transition Rate [°C/s]	20.0
Secondary Target Temperature [°C]	0
Step Size [°C]	0.0
Step Delay [Cycles]	0
Acquisition Mode	None

Result Reading:

- Select the fluorescence channels F1 and F2 or Channel 1 and 3 at the LC 2.0
- Click on *Quantification* to generate the amplification plots and the specific *ct*-values
- The threshold will be generated automatically.
- Samples showing no significant increase in the amplification plot can be considered as negative.

	<i>Chl. trachomatis</i>	Internal control
channel for LC 1.2, 1.5	F1	F2
channel for LC 2.0	channel 1	channel 3

2. Rotor-Gene™ 6000 (5-plex)

Program Step 1: Pre-incubation

Setting	Hold
Hold Temperature	95 °C
Hold Time	3 min 0 sec

Please check the correct settings for the filter combination:

green filter (470-510): *Chl. trachomatis*

filter orange (585-610): **internal control**

Program Step 2: Amplification

Setting	Cycling
Cycles	45
Denaturation	95°C for 5 sec
Annealing	55°C for 10 sec —> acquiring to Cycling A (green and orange)
Elongation	72°C for 7 sec
Gain setting	automatic (auto Gain)
Slope Correct	activated
Ignore First	deactivated

Result Reading:

- Open the menu *Analysis*
- Select *Quantitation*
- Check the required filter set (green and orange) according to the following table and start data analysis by double click.
- The following windows will appear:
 - Quantitation Analysis - Cycling A* (green or orange)
 - Quant. Results - Cycling A* (green or orange)
 - Standard Curve - Cycling A* (green or orange)
- In window *Quantitation Analysis*, select first *linear scale* and then *slope correct*
- Threshold setup (not applicable if a standard curve was carried with the samples and auto threshold was selected)
 - In window *CT Calculation* set the threshold value to 0-1
 - Pull the threshold line into the graph. Adapt the threshold line to the initial linear section of the positive control reaction.
- The *ct*-values can be taken from the window *Quant. Results*.
- Samples showing no *ct*-value can be considered as negative.

3. ABI Prism® 7500

Detector Settings:

<i>Chl. trachomatis</i> Target Probe:	Reporter - FAM™	Quencher - none
Internal Control Probe:	Reporter - ROX™	Quencher - none

**The ROX™ Reference needs to be disabled. Activate both detectors for each well.
Measurement of fluorescence during extension.**

Program Step 1: Pre-incubation

Setting	Hold
Temperature	95 °C
Incubation time	3:00 min

Program Step 2: Amplification

Cycles	45
Setting	Cycle
Denaturing	95 °C for 30 sek
Annealing	55 °C for 30 sec
Extension	60 °C for 45 sec

Result Reading:

- Enter the following basic setting at the right task bar:
Data: Delta RN vs. Cycle
Detector: FAM and ROX

	<i>Chl. trachomatis</i>	Internal control
channel	FAM™	ROX™

Line Color: Well Color

- Open a new window with for the *Graph settings* by clicking the right mouse button
Select the following setting and confirm with *ok*:
Real Time Settings: Linear
Y-Axis Post Run Settings: Linear and Auto Scale
X-Axis Post Run Settings: Auto Scale
Display Options: 2
- Initiate the calculation of the ct-values and the graph generation by clicking on *Analyze* within the report window.
- Pull the threshold line into the graph. Adapt the threshold line to the initial linear section of the positive control reaction.
- Samples showing no ct-value can be considered as negative.

4. Mx3005P®

- Mark the wells in use in the setup menu, tab "Plate Setup"
- Select FAM™ und ROX™ for "Collect Fluorescence Data"
- Check that the Reference Dye function is set to "none" according to the ground settings of the device
- Specify the type of sample (negative control, positive control, standards, sample) by using the tab "well type"
- Program the temperature profile by using the tab "Thermal Profile Design":
 - Segment 1 (Pre-Melt): 3 min, 95 °C
 - Segment 2: 30 sec, 95 °C
30 sec, 55 °C - data collection end,
45 sec, 60 °C
45 cycles
- Start the program in the menu "Run" by clicking on "Start" in the window "Run Status"

Interpretation of the raw data:

- Select the "Analysis" menu, click the tab "Analysis Selection/Setup" and mark all positions for analysis
- Activate the following options in window "algorithm enhancement":
 - Amplification-based threshold
 - Adaptive baseline
 - Moving average
- Set the threshold to automatic in tab "Results" and "Amplification Plots"
- Read the *Ct values* in "Text Report"

APPENDIX 2: USE OF DIFFERENT REHYDRATION BUFFERS

Rehydration Buffer A

LightCycler® 1.2

LightCycler® 1.5

LightCycler® 2.0

LightCycler® 480

Rotor-Gene™ 3000

Rotor-Gene™ 6000

iCycler iQ®

iQ™ 5

Opticon 2

Chromo 4

MX3000P®

MX3005P®

MX4000®

Rehydration Buffer B

ABI Prism® 7000

ABI Prism® 7300+

ABI Prism® 7500

ABI Prism® 7700

ABI Prism® 7900

APPENDIX 3: LEGALS

Limited Product Warranty

This warranty limits our liability for replacement of this product. No warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided. Minerva Biolabs shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

Limited License

The use of this product for clinical diagnostics of pathogens is covered by a Scorpions license from DxS Ltd. Further information on Scorpions licenses can be obtained from DxS Ltd, 48 Grafton St, Manchester, UK M13 9XX.

Trademarks

LightCycler is a registered trademark of a member of the Roche Group. ABI Prism is a registered trademark of Applied Biosystems Corporation or its subsidiaries in the US and certain other countries. FAM™ and ROX™ are trademarks of Applied Biosystems Corporation or its subsidiaries in the US and certain other countries. Rotor-Gene is a trademark of Corbett Life Technologies. Mx3005P is a trademark of Agilent Technologies. Dacron® is a registered trademark of DuPont. Onar is registered trademarks of Minerva Biolabs GmbH.

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Tel. +49 (0)30 2000 437-0
support@minerva-biolabs.com

Made in Germany

Minerva Biolabs GmbH develops and manufactures products in accordance with EN ISO 9001:2008 and EN ISO 13485:2003 quality system requirement.
Reg.No. SY 60026567 0001 & SX 60025009 0001

...as precise as diagnostics should be™