**Onar®Ct**  
**Chlamydia trachomatis Diagnosis Kit for qPCR**  
**- Type 1 and Type 2 -**

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**Contents**

1. General Advice	13
1.1 Test Kit Components	13
1.2 Stability and Storage	13
1.3 Supplemental Requirements	13
2. Application and Test Principle	14
3. Test Protocol	14
3.1 Specimen Collection and Sample Storage	14
3.2 Preparation of Sample Material	16
3.3 Rehydration of the Reagents	16
3.4 Experimental Protocol and Result Reading	17
3.4.1 LightCycler® 1.2, 1.5 and 2.0 (for type 1 kits)	17
3.4.2 Rotorgene 6000 (5-plex) (for type 1 kits)	18
3.4.3 ABI 7500 (for type 2 kits)	19
3.5 PCR Master Mix Setup	20
3.6 Result Interpretation	20
4. Trouble Shooting	21
5. Performance Data	22

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# 1. Reagents and Materials

**Please do not use the product after the date of expiry. Discard remaining reagents.**

## 1.1 Test Kit Components

<i>Primer/Probe/Nucleotide Mix (PPN Mix)</i> , with type 1 kits only Primer set, probes and deoxynucleotide triphosphates dATP, dCTP, dGTP and dUTP, aliquoted for 25 reactions, lyophilized	red cap
<i>Primer/Probe/Nucleotide Mix (PPN Mix)</i> , with type 2 kits only Primer set, probes and deoxynucleotide triphosphates dATP, dCTP, dGTP and dUTP, aliquoted for 25 reactions, lyophilized	orange cap
<i>Rehydration Buffer</i> 1.8 ml	blue cap
<i>Positive Control DNA</i> Fragment of <i>Chlamydia trachomatis</i> DNA, non-infectious, lyophilized	green cap
<i>Internal Control DNA</i> Plasmid DNA, non-infectious, lyophilized	yellow cap
<i>PCR grade Water</i> 2 ml	white cap
<i>MB Taq DNA Polymerase</i> 1 U/ $\mu$ l	brown cap

## 1.2 Stability and Storage

Kit components are stable during shipping. Upon receipt, store at +2 °C to +8 °C. After receive the product take immediately the polymerase and store below – 20 °C. After rehydration of the *Primer/Probe/Nucleotide Mix* store below -18 °C and avoid repeated freezing and thawing. Protect the *Primer/Probe/Nucleotide Mix* from light. For repeated testing of low sample numbers, the *Primer/Probe/Nucleotide Mix* and the controls should be aliquoted. By following these recommendations, the kit is stable until the expiration date stated on the box label.

## 1.3 Supplemental Requirements

- qPCR machine
- corresponding PCR reaction tubes
- micro centrifuge, micropipettes and filtered tips (1 - 1000  $\mu$ l)

## 2. Application and Test Principle

Onar<sup>®</sup>Ct is an *in vitro* test for qualitative and semi-quantitative PCR-based diagnostics of *Chlamydia trachomatis* in clinical samples. The detection procedure can be performed within 2 to 3 hours. The primer set and probe are specific for a segment of the chromosomal *momp* gene (*major outer-membrane protein*). The genome of *Chlamydia trachomatis* contains only one copy of the *momp* gene locus. In contrast to other test methods amplifying the multicopy *cryptic plasmid*, Onar<sup>®</sup>Ct allows the relative quantification of the pathogen for significant clinical valuation and infection process monitoring.

The target probe emits fluorescent light at ~520 nm (type 1 and type 2). Onar<sup>®</sup>Ct also provides internal control DNA, which can be added to the reaction. When running the PCR with the internal control DNA, a successful performed reaction with a negative result is indicated by a distinct fluorescent signal. The internal control is detected by another probe, which emits light at ~610 nm (Type 1) and ~ 560 nm (Type 2).

Onar<sup>®</sup>Ct contains the nucleotide dUTP instead of dTTP. The heat-labile Uracil-DNA Glycosylase (UNG) is suitable to prevent carry-over contamination between PCRs. This technique relies on the incorporation of deoxyuridine triphosphate (dUTP) during all amplified reactions and the pretreatment of all successive PCR mixtures with the heat-labile UNG. UNG is not provided with Onar<sup>®</sup>Ct.

## 3. Test Protocol

### 3.1 Specimen Collection and Sample Storage

#### 1. Urine Specimens

- The patient should be advised not to urinate for at least two hours prior to specimen collection
- Collect 10 ml of first-catch urine in a clean container of polypropylene without preservatives.
- The urine specimens may be stored at room temperature (19 - 23 °C) if testing of samples can be completed within 24 hours.
- Samples have to be stored at 2 - 8 °C in the fridge, if the testing is done within 7 days.
- Stored at < - 20 °C for up to 30 days.

#### 2. Ejaculate Specimens

- The ejaculates must be collected the same day to be going to be shipped to the laboratory.
- Sample must be stored at < - 20 °C
- If a ejaculates cannot be processed directly after receipt in the laboratory, it must be stored at < - 20 °C and tested within 30 days.

### 3. Swab Specimens

- Samples must be shipped refrigerated.
- If the swab specimens are not processed directly after receipt in the laboratory, they have to be stored at 2 - 8 °C and have to be processed within seven days. If they are not tested within seven days after collection, the samples have to be stored at < -20 °C and tested within up to 30 days.

For the collection of eye, end cervical and urethral swabs, please use Darcon® or rayon swabs:

- ⇒ MASTASWAB, MAST DIAGNOSTICA, Cat. No. 800108 MD 508
- ⇒ Sterile Swab Applicators, micro Rheologics/Copan, Cat. No. 502CS01
- ⇒ Forensic swabs predetermined break zone, sterile wrapped, Copan, Cat. No. 3520CS01-4N6
- ⇒ Flocked swabs, microRheologics/Copan, Cat. No. 552C nylon

### 3.2 Preparation of Sample Material

Clinical samples should be extracted from infected or routine-screening patients. This diagnostic assay can be performed with various specimens, especially urine. It is necessary that templates for PCR analysis are prepared by DNA extraction using a suitable DNA extraction kit. The kit was validated using Minerva Biolabs DNA Extraction Kit (Cat. No. 56-1100) with urine and pool urine.

Alternatively, the following manual protocol has been proved of value and used for determination of analytical sensitivity:

1. centrifuge 1 ml urine or 1 ml transport medium (swabs) at 18.000 x g for 20 min.
2. removed supernatant
3. washed sediment in 500  $\mu$ l 1 x PBS pH 7.3 (phosphate-buffered saline)
4. centrifuge (18.000 x g for 20 min) and resuspend the sediment in 25  $\mu$ l 10 mM Tris-HCl, pH 7.2
5. add 2,5  $\mu$ l proteinase K (1 mg/ml) and incubate at 56 °C for 30 min
6. incubate at 95 °C for 10 min to lyse the cells and inactivate the proteinase K
7. after a short centrifugation, the sample is ready to us for PCR
8. if not tested immediately, samples are stored at -20 °C.

**Using extraction kits from other manufacturers please note that the extract must be eluted in water or 10 mM Tris-HCl buffer, pH 8.0. Elution buffers provided with the extraction kits from other manufacturers may cause inhibition of the PCR.**

10  $\mu$ l of the DNA extract can be used directly as PCR template. The extracts can be stored at a temperature of below -20 °C for a period of one year. Repeated freezing and thawing or storage in the refrigerator for more than 12 hours should be avoided. The sample should not contain more than 100  $\mu$ g/ml DNA.

### 3.3 Rehydration of the Reagents

1. centrifuge tubes with lyophilized components (5 sec at maximum speed)
2. add 365  $\mu$ l *Rehydration Buffer* to the *Primer/Probe/Nucleotide Mix*
4. add 300  $\mu$ l of *PCR grade Water* to the *Internal Control DNA* and *Positive Control DNA*
3. incubate for 5 minutes at room temperature
4. vortex and centrifuge again

### 3.4 Experimental Protocols and Result Reading

#### 3.4.1 LightCycler® 1.2, 1.5 and 2.0 (with type 1 kits)

##### Program 1: Pre-incubation

Cycles 1  
Analysis Mode None

##### **Temperature Targets [°C] Segment 1**

Target Temperature [°C] 95  
Incubation time [min] 2:00  
Temperature Transition Rate [°C/s] 20.0  
Secondary Target Temperature [°C] 0  
Step Size [°C] 0.0  
Step Delay [Cycles] 0  
Acquisition Mode None

##### **LightCycler® 2.0:**

**Please check the correct settings for the „seek temperature“, this should be adjusted on at least 90 °C.**

##### Program 2: Amplification

Cycles 45  
Analysis Mode Quantification

##### **Temperature Targets [°C] Segment 1 Segment 2 Segment 3 Segment 4**

Target Temperature [°C]	95	55	60	72
Incubation time [s]	0	5	7	5
Temperature Transition Rate [°C/s]	20.0	20.0	20.0	20.0
Secondary Target Temperature [°C]	0	0	0	0
Step Size [°C]	0.0	0.0	0.0	0.0
Step Delay [Cycles]	0	0	0	0
Acquisition Mode	None	None	Single	None

**Segment 3 is essential for the performance of the kit. Do not skip or modify this step.**

##### Program 3: Cooling

Cycles 1  
Analysis Mode None

##### **Temperature Targets [°C] Segment 1**

Target Temperature [°C] 40  
Incubation time [s] 30  
Temperature Transition Rate [°C/s] 20.0  
Secondary Target Temperature [°C] 0  
Step Size [°C] 0.0  
Step Delay [Cycles] 0  
Acquisition Mode None

Result Reading:

- Select the fluorescence channels F1 and F2 or Channel 1 and Channel 3 at the LC 2.0
- Click on *Quantification* to generate the amplification plots and the specific ct-values
- The threshold will be generated automatically.
- Samples showing no significant increase in the amplification plot can be considered as negative.

target	<b>Ct-momp</b>	<b>Internal control</b>
channel for LC 1.2, 1.5	F1 (530)	F2 (640)
channel for LC 2.0	channel 1 (530)	channel 3 (610)

3.4.2 Rotorgene 6000 (5-plex) (with type 1 kits)

Program Step 1: Pre-incubation

Setting	Hold
Hold Temperature	95 °C
Hold Time	3 min 0 sec

**Please check the correct settings for the filter combination:**

**green filter (470-510): C. trachomatis (Ct-momp)**  
**filter orange (585-610): internal control**

Program Step 2: Amplification

Setting	Cycling
Cycles	45
Denaturation	95°C for 5 sec
Annealing	55°C for 10 sec —> <b>acquiring to Cycling A (green and orange)</b>
Elongation	72°C for 7 sec
Gain setting	automatic (auto Gain)
Slope Correct	activated
Ignore First	deactivated

Result Reading:

- Open the menu *Analysis*
- Select *Quantitation*
- Check the required filter set (green or orange) according to the following table and start data analysis by double click.
- The following windows will appear:
  - Quantitation Analysis - Cycling A* (green or orange)
  - Quant. Results - Cycling A* (green or orange)
  - Standard Curve - Cycling A* (green or orange)
- In window *Quantitation Analysis*, select first *inear scale* and than *slope correct*

Threshold setup (not applicable if a standard curve was carried with the samples and auto threshold was selected)

- In window *CT Calculation* set the threshold value to 0-1
- Pull the threshold line into the graph. Adapt the threshold line to the initial linear section of the positive control reaction.
- The ct-values can be taken from the window *Quant. Results*.
- Samples showing no ct-value can be considered as negative.

target	<b>Ct-momp</b>	<b>Internal control</b>
channel	green	orange

### 3.4.3 ABI 7500 (with type 2 kits)

#### Program Step 1: Pre-incubation

Setting Hold  
Temperature 95 °C  
Incubation time 3:00 min

**Please check the correct settings of the detectors!**

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**The ROX reference needs to be disabled. Activate both detectors for each well. Measurement of fluorescence during extension.**

#### Result Reading:

- Enter the following basic setting at the right task bar:  
Data: Delta RN vs. Cycle  
Detector: FAM and VIC  
Line Color: Well Color
- Open a new window with for the *Graph settings* by clicking the right mouse button

target	<b>Ct-momp</b>	<b>Internal control</b>
channel	FAM	VIC

Select the following setting and confirm with *ok*:

Real Time Settings: Linear  
Y-Axis Post Run Settings: Linear and Auto Scale  
X-Axis Post Run Settings: Auto Scale  
Display Options: 2

- Initiate the calculation of the ct-values and the graph generation by clicking on *Analyze* within the report window.
- Pull the threshold line into the graph. Adapt the threshold line to the initial linear section of the positive control reaction.
- Samples showing no ct-value can be considered as negative.

### 3.5 PCR Master Mix Setup

Total volume per reaction is 25  $\mu\text{l}$ . When setting up reactions, calculations should also include positive (PC) and negative controls (NC). Pipette master mix on ice into a 1.5 ml reaction tube and mix gently.

#### Pipetting scheme:

	for 1 reaction	for 25 reactions*
<i>Primer/Probe/Nucleotide Mix</i>	14,0 $\mu\text{l}$	350,0 $\mu\text{l}$
<i>Internal Control DNA</i>	1,0 $\mu\text{l}$	25,0 $\mu\text{l}$
<i>MB DNA Polymerase</i>	1,0 $\mu\text{l}$	25,0 $\mu\text{l}$
total	16,0 $\mu\text{l}$	

\* the equivalent of the content of one red-capped vial

**The total duration from master mix preparation to PCR cycling must not exceed 45 minutes to avoid a decrease in the fluorescent signal.**

Aliquot 15  $\mu\text{l}$  of master mix into each PCR reaction tube. Discard the excess of master mix (1  $\mu\text{l}$ ).

Pipette 10  $\mu\text{l}$  of sample elution buffer or DNA-free water per tube as negative control. Seal the tubes. Add 10  $\mu\text{l}$  of sample per PCR reaction tube. Sealed these tubes as well before proceeding with the positive controls or standards. Add 10  $\mu\text{l}$  of the positive control to the remaining tubes as the last step in order to avoid cross contamination.

### 3.6 Result Interpretation

The presence of *C. trachomatis* in the sample is indicated by an increasing fluorescence signal in the *Ct-momp* channel during PCR. The concentration of the contaminant can be calculated by a software comparing the crossing cycle number of the sample with a standard curve created in the same run.

A successfully performed PCR without inhibition is indicated by an increasing fluorescence signal in the internal control channel, provided the Internal Control was added to the master mix. *C. trachomatis* DNA and Internal Control DNA are competitors in PCR. Because of the very low concentration of Internal Control in the PCR mix, the signal strength in this channel is reduced with increasing *C. trachomatis* DNA loads in the sample.

<b><i>Ct-momp</i> PCR</b>	<b>Internal Control</b>	<b>Interpretation</b>
positive	irrelevant	<i>C. trachomatis</i> positive
negative	negative	PCR inhibition
negative	positive	<i>C. trachomatis</i> negative

#### **4. Trouble shooting**

Before repeating a negative and a positive control run please check the cyclor protocol and pipetting scheme.

Insufficient (inhibition) or no amplification of the internal control may be due to the following reasons:

- programming mistake
- pipetting mistake
- sample elution buffer and sample used to complete the negative control are different

## 5. Performance Data

Performance Parameter	Performance Characteristics
Analytical sensitivity including purification	0.9 IFU/ $\mu$ l, corresponding to 2.4 GU/ $\mu$ l
Analytical sensitivity without purification	2.5 IFU/ $\mu$ l, corresponding to 8.1 GU/ $\mu$ l
Analytical specificity	100 % for serovars A - K and L1 - L3 no cross reactivity with relevant commensals
Diagnostic sensitivity urines: n = 61 ejaculates n = 21 swabs n = 61	urines: 93.6 % ejaculates: 100 % swabs: 100 %
Diagnostic specificity n = 143	100 %
Precision (intra assay)	4 % at 2000 GU/PCR
Precision (inter assay)	6 % at 2000 GU/PCR
Robustness	cross contamination rate: 0 % fatal error rate: 0 % matrix effects: 0 %
Additional Parameters	lot-to-lot variations: 0.9 % linearity: 200 to $2.0 \times 10^7$ GU/PCR pipetting time: < 60 minutes without negative effect inter laboratory comparison: Instand RV531 completed

## Related Products

### MB Taq DNA Polymerase

53-0050/-0100/-0200/-0250 MB Taq DNA Polymerase 50/100/200/250 units

### Clinical Diagnostic Kits for qPCR

20-2025/-2100/-2250	Venor <sup>®</sup> Mp-QP <i>Mycoplasma pneumoniae</i>	25/100/250 tests
21-2025/-2100/-2250	Onar <sup>®</sup> Ls-QP <i>Legionella species</i>	25/100/250 tests
21-3025/-3100/-3250	Onar <sup>®</sup> Lp-QP <i>Legionella pneumophila</i>	25/100/250 tests
22-2025/-2100/-2250	Onar <sup>®</sup> MRSA Methicillin-resistant <i>S. aureus</i>	25/100/250 tests
23-2025/-2100	Onar <sup>®</sup> Tp <i>Treponema pallidum</i>	25/100 tests
24-2025/-2100/-2250	Onar <sup>®</sup> Ct <i>Chlamydia trachomatis</i>	25/100/250 tests
25-2025/-2100/-2250	Onar <sup>®</sup> Pertussis <i>Bordetella pertussis/parapertussis</i>	25/100/250 tests

### Quantification Standards, 100 µl each, 1x10<sup>6</sup> genomes/µl

52-0101	<i>Legionella pneumophila</i> DNA Standard
52-0119	<i>Mycoplasma pneumoniae</i> DNA Standard

### Genomic DNA Extracts, 100 µl each, +/- 10 ng / 100 µl

51-0566	<i>Streptococcus pneumoniae</i> , DSMZ 20566
51-5571	<i>Bordetella pertussis</i> , DSMZ 5571
51-3415	<i>Bordetella parapertussis</i> , DSMZ 13415
51-0177	<i>Ureaplasma urealyticum</i> , NC 10177
51-0111	<i>Mycoplasma hominis</i> , NC 010111
51-0112	<i>Mycoplasma orale</i> , NC 010112
51-0117	<i>Mycoplasma fermentans</i> PG19, NC 010117
51-0119	<i>Mycoplasma pneumoniae</i> , NC 010119
51-0129	<i>Mycoplasma arginini</i> , NC 010129
51-0162	<i>Mycoplasma arthritidis</i> , NC 010162
51-0195	<i>Mycoplasma genitalium</i> , NC 010195
51-1746	<i>Mycoplasma penetrans</i> , NC 11746

### DNA Remover™

15-2025	DNA Decontamination Reagent, spray bottle	250 ml
15-2200	DNA Decontamination Reagent, refill bottles	4 x 500 ml

### Extraktion Kit

56-1100	MB DNA Extraction Kit	100 extractions
---------	-----------------------	-----------------