**Art. Nr. / Ordering No. 21-3025/-3100/-3250**

# **Onar<sup>®</sup>Lp**

***Legionella pneumophila Diagnostic Kit for qPCR***  
**version 2.0**

## **Gebrauchsinformation / Instructions for Use**

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Reagents for 25/100/250 reactions

## **Hersteller / Manufacturer:**

Minerva Biolabs GmbH, Koepenicker Strasse 325, 12555 Berlin, Germany

**ZUR IN VITRO DIAGNOSTIK**  
**FOR IN VITRO USE ONLY**

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## INDICATION

The Onar<sup>®</sup>Lp test system is an in vitro test for quantitative diagnosis of *Legionella pneumophila* in clinical samples.

## EXPLANATION OF THE TEST

*Legionella pneumophila* is ubiquitous in surface water and moist soil and is known as an opportunistic pathogen of Legionnaires' disease and Pontiac fever. The incubation period is 20-10 days. The airborne transmission is unlikely. Typical symptoms of Legionnaires' disease include high fever, cough, chest pain, diarrhea, confusion, some with severe cases.

Onar<sup>®</sup>Lp is an in vitro assay for the quantitative detection of *Legionella pneumophila* in clinical specimens. Compared to immunological methods that can detect a primary atypical pneumonia first 2 to 3 weeks after infection, the PCR-based detection of pathogens has proven especially helpful in the early diagnosis. Compared with the detection by culture, PCR is characterized by a shorter analysis time and higher sensitivity.

The direct detection of *Legionella pneumophila* is possible within 1-2 hours. The kits were designed under the aspect of high practicality and ease of use. This ensures high accuracy and reproducibility of the assay. Extensive studies demonstrate the high specificity and sensitivity of Onar<sup>®</sup>Lp, required for the early identification of *Legionella pneumophila* infections.

## TEST PRINCIPLE

The amplified *mip* gene is highly conserved within the species *Legionella*. It allows detection of all *L. pneumophila* serogroups 1-14 and 16. Cross activities to commensals of the oropharynx and other pathogens in the respiratory tract, eg bordetella, chlamydia, mycoplasma, mycobacteria, pseudomonas and streptococci are not known. The detection limit is 10 genome equivalents per PCR. Because of the wide linear range reliable results without a time-consuming conditioning of the sample material are possible with Onar<sup>®</sup>Lp. The samples are compared to the antigen test significantly more stable.

Onar<sup>®</sup>Lp contains the nucleotide dUTP instead of dTTP. The heat-labile Uracil-DNA Glycosylase (UNG) is suitable to prevent carry-over contamination between PCRs. This technique relies on the incorporation of deoxyuridine triphosphate (dUTP) during all amplified reactions and the pretreatment of all successive PCR mixtures with the heat-labile UNG. UNG is not provided with Onar<sup>®</sup>Lp.

## REAGENTS

Each kit contains reagents for 25, 100 or 250 reactions. The expiry date of the unopened package is marked on the package label. The kit components are stored until use at +2 to +8 °C. The lot specific Certificate of Analysis can be downloaded from our website ([www.minerva-biolabs.com](http://www.minerva-biolabs.com)).

Kit Component Label Information	Quantity			Cap Color
	25 reactions Order No. 21-3025	100 reactions Order No. 21-3100	250 reactions Order No. 21-3250	
Lp Mix	1 x lyophilized	4 x lyophilized	10 x lyophilized	red
Rehydration Buffer	1 x 1.8 ml	1 x 1.8 ml	2 x 1.8 ml	blue
Positive Control DNA	1 x lyophilized	1 x lyophilized	1 x lyophilized	green
Internal Control DNA	1 x lyophilized	1 x lyophilized	2 x lyophilized	yellow
PCR grade Water	1 x 2.0 ml	1 x 2.0 ml	1 x 2.0 ml	white

## NEEDED, BUT NOT INCLUDED IN THE KIT

The Onar<sup>®</sup>Lp kit contains the reagents for the specific detection of *Legionella pneumophila*. General industrial supplies and reagents, usually available PCR laboratories are not included:

- qPCR device with filter sets for the detection of the fluorescence dyes FAM<sup>™</sup> and ROX<sup>™</sup>
- PCR reaction tubes for the specific qPCR device
- 1,5 ml reaction tubes, DNA- and RNA-free
- Micro centrifuge for 1.5 ml PCR reaction tubes
- Pipettes with corresponding filter tips to prepare and dispense the reaction mix (10, 100 und 1000 µl)
- Optional: For calibration we recommend our *Legionella pneumophila* DNA Quantification Standard (Cat.-No. 52-0101).

## SPECIMEN

To avoid inhibition of PCR, samples should be appropriately prepared. The assay was validated using the following extraction procedures for swab and sputum:

- Minerva Biolabs, MB DNA Extraction Kit (No. 56-1100)
- Qiagen, QIAamp DNA Mini Kit (No. 51 304), Protocol for Swab Samples

Good experiences were made with the direct use of swabs squeezed in an appropriate buffer solution. The swab should not be stored for longer than 6 hours without and not longer than 48 hours with refrigeration. Particularly suitable are dry swab from Dacron® or synthetic fibers

- |    |   |
|----|---|
| 1. | Pipette 200 µl wash buffer (10 mM Tris.HCl pH 8.5, 1 mM EDTA, 0.5 % Tween20) or Swab Wash Buffer (Minerva Biolabs, Order No. 57-1100) into a 1.5 ml reaction tube |
| 2. | Stick the swab into the buffer, incubate for 5 min at room temperature  |
| 3. | Squeeze out the buffer by pressing the swab against the inner wall of the reaction tube   |
| 4. | Use the eluate directly for PCR   |

The quality of the smear taken by the sender and the type of swab used may have uncontrollable influences on the amount of accompanied substances and thus on the quality of the sample material. It is not possible to validate this practicable method universally and in appropriate detail. DNA extracts and eluates can be kept for 6 days at +2 to +8 °C. Longer storage times require a temperature of -18 °C. Repeated freezing and thawing should be avoided.

#### PRECAUTIONS

For in vitro diagnostic use. This kit should be used only by trained persons. All samples should be considered potentially infectious and handled at the local or national regulations. This kit does not contain hazardous substances and may be disposed of according to local regulations.

#### TEST PROCEDURE

The test should be carried out with negative and positive controls and patient samples in duplicate. For quantification, a dilution series of an appropriate standard (Cat. No. 52-0101) is kept. All reagents and samples must be equilibrated to 2 to 8 °C prior use.

#### 1. Rehydration of the Reagents

After reconstitution, the reagents should be stored at 2 to 8 °C for up to 4 days or below -18 °C till the expiration date stated on the box label. Repeated freezing and thawing should be avoided and reconstituted controls (internal amplification control and the positive control) stored in aliquots.

1.	<i>Lp</i> Mix Internal Control Positive Control	red cap yellow cap green cap	Spin all lyophilized components for 5 sec at maximum speed of the mini centrifuge
2.	<i>Lp</i> Mix	red cap	Add 365 µl rehydration buffer (blue cap)
3.	Internal Control	yellow cap	Add 300 µl PCR grade Water (white cap)
4.	Positive Control	green cap	Add 300 µl PCR grade Water (white cap)
5.	<i>Lp</i> Mix Internal Control Positive Control	red cap yellow cap green cap	Incubate 5 min at room temperature
6.	<i>Lp</i> Mix Internal Control Positive Control	red cap yellow cap green cap	Vortex DNA briefly and spin for 5 sec

#### 2. Preparation of the Reaction Mix

Prepare the required amount of master mix at room temperature in a 1.5 ml reaction tube for all control and test reactions.			
1.		for 1 reaction	for 25 reactions
	<i>Lp</i> Mix	15.0 µl	350.0 µl
	Internal Control DNA	1.0 µl	25.0 µl
2.	Homogenize the reaction mix by carefully snapping the tube		
3.	Add 15 µl to each PCR tube, discard remaining material.		

### 3. Loading the test tubes

- Negative Control:
1. Add 10 µl elution buffer from DNA extraction (ref. Chapter sample material) or PCR Grade Water (white cap).
  2. Sample: Add 10 µl of each sample.
  3. Positive Control: Add 10 µl Positive Control (green cap).
  4. Spin all PCR tubes briefly, load the qPCR cyclers and start the program.

Preparation of the master mix and sample loading should not take more than 45 minutes to avoid a reduction in the fluorescent signal. The pipetting sequence should be respected and the tubes closed after each sample load.

### 4. Starting the reaction

1. Load the cycler, check each PCR tube and the cycler lid for tight fit.
2. Program the qPCR cycler or check stored temperature profiles. See Appendix for temperature profiles of selected qPCR cyclers. Programs for additional cyclers might be available on request.
3. Start the program and data reading.

### 5. Result reading

1. Save the data at the end of the run.
2. Read the channels for the fluorescence dyes FAM™ and ROX™ and show the 2nd deviation of the data.
3. Read the calculation of the Ct values for the negative controls, the positive controls and the samples.
4. Evaluate in accordance with the laboratory and instrument-specific reference ranges.

#### Reference range for selected qPCR devices:

qPCR device	Manufacturer	Evaluating the threshold cycle		
		positive	borderline	negative
LightCycler® 1.2	Roche Diagnostics	< 40	≥ 40	no Ct
ABI Prism® 7500	Applied Biosystems	< 40	≥ 40	no Ct
RotorGene® 6000	Corbett Research	< 40	≥ 40	no Ct
Mx3005P®	Agilent Technologies	< 40	≥ 40	no Ct

### NOTES ON THE TEST PROCEDURE

1. This leaflet must be widely understood for a successful use of Onar®Lp. The reagents supplied should not be mixed with reagents from different lots and used as an integral unit. The reagents of the kit should not be used beyond its shelf life.
2. Any deviation from the test method can affect the results.
3. Inhibition may be caused by the sample matrix, but also by sample elution buffer of DNA extraction kits not recommended or validated. Negative controls should always be completed with the use sample elution buffer.
4. For each test setup, at least one negative control may be included also including sample preparation. Positive controls facilitate the evaluation of the test. Typical Ct values for the internal control and positive control are shown on the Certificate of Analysis and can be used for orienting quality control.
5. The use of control samples is advised to secure the day-to-day validity of results. The controls should be carried out in the same manner as the samples. It is recommended to run laboratory specific control samples with a high, medial and low (e.g. 3x LOD<sub>95</sub>) level or establish commercial controls, for example, Minerva Biolabs Quantification Standard *Legionella pneumophila* (Order No. 52-0101). Participation in external quality control programs, such as offered by the INSTAND e.V., the Society for the promotion of quality assurance in medical laboratories, is recommended.

### INTERPRETATION OF RESULTS

As with all diagnostic tests the final clinical assessment should not be based on the results of a single test, but should be done by doctors only after the evaluation of all clinical and laboratory findings. The values given here are only a guide. Each laboratory should establish its own reference range for the Ct values of controls and samples.

The presence of *Legionella pneumophila* in the sample is indicated by an increasing fluorescence signal in the target channel. When using the internal control, an error-free PCR is indicated by an increasing fluorescence signal in the channel for internal control. The Internal control is present at a very low concentration that it is no longer amplified at high concentrations of target DNA. A successful PCR is then detected by the signal in the target channel.

Detection of <i>Legionella pneumophila</i> FAM™ channel	Internal control ROX™ channel	interpretation
positive	irrelevant	<i>Legionella pneumophila</i> positive
negative	negative	PCR inhibition
negative	positive	<i>Legionella pneumophila</i> negative


## ANALYTICAL CHARACTERISTICS OF THE TEST

The precision of the Onar<sup>®</sup>Lp was determined using genomic DNA preparations (*Legionella pneumophila*) at different concentrations generated on a RotorGene<sup>®</sup> 6000 instrument.

The intra-assay variation shown is based on the scattering of the Ct values of samples of the same concentrations within a test batch, each with four replicates per concentration. Given is the mean coefficient of variation (CV) from five repetitions. The inter-assay variation shown is based on the scattering of the Ct values from samples of equal concentrations of two independent assays, each with four replicates per concentration. Given is the average of the retail price of 3 reps.

Concentration [GU/PCR]	Intraassay Variation VC [%]	Interassay Variation VC [%]	Lot Variation VC [%]
10	2,51	5,27	4,08
100	0,90	1,97	2,97
1000	0,37	1,73	2,80
10000	0,46	1,86	2,93
100000	0,45	2,12	3,50
1000000	0,29	2,14	4,87
10000000	0,33	1,21	2,07

A typical standard curve for Onar<sup>®</sup>Lp is shown below. The linear regression was determined from four replicates per dilution level (1 to 10<sup>6</sup> GU/PCR) with the RotorGene<sup>®</sup> 6000.

Analytical Sensitivity: 2 GU/μl

Analytical Specificity: 100 % for *Legionella pneumophila* sero group 1-14 and 16

No cross reactivity with other *Legionella* species, relevant commensals or DNA of human origin

Linearity 50 to 10<sup>6</sup> GU/PCR

Pipetting time < 60 minutes without influence

Storage of components after rehydration at 2 to 8 °C up to 4 days without effect

## Appendix: Programming and data recording devices of different qPCR

### 1. LightCycler<sup>®</sup> 1.2, 1.5, 2.0 and 480

#### Program 1: Pre-incubation

Cycles	1
Analysis Mode	None
<b>Temperature Targets</b>	<b>Segment 1</b>
Target Temperature [°C]	95
Incubation time [min]	2:00
Temperature Transition Rate [°C/s]	20.0
Secondary Target Temperature [°C]	0
Step Size [°C]	0.0
Step Delay [Cycles]	0
Acquisition Mode	None

#### LightCycler<sup>®</sup> 2.0:

**Please check the correct settings for the „seek temperature“, this should be adjusted on at least 90°C.**

#### Program 2: Amplification

Cycles	45			
Analysis Mode	Quantification			
<b>Temperature Targets</b>	<b>Segment 1</b>	<b>Segment 2</b>	<b>Segment 3</b>	<b>Segment 4</b>
Target Temperature [°C]	95	55	60	72
Incubation time [s]	0	5	7	5
Temperature Transition Rate [°C/s]	20.0	20.0	20.0	20.0
Secondary Target Temperature [°C]	0	0	0	0
Step Size [°C]	0.0	0.0	0.0	0.0
Step Delay [Cycles]	0	0	0	0
Acquisition Mode	None	None	Single	None

#### Program 3: Cooling

Cycles	1
Analysis Mode	None
<b>Temperature Targets</b>	<b>Segment 1</b>
Target Temperature [°C]	40
Incubation time [s]	30
Temperature Transition Rate [°C/s]	20.0
Secondary Target Temperature [°C]	0
Step Size [°C]	0.0
Step Delay [Cycles]	0
Acquisition Mode	None

#### Result Reading:

- Select the fluorescence channels F1 and F2 or Channel 1 and 3 at the LC 2.0
- Click on *Quantification* to generate the amplification plots and the specific ct-values
- The threshold will be generated automatically.
- Samples showing no significant increase in the amplification plot can be considered as negative.

target	<i>L. pneumophila</i>	Internal control
channel for LC 1.2, 1.5	F1 (530)	F2 (640)
channel for LC 2.0, 480	channel 1 (530)	channel 3 (610)

## 2. RotorGene® 6000 (5-plex)

#### Program Step 1: Pre-incubation

Setting	Hold
Hold Temperature	95 °C
Hold Time	3 min 0 sec

**Please check the correct settings for the filter combination:**  
**green filter (470-510): *Legionella pneumophila***  
**filter orange (585-610): internal control**

#### Program Step 2: Amplification

Setting	Cycling
Cycles	45
Denaturation	95°C for 5 sec
Annealing	55°C for 10 sec —> <b>acquiring to Cycling A (green and orange)</b>
Elongation	72°C for 7 sec
Gain setting	automatic (auto Gain)
Slope Correct	activated
Ignore First	deactivated

#### Result Reading:

- Open the menu *Analysis*
- Select *Quantitation*
- Check the required filter set (green and orange) according to the following table and start data analysis by double click.
- The following windows will appear:  
*Quantitation Analysis - Cycling A (green or orange)*  
*Quant. Results - Cycling A (green or orange)*  
*Standard Curve - Cycling A (green or orange)*
- In window *Quantitation Analysis*, select first *linear scale* and than *slope correct*  
Threshold setup (not applicable if a standard curve was carried with the samples and auto threshold was selected)
  - In window *CT Calculation* set the threshold value to 0-1
  - Pull the threshold line into the graph. Adapt the threshold line to the initial linear section of the positive control reaction.
- The ct-values can be taken from the window *Quant. Results*.
- Samples showing no ct-value can be considered as negative.

## 3. ABI Prism® 7500

#### Detector Settings:

*Legionella* Target Probe: Reporter - FAM™, Quencher - none  
Internal Control Probe: Reporter - ROX™, Quencher - none

**The ROX™ Reference needs to be disabled. Activate both detectors for each well.**  
**Measurement of fluorescence during extension.**

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#### Result Reading:

- Enter the following basic setting at the right task bar:

Data:	Delta RN vs. Cycle	Target	<i>L. pneumophila</i>	Internal Control
Detector:	FAM™ and ROX™	Channel	FAM™	ROX™
Line Color:	Well Color			

- Open a new window with for the *Graph settings* by clicking the right mouse button  
Select the following setting and confirm with *ok*:  
Real Time Settings: Linear  
Y-Axis Post Run Settings: Linear and Auto Scale  
X-Axis Post Run Settings: Auto Scale  
Display Options: 2
- Initiate the calculation of the ct-values and the graph generation by clicking on *Analyze* within the report window.
- Pull the threshold line into the graph. Adapt the threshold line to the initial linear section of the positive control reaction.
- Samples showing no ct-value can be considered as negative.

#### 4. Mx3005P®

- Mark the wells in use in the setup menu, tab “Plate Setup”
- Select FAM™ und ROX™ for “Collect Fluorescence Data”
- Check that the Reference Dye function is set to “none” according to the ground settings of the device
- Specify the type of sample (negative control, positive control, standards, sample) by using the tab “well type”
- Program the temperature profile by using the tab “Thermal Profile Design”:
  - Segment 1 (Pre-Melt): 2 min, 95 °C
  - Segment 2:
    - 30 sec, 95 °C
    - 30 sec, 55 °C - data collection end
    - 45 sec, 60 °C - data collection end
    - 45 cycles
- Start the program in the menu “Run” by clicking on “Start” in the window “Run Status”

#### Interpretation of the raw data:

- Select the “Analysis” menu, click the tab “Analysis Selection/Setup” and mark all positions for analysis
- Activate the following options in window “algorithm enhancement”:
  - Amplification-based threshold
  - Adaptive baseline
  - Moving average
- Set the threshold to automatic in tab “Results” and “Amplification Plots”
- Read the *Ct values* in “Text Report”

## APPENDIX

### Limited Product Warranty

This warranty limits our liability for replacement of this product. No warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided. Minerva Biolabs shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

### Limited License

The use of this product for clinical diagnostics of pathogens is covered by a Scorpions license from DxS Ltd. Further information on Scorpions licenses can be obtained from DxS Ltd, 48 Grafton St, Manchester, UK M13 9XX.

### Trademarks

LightCycler is a registered trademark of a member of the Roche Group. ABI Prism is a registered trademark of Applied Biosystems Corporation or its subsidiaries in the US and certain other countries. FAM™ and ROX™ are trademarks of Applied Biosystems Corporation or its subsidiaries in the US and certain other countries. Mx3005P is a trademark of Agilent Technologies. RotorGene is a registered trademark of Qiagen GmbH. Dacron® is a registered trademark of DuPont. Onar is registered trademarks of Minerva Biolabs GmbH.

## Related Products

### MB Taq DNA Polymerase

53-0050/-0100/-0200/-0250	MB Taq DNA Polymerase (5 U/ $\mu$ l)	50/100/200/250 units
53-1050/-1100/-1200/-1250	MB Taq DNA Polymerase (1 U/ $\mu$ l)	50/100/200/250 units

### Clinical Diagnostic Kits for qPCR

20-2025/-2100/-2250	Onar <sup>®</sup> Mp <i>Mycoplasma pneumoniae</i>	25/100/250 tests
21-3025/-3100/-3250	Onar <sup>®</sup> Lp <i>Legionella pneumophila</i>	25/100/250 tests
22-2025/-2100/-2250	Onar <sup>®</sup> MRSA Methicillin-resistant <i>S. aureus</i>	25/100/250 tests
23-2025/-2100	Onar <sup>®</sup> Syphilis <i>Treponema pallidum</i>	25/100 tests
24-2025/-2100/-2250	Onar <sup>®</sup> Ct <i>Chlamydia trachomatis</i>	25/100/250 tests
25-2025/-2100/-2250	Onar <sup>®</sup> Pertussis <i>Bordetella pertussis/parapertussis</i>	25/100/250 tests

### Quantification Standards, 100 $\mu$ l each, 1x10<sup>8</sup> genomes/ $\mu$ l

52-0101	<i>Legionella pneumophila</i> DNA Standard
52-0119	<i>Legionella pneumophila</i> DNA Standard

### Genomic DNA Extracts, 100 $\mu$ l each, +/- 10 ng / 100 $\mu$ l

51-0566	<i>Streptococcus pneumoniae</i> , DSMZ 20566
51-5571	<i>Bordetella pertussis</i> , DSMZ 5571
51-3415	<i>Bordetella parapertussis</i> , DSMZ 13415
51-0177	<i>Ureaplasma urealyticum</i> , NC 10177
51-0111	<i>Mycoplasma hominis</i> , NC 010111
51-0112	<i>Mycoplasma orale</i> , NC 010112
51-0117	<i>Mycoplasma fermentans</i> PG19, NC 010117
51-0119	<i>Mycoplasma pneumoniae</i> , NC 010119
51-0129	<i>Mycoplasma arginini</i> , NC 010129
51-0162	<i>Mycoplasma arthritidis</i> , NC 010162
51-0195	<i>Mycoplasma genitalium</i> , NC 010195
51-1746	<i>Mycoplasma penetrans</i> , NC 11746

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Tel. +49 (0)30 2000 437-0  
support@minerva-biolabs.com

## Made in Germany

Minerva Biolabs GmbH develops and manufactures products in accordance with EN ISO 9001:2008 and EN ISO 13485:2003 quality system requirement. Reg.No. SY 60026567 0001 & SX 60025009 0001

...as precise as diagnostics should be™