**Art. Nr. / Order No. 11-4025/-4100/-4250 (Type 1)**

**Art. Nr. / Order No. 11-5025/-5100/-5250 (Type 2)**

# Venor<sup>®</sup>GeM-qEP

***Mycoplasma Detection Kit for qPCR***

**version 1.2**

## **Gebrauchsinformation / Instructions for Use**

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Reagents for 25/100/250 reactions

## **Hersteller / Manufacturer:**

Minerva Biolabs GmbH, Koepenicker Strasse 325, 12555 Berlin, Germany

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#### INDICATION

Venor®GeM-qEP is used for the detection of *Mollicutes* (*Mycoplasma*, *Acholeplasma*, *Spiroplasma*) contamination in cell cultures and other cell culture derived biologicals.

#### EXPLANATION OF THE TEST

The kit was validated according to the *European Pharmacopoeia* (EP) with different kinds of sample material (chondrocytes, serum, cell culture supernatant, etc.). The test has been validated for both methods described in the *European Pharmacopoeia* 2.6.7: direct detection and cell-culture enrichment followed by PCR. The kit complies fully with the requirements of EP 2.6.7.

#### TEST PRINCIPLE

Venor®GeM-qEP utilizes the polymerase chain reaction (PCR), which was established as the method of choice for high sensitivity. The kit includes three different Primer/Probe/Nucleotide mixes. These mixes contain FAM labeled Scorpions-probes specific for different *mycoplasma*. The *Species Mix* detects a broad range of *Mycoplasma* (see detailed list), the *Mp Mix* detects *Mycoplasma pneumoniae* and the *Al Mix* is specific for *Acholeplasma laidlawii*.

False negative results due to PCR inhibitors or improper DNA extraction are detected by the internal amplification control. The Internal Control DNA can be added to the sample prior to DNA extraction and analysis for verification of the complete process (DNA extraction and PCR reaction). The Internal Control DNA can also be added directly to the PCR master mix to act as a PCR control only. The amplification of the control reaction is detected at 610 nm (ROX channel, type 1) or at 560 nm (VIC channel, type 2) and the pathogen-specific sequence at 520 nm (FAM channel).

The kit contains dUTP instead of dTTP, so the option is available to degrade amplicons from previous analysis by use of uracil-DNA glycosylase (UNG). Thus the occurrence of false-positive result can be minimized. UNG is not part of Venor®GeM-qEP.

## REAGENTS

Each kit contains reagents for 25, 100 or 250 reactions. The expiry date of the unopened package is marked on the package label. The kit components are stored until use at +2 to +8 °C. The lot specific *Certificate of Analysis* can be downloaded from our website ([www.minerva-biolabs.com](http://www.minerva-biolabs.com)).

Kit Component Label Information	Quantity			Cap Color
	25 reactions Order No. 11-4025 Order No. 11-5025	100 reactions Order No. 11-4100 Order No. 11-5100	250 reactions Order No. 11-4250 Order No. 11-5250	
Species Mix	1 x lyophilized	4 x lyophilized	10 x lyophilized	red
Mp Mix	1 x lyophilized	4 x lyophilized	10 x lyophilized	black
AI Mix	1 x lyophilized	4 x lyophilized	10 x lyophilized	brown
Rehydration Buffer	1 x 1.8 ml	3 x 1.8 ml	7 x 1.8 ml	blue
Positive Control DNA	1 x lyophilized	1 x lyophilized	2 x lyophilized	green
Internal Control DNA	1 x lyophilized	1 x lyophilized	2 x lyophilized	yellow
PCR grade Water	1 x 2.0 ml	1 x 2.0 ml	1 x 2.0 ml	white

## NEEDED, BUT NOT INCLUDED IN THE KIT

The Venor®GeM-qEP kit contains reagents for the specific detection of *Mollicutes*. General industrial supplies and reagents, usually available PCR laboratories are not included:

- qPCR device with filter sets for the detection of the fluorescence dyes FAM and ROX
  - PCR reaction tubes for the specific qPCR device
  - 1,5 ml reaction tubes, DNA- and RNA-free
  - Micro centrifuge for 1.5 ml PCR reaction tubes
  - Pipettes with corresponding filter tips to prepare and dispense the reaction mix (10, 100 und 1000 µl)
  - hot-start DNA polymerase (1 unit/reaction; 3 units/sample) with a concentration of 1 U/µl or 5 U/µl
- optional:
- DNA Quantification Standard available for all EP listed species (Minerva Biolabs)

This test provides excellent results with MB Taq DNA polymerase (Cat # 53-0050/100/200/250). The polymerase is not part of the kit. A suitable performance of other polymerases can not be guaranteed. Please request your free samples of MB Taq polymerase (10 units).

## SPECIMEN

Samples should be derived from cultures which are at 90-100 % confluence. PCR inhibiting substances may accumulate in the medium of older cultures. For these sample materials a DNA extraction (e.g. MB DNA Extraction Kit Cat # 56-1100) is strictly recommended prior testing. Penicillin and streptomycin in the culture media do not inhibit mycoplasma or affect test sensitivity. Only cell culture supernatant should be applied to test for mycoplasma. Cell pellets should only be tested after suitable DNA extraction, since debris will interfere with the PCR reaction. With average titers at 10<sup>6</sup> and a maximum titer at 10<sup>8</sup> you will find sufficient mycoplasma in the supernatant to guarantee a sensitive PCR. However, other materials that can be tested are Fetal Calf Serum, vaccines, and paraffin-embedded samples following DNA extraction. If necessary, templates

for PCR analysis are prepared by DNA extraction using commercially available extraction kits (MB DNA Extraction Kit Cat. No. 56-1100). Please be sure to remove any alcohol containing wash buffer from the preparation to avoid co-elution of alcohol and sample material. Any remaining alcohol may inhibit the PCR. 2 µl of the extract can be used directly as PCR template. To avoid false positive results, we recommend the use of the PCR grade water delivered with the kit, aerosol-preventive filter tips and gloves.

The preparation of sample material could be performed by one of the following methods:

### Heat-inactivation of the sample material

The templates for the PCR analysis are prepared by direct heating of the cell culture supernatant or the biological sample material:

- 100 µl liquid supernatant of the sample material is transferred into a sterile reaction tube;
- the supernatant is incubated at 95 °C for 10 minutes;
- the supernatant is centrifuged briefly (5 seconds, 1000 x g) to remove cellular debris.
- The supernatant is used in the PCR. Alternatively, the DNA can be purified with a commercial extraction kit (e.g. MB DNA Extraction Kit Cat # 56-1100).

### Enrichment of mycoplasma by centrifugation

- 1 ml liquid supernatant of the sample material is transferred into a sterile reaction tube;
- the supernatant is centrifuged (15 minutes, 10.000 x g) to sediment mycoplasma particles. Alternatively: centrifuge the supernatant 6 min with 13.000 x g.
- The supernatant is discarded and the pellet is suspended into 50 µl buffer (10 mm Tris, pH 8.4).
- The sample should be vortexed and finally heated up to 95 °C for 10 min.

The extracts can be stored at a temperature of at least -18 °C for a period of one year. Repeated freezing and defrosting, or storage in the refrigerator for longer than 12 hours should be avoided. The sample should not contain more than 100 µg/ml DNA.

## PRECAUTIONS

For in vitro diagnostic use. This kit should be used only by trained persons. All samples should be considered potentially infectious and handled at the local or national regulations. This kit does not contain hazardous substances and may be disposed of according to local regulations.

## TEST PROCEDURE

The test should be carried out with negative and positive controls and samples in duplicate. For quantification, a dilution series of an appropriate standard is kept. All reagents and samples must be equilibrated to +2 to +8 °C prior use.

## 1. Rehydration of the reagents

After reconstitution, the reagents should be stored at below -18 °C. Repeated freezing and thawing should be avoided and reconstituted controls (internal amplification control and the positive control) stored in aliquots.

1.	Species Mix Mp Mix AI Mix Internal Control Positive Control	red cap black cap brown cap yellow cap green cap	Spin all lyophilized components for 5 sec at maximum speed of the centrifuge
2.	Species Mix Mp Mix AI Mix	red cap black cap brown cap	Add 365 µl rehydration buffer (blue cap)
3.	Internal Control	yellow cap	Add 300 µl PCR grade Water (white cap)
4.	Positive Control	green cap	Add 300 µl PCR grade Water (white cap)
5.	Species Mix Mp Mix AI Mix Internal Control Positive Control	red cap black cap brown cap yellow cap green cap	Incubate 5 min at room temperature
6.	Species Mix Mp Mix AI Mix Internal Control Positive Control	red cap black cap brown cap yellow cap green cap	Vortex DNA briefly and spin for 5 sec

## 2. Preparation of the Reaction Mix

- Prepare the required amount of master mix individually for each species probe. Mix the kit components at room temperature in a 1.5 ml reaction tube for all control and test reactions.

Pipetting scheme using a polymerase with a concentration of 5 U/µl:

	Color code	for 1 reaction	for 25 reactions
Species or <i>Mp</i> or <i>AI</i> Mix	red/black/brown	14.0 µl	350.0 µl
PCR grade water	white	7.8 µl	195.0 µl
Internal Control DNA	yellow	1.0 µl	25.0 µl
Polymerase (5 U/µl)		0.2 µl	5.0 µl

Pipetting scheme using a polymerase with a concentration of 1 U/µl:

	Color code	for 1 reaction	for 25 reactions
Species or <i>Mp</i> or <i>AI</i> Mix	red/black/brown	14.0 µl	350.0 µl
PCR grade water	white	7.0 µl	180.0 µl
Internal Control DNA	yellow	1.0 µl	25.0 µl
Polymerase (5 U/µl)		1.0 µl	25.0 µl

- Homogenize the reaction mix by carefully snapping the tube
- Add 23 µl to each PCR tube, discard remaining material.

## 3. Loading the test tubes

- Negative Controls: add 2 µl elution buffer from DNA extraction (ref. chapter Sample Material) or PCR Grade Water (white cap)
- Sample: add 2 µl of sample
- Positive Control: add 2 µl Positive Control DNA (green cap)
- Close and spin all PCR tubes briefly, load the qPCR cyclers and start the program

Preparation of the master mix and sample loading should not take more than 45 minutes to avoid a reduction in the fluorescent signal. The pipetting sequence should be respected and the tubes closed after each sample load. If purified DNA in aqueous solution is tested sensitivity of the test system can be further increased by using a volume of up to 10 µl for the sample and the controls. The master mix volume must be decreased accordingly by reducing the amount of PCR grade Water.

## 4. Starting the reaction

- Load the cycler, check each PCR tube and the cycler lid for tight fit
- Program the qPCR cycler or check stored temperature profiles. See Appendix for temperature profiles of selected qPCR cyclers. Programs for additional cyclers might be available on request.
- Start the program and data reading.

## 5. Result reading

- Save the data at the end of the run.
- Read the channels for the wavelengths 510 and 560/610 nm and show the 2nd deviation of the data.
- Read the calculation of the Ct values for the negative controls, the positive controls and the samples
- Evaluate in accordance with the laboratory and instrument-specific reference ranges

### Reference range for selected qPCR devices:

qPCR device	Manufacturer	Evaluating the threshold cycle		
		positive	borderline	negative
LightCycler® 1.2	Roche Diagnostics	< 40	≥ 40	no Ct
ABI Prism® 7500	Applied Biosystems	< 40	≥ 40	no Ct
RotorGene® 6000	Corbett	< 40	≥ 40	no Ct
Mx3005P®	Agilent Technologies	< 40	≥ 40	no Ct

## NOTES ON THE TEST PROCEDURE

1. This leaflet must be widely understood for a successful use of Venor®GeM-qEP. The reagents supplied should not be mixed with reagents from different lots and used as an integral unit. The reagents of the kit should not be used beyond its shelf life.
2. Any deviation from the test method can affect the results.
3. When using other than the validated MB *Taq* DNA polymerase please note the instructions outlined in the section "Required but not provided material". The amount of polymerase can be increased to 2.5 U/ reaction. Please note that the pipetting scheme of the reaction mix needs to be recalculated.
4. Inhibition may be caused by the sample matrix, but also by sample elution buffer of DNA extraction kits not recommended or validated. Negative controls should always be completed with the use sample elution buffer.
5. For each test setup, at least one negative control may be included reflecting sample preparation. Positive controls facilitate the evaluation of the test. Typical Ct values for the internal control and positive control are shown on the Certificate of Analysis and can be used for orienting quality control.
6. The use of control samples is advised to secure the day-to-day validity of results. The controls should be carried out in the same manner as the samples. It is recommended to run laboratory specific control samples with a high, medial and low (e.g. 3x LOD<sub>95</sub>) level or establish commercial controls, for example, MB Quantification Standards for all mycoplasma species listed in the EP. Participation in external quality control programs, such as offered by Minerva Biolabs on an annual base, is recommended.

## INTERPRETATION OF RESULTS

The presence of *Mycoplasma* in the sample is indicated by an increasing fluorescence signal in the *Mycoplasma* channel during PCR. The concentration of the contaminant can be calculated by a software comparing the crossing cycle number of the sample with a standard curve created in the same run.

A successfully performed PCR without inhibition is indicated by an increasing fluorescence signal in the internal control channel, provided the Internal Control was added to the master mix. *Mycoplasma* DNA and Internal Control DNA are competitors in PCR. Because of the very low concentration of Internal Control in the PCR mix, the signal strength in this channel is reduced with increasing *Mycoplasma* DNA loads in the sample.

Detection of <i>Mollicutes</i> FAM channel	Internal control ROX/VIC channel	Interpretation
positive	irrelevant	<i>Mollicutes</i> positive
negative	negative	PCR inhibition
negative	positive	<i>Mollicutes</i> negative

## ANALYTICAL CHARACTERISTICS OF THE TEST

The positive Cut-Off point was determined for the three Primer/Probe/Nucleotide mixes by using the mycoplasma species listed in the EP 2.6.7 as template either directly from broth culture (in parallel titrated by the traditional culture method), EDQM standard, or by spiking quantified genomic DNA into the sample matrix.

Species	Cut-Off limit		
	DNA extract [GU/ml]	EDQM Standard [CFU/ml]	Culture Broth [CFU/ml]
<i>A. laidlawii</i>	1000	n.a.	1.2
<i>M. hyorhinis</i>	1562	1.2	1.2
<i>M. fermentans</i>	1562	9.6	1.8
<i>M. orale</i>	1562	0.5	1.2
<i>M. synoviae</i>	1562	1.8	1.2
<i>M. pneumoniae</i>	4500	n.a.	9.6
<i>M. arginini</i>	1562	n.a.	1.8

### 2.1 Specificity

Cross-detection of bacteria with close phylogenetic relation to *Mycoplasma* is not monitored. The „European Pharmacopoeia“ recommends to check for unspecific detection of *Clostridium*, *Lactobacillus*, and *Streptococcus*. None of the following species is detected with the Venor®GeM-qEP detection kit: *Clostridium acetobutylicum*, *Lactobacillus acidophilus* and *Streptococcus pneumoniae*. Furthermore there is no positive signal with human, murine or other bacterial DNA (e.g.: *Legionella sp.*, *Chlamydomphila sp.*, *Bordetella sp.*) as template. The following list shows the *Mycoplasma* detectable with Venor®GeM-qEP.

#### Detectable species:

<i>A. laidlawii</i> *	<i>M. cloacale</i>	<i>M. glycyphilum</i>	<i>M. pneumoniae</i> *
<i>M. agalactiae</i>	<i>M. collis</i>	<i>M. gypis</i>	<i>M. pulmonis</i>
<i>M. agassizii</i>	<i>M. columbinasale</i>	<i>M. hominis</i>	<i>M. salivarium</i>
<i>M. alkalescens</i>	<i>M. columbinum</i>	<i>M. hyopharyngis</i>	<i>M. simbae</i>
<i>M. anseris</i>	<i>M. columborale</i>	<i>M. hyorhinis</i>	<i>M. sp.ovine/caprine</i>
<i>M. arginini</i>	<i>M. cricetuli</i>	<i>M. hyosynoviae</i>	<i>M. spermatophilum</i>
<i>M. arthritis</i>	<i>M. cynos</i>	<i>M. iguanae</i>	<i>M. sphenisci</i>
<i>M. bovinegenitalium</i>	<i>M. edwardii</i>	<i>M. indienne</i>	<i>M. spumans</i>
<i>M. bovirhinis</i>	<i>M. equirhinis</i>	<i>M. iners</i>	<i>M. sualvi</i>
<i>M. bovis</i>	<i>M. falconis</i>	<i>M. lagogenitalium</i>	<i>M. subdolum</i>
<i>M. buccale</i>	<i>M. faucium</i>	<i>M. lipofaciens</i>	<i>M. synoviae</i>
<i>M. buteonis</i>	<i>M. felifaucium</i>	<i>M. lipophilum</i>	<i>M. testudineum</i>
<i>M. californicum</i>	<i>M. fermentans</i>	<i>M. maculosum</i>	<i>M. timone</i>
<i>M. canadense</i>	<i>M. gallinaceum</i>	<i>M. meleagridis</i>	<i>M. tumidae</i>
<i>M. capricolum</i>	<i>M. gallinarum</i>	<i>M. moatsii</i>	<i>M. verecundum</i>
<i>M. caviae</i>	<i>M. gallopavonis</i>	<i>M. opalescens</i>	<i>M. zalophi</i>
<i>M. citelli</i>	<i>M. gateae</i>	<i>M. orale</i>	

\*Detection with A.I Mix / Mp Mix

**Appendix: Programming and data recording devices of different qPCR**

**1. LightCycler® 1.2, 1.5 and 2.0**

Program 1: Pre-incubation

Cycles 1  
Analysis Mode None

**Temperature Targets**

Target Temperature [°C] 95  
Incubation time [min] 2:00  
Temperature Transition Rate [°C/s] 20.0  
Secondary Target Temperature [°C] 0  
Step Size [°C] 0.0  
Step Delay [Cycles] 0  
Acquisition Mode None

**LightCycler® 2.0:**

**Please check the correct settings for the „seek temperature“ of at least 90 °C.**

**The pre-incubation time depends on the polymerase used. Hot start enzymes need to be activated at 95°C. Please see polymerase data sheet for duration**

Program 2: Amplification

Cycles 45  
Analysis Mode Quantification

**Temperature Targets**

	Segment 1	Segment 2	Segment 3	Segment 4
Target Temperature [°C]	95	55	60	72
Incubation time [s]	0	5	7	5
Temperature Transition Rate [°C/s]	20.0	20.0	20.0	20.0
Secondary Target Temperature [°C]	0	0	0	0
Step Size [°C]	0.0	0.0	0.0	0.0
Step Delay [Cycles]	0	0	0	0
Acquisition Mode	None	None	Single	None

Program 3: Cooling

Cycles 1  
Analysis Mode None

**Temperature Targets**

Target Temperature [°C] 40  
Incubation time [s] 30  
Temperature Transition Rate [°C/s] 20.0  
Secondary Target Temperature [°C] 0  
Step Size [°C] 0.0  
Step Delay [Cycles] 0  
Acquisition Mode None

Result Reading:

- Select the fluorescence channels F1 and F2 or Channel 1 and 3 at the LC 2.0
- Click on *Quantification* to generate the amplification plots and the specific Ct-values
- The threshold will be generated automatically.
- Samples showing no significant increase in the amplification plot can be considered as negative.

	Mollicutes	Internal Control
channel for LC 1.2, 1.5	F1 (520)	F2 (610)
channel for LC 2.0	channel 1 (520)	channel 3 (610)

**2. RotorGene® 6000 (5-plex)**

Program Step 1: Pre-incubation

Setting Hold  
Hold Temperature 95 °C  
Hold Time 3 min 0 sec

**Please check the correct settings for the filter combination:**

**green filter (470-510): Mollicutes**

**filter orange (585-610): Internal Control**

**The pre-incubation time depends on the polymerase used. Hot start enzymes need to be activated at 95°C. Please see polymerase data sheet for duration**

Program Step 2: Amplification

Setting Cycling  
Cycles 45  
Denaturation 95°C for 5 sec  
Annealing 55°C for 10 sec —> **acquiring to Cycling A (green and orange)**  
Elongation 72°C for 7 sec  
Gain setting automatic (auto Gain)  
Slope Correct activated  
Ignore First deactivated

Result Reading:

- Open the menu *Analysis*
- Select *Quantitation*
- Check the required filter set (green and orange) according to the following table and start data analysis by double click.
- The following windows will appear:  
*Quantitation Analysis - Cycling A* (green or orange)  
*Quant. Results - Cycling A* (green or orange)  
*Standard Curve - Cycling A* (green or orange)
- In window *Quantitation Analysis*, select first *linear scale* and then *slope correct*  
Threshold setup (not applicable if a standard curve was carried with the samples and auto threshold was selected)
  - In window *CT Calculation* set the threshold value to 0-1
  - Pull the threshold line into the graph. Adapt the threshold line to the initial linear section of the positive control reaction.
- The ct-values can be taken from the window *Quant. Results*.
- Samples showing no ct-value can be considered as negative.

Target	Mollicutes	Internal Control
channel	green	orange

### 3. ABI Prism® 7500

#### Detector Settings:

*Bordetella* Target Probe: Reporter - FAM      Quencher - none  
Internal Control Probe: Reporter - ROX      Quencher - none

**The ROX Reference needs to be disabled. Activate both detectors for each well.**

#### Measurement of fluorescence during extension.

##### Program Step 1: Pre-incubation

Setting                      Hold  
Temperature                95 °C  
Incubation time            3:00 min

**The pre-incubation time depends on the polymerase used. Hot start enzymes need to be activated at 95 °C. Please see polymerase data sheet for duration.**

##### Program Step 2: Amplification

Cycles                      45  
Setting                      Cycle  
Denaturing                 95 °C for 30 sek  
Annealing                 55 °C for 30 sek  
Extension                 60 °C for 45 sek

#### Result Reading:

- Enter the following basic setting at the right task bar:  
Data:                      Delta RN vs. Cycle  
Detector:                 FAM and ROX

	<i>Mollicutes</i>	Internal control
channel	FAM	ROX

Line Color:                Well Color

- Open a new window with for the *Graph settings* by clicking the right mouse button  
Select the following setting and confirm with *ok*:  
Real Time Settings:      Linear  
Y-Axis Post Run Settings: Linear and Auto Scale  
X-Axis Post Run Settings: Auto Scale  
Display Options:         2
- Initiate the calculation of the ct-values and the graph generation by clicking on *Analyze* within the report window.
- Pull the threshold line into the graph. Adapt the threshold line to the initial linear section of the positive control reaction.
- Samples showing no ct-value can be considered as negative.

### 4. Mx3005P®

- Go to the setup menu, click on „Plate Setup“, check all positions which apply
- Click on „Collect Fluorescence Data“ and check FAM and ROX
- Corresponding to the basic settings the „Reference Dye“ function should be deactivated
- Specify the type of sample (negative or positive control, sample, standard) at „well type“
- Edit the temperature profile at „Thermal Profile Design“:  
Segment 1 (Pre-Melt): 3 min, 95 °C  
Segment 2: 15 sec, 95 °C  
1 min, 60 °C  
data collection end, 45 cycles
- at menu „Run Status“ select „Run“ and start the cycler by pushing „Start“

#### Analysis of raw data:

- In the window "Analysis" tab on "Analysis Selection / Setup" to analyze the marked positions
- Ensure that in window „algorithm enhancement“ all options are activated:  
Amplification-based threshold  
Adaptive baseline  
Moving average
- Click on „Results“ and „Amplification Plots“ for an automatic threshold
- Read the Ct values at „Text Report“

#### Appendix

##### Limited Product Warranty

This warranty limits our liability for replacement of this product. No warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided. Minerva Biolabs shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

##### Notice to Purchaser

This product is optimized for use in the Polymerase Chain Reaction („PCR“) covered by patents owned by Hoffmann-La Roche, Inc., and F. Hoffmann-La Roche Ltd. („Roche“). No license under these patents to use the PCR process is conveyed expressly or by implication to the purchaser of this product. Minerva Biolabs does not encourage or support the unauthorized or unlicensed use of the PCR process. Use of this product is restricted to persons that either have a license to perform PCR or are not required to obtain a license.

##### Limited License

The use of this product for the detection of Mycoplasma contamination is covered by a Scorpions license from DxS Ltd. Further information on Scorpions licenses can be obtained from DxS Ltd, 48 Grafton St, Manchester, UK M13 9XX.

##### Trademarks

*LightCycler* is a registered trademark of a member of the Roche Group. *ABI Prism* is a registered trademark of Applied Biosystems Corporation or its subsidiaries in the US and certain other countries. *Venor*, *Onar*, *Mynox* and *Mycoplasma Off* are registered trademarks of Minerva Biolabs.

## Related Products

### MB Taq DNA Polymerase

53-0050/-0100/-0200/-0250	MB Taq DNA Polymerase (5 U/ $\mu$ l)	50/100/200/250 units
53-1050/-1100/-1200/-1250	MB Taq DNA Polymerase (1 U/ $\mu$ l)	50/100/200/250 units

### Detection Kits for endpoint PCR

11-1025/-1050/-1100/-1250	Venor <sup>®</sup> GeM Mycoplasma Detection Kit	25/50/100/250 tests
11-7024/-7048/-7096/-7240	Venor <sup>®</sup> GeM Advance Mycoplasma Detection Kit	24/48/96/240 tests
11-2025/-2050/-2100/-2250	Venor <sup>®</sup> GeM OneStep Mycoplasma Detection Kit	25/50/100/250 tests
12-1025/-1050/-1100/-1250	Onar <sup>®</sup> EUB Bacteria Detection Kit	25/50/100/250 tests

### Diagnostic Kits for qPCR

11-4025/100/250	Venor <sup>®</sup> GeM-qEP Mycoplasma Detection Kit	25/100/250 tests
12-2025/100/250	Onar <sup>®</sup> qEUB Bacteria Detection Kit	25/100/250 tests

### Mycoplasma Elimination

10-0200/0500/1000	Mynox <sup>®</sup> Mycoplasma Elimination Reagent	2/5/10 treatments
10-0201/0501/1001	Mynox <sup>®</sup> Gold Mycoplasma Elimination Reagent	2/5/10 treatments

### Quantification Standards, 100 $\mu$ l each, 1x10<sup>6</sup> genomes/ $\mu$ l

52-0116	<i>Acholeplasma laidlawii</i>
52-0129	<i>Mycoplasma arginini</i>
52-0117	<i>Mycoplasma fermentans</i>
52-0130	<i>Mycoplasma hyorhinis</i>
52-0112	<i>Mycoplasma orale</i>
52-0119	<i>Mycoplasma pneumoniae</i>
52-0124	<i>Mycoplasma synoviae</i>
52-0164	<i>Spiroplasma citri</i>

### Genomic DNA Extracts, 100 $\mu$ l each, +/- 10 ng / 100 $\mu$ l (selection)

51-0177	<i>Ureaplasma urealyticum</i> , NC 10177
51-0111	<i>Mycoplasma hominis</i> , NC 010111
51-0112	<i>Mycoplasma orale</i> , NC 010112
51-0117	<i>Mycoplasma fermentans</i> PG19, NC 010117
51-0119	<i>Mycoplasma pneumoniae</i> , NC 010119
51-0129	<i>Mycoplasma arginini</i> , NC 010129
51-0162	<i>Mycoplasma arthritidis</i> , NC 010162
51-0195	<i>Mycoplasma genitalium</i> , NC 010195
51-1746	<i>Mycoplasma penetrans</i> , NC 11746

### DNA Remover™

15-2025	DNA Decontamination Reagent, spray bottle	250 ml
15-2200	DNA Decontamination Reagent, refill bottles	4x 500 ml

### Extraction Kit

56-1100	MB DNA Extraction Kit	100 extractions
---------	-----------------------	-----------------

### Mycoplasma Off®

15-1000	Surface Disinfectant Spray, spray bottle	1000 ml
15-5000	Surface Disinfectant Spray, refill bottles	5x 1000 ml

### Mycoplasma Off® Wipes

15-1001	Surface Disinfectant Wipes	120 wipes
15-5001	Surface Disinfectant Wipes, refill sachets	5 x 120 wipes

### ZellShield®

13-0050/-0150	Cell Culture Contamination Prevention, 100x concentrate	50 ml/ 3x 50 ml
---------------	---	-----------------

### WaterShield™

15-3025/-3075	Water Bath Contamination Prevention, 200x concentrate	5x 5 ml/ 15x 5 ml
---------------	---	-------------------

## Manufacturer

Minerva Biolabs GmbH  
Koeppenicker Str. 325  
D-12555 Berlin  
Germany

## Ordering

Tel. +49 (0)30 2000 437-0  
Fax +49 (0)30 2000 437-9  
order@minerva-biolabs.com

## Product Information

www.minerva-biolabs.com  
info@minerva-biolabs.com

## Technical Service

Tel. +49 (0)30 2000 437-0  
support@minerva-biolabs.com

## Made in Germany

Minerva Biolabs GmbH develops and manufactures products in accordance with EN ISO 9001:2008 and EN ISO 13485:2003 quality system requirement. Reg.No. SY 60026567 0001 & SX 60025009 0001

...as precise as diagnostics should be™