# **Onar®MRSA**

## **Methicillin resistant *Staphylococcus aureus***

### **Diagnostic Kit for qPCR**

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# 1. Reagents and Materials

Please do not use the product after the date of expiry. Discard remaining reagents.

The test provides excellent results with MB *Taq* DNA Polymerase (Cat # 53-0050/0100/0200/0250). We can neither guarantee a high level of sensitivity nor compatibility with other polymerases. If you intend to test our MB *Taq* DNA Polymerase in parallel with your in-house polymerase, please feel free to contact us and get a gratis MB *Taq* sample (10 units).

Onar<sup>®</sup>MRSA is available in tailored package sizes:

Cat.-No. 22-2025	25 tests	1 x 25 reactions
Cat.-No. 22-2100	100 tests	4 x 25 reactions each
Cat.-No. 22-2250	250 tests	10 x 25 reactions each
Cat.-No. 22-2101	100 tests	high throughput (HT) version with 100 reactions per tube

## 1.1 Test Kit Components

*Primer/Probe/Nucleotide Mix (PPN mix)*

Primer set, probes and deoxynucleotide triphosphates dATP, dCTP, dGTP and dUTP, aliquoted for 25 reactions, lyophilized

red cap

*Rehydration Buffer*

1.8 ml


blue cap

*Positive Control DNA*

Fragment of methicillin resistant *Staphylococcus aureus* DNA, prepared by PCR, non-infectious, 300 µl in Tris/HCl


green cap

*Internal Control DNA*

Plasmid DNA, non-infectious, 300 µl in Tris/HCl


yellow cap

*PCR Grade Water*

water for no-template controls, 2 ml

white cap

Please find the guarantee certificate and the conformity declaration on our website.

## 1.2 Stability and Storage

Kit components are stable during shipping. Upon receipt, store at +2°C to +8°C. After rehydration of the *Primer/Probe/Nucleotide Mix* store below -18°C and avoid repeated freezing and thawing. Protect the *Primer/Probe/Nucleotide Mix* from light. For repeated testing of low sample numbers, the *Primer/Probe/Nucleotide Mix* and the controls should be aliquoted. By following these recommendations, the kit is stable until the expiration date stated on the box label.

## 1.3 Supplemental Requirements

- qPCR machine
- corresponding PCR reaction tubes
- microcentrifuge, micropipettes and filtered tips (1 - 1000 µl)
- polymerase

## 2. Application and Test Principle

The Onar<sup>®</sup>MRSA is an in vitro test for semi-quantitative and rapid PCR-based detection of methicillin resistant *Staphylococcus aureus* in clinical samples. This one-step triplex test system is a CE-marked IVD and approved for in vitro clinical diagnostics.

The test can basically be used with all relevant PCR machines with filters for the detection of 520, 560 and 610 nm. The supplied primer, probes and the positive control are specific for a segment of the *mecA* gene of the methicillin resistance encoding gene and the *Staphylococcus aureus* verification gene *femA* as *mecA* is also present in other non-pathogenic *Staphylococcus* species like *S. epidermidis*, *S. intermedius* and *S. pseudointermedius*. The internal control DNA can be added separately directly to the sample to monitor DNA extraction efficiency or later on to the PCR master mix to control the PCR setup individually. The internal control DNA is detected by a homologous primer/probe system.

Onar<sup>®</sup>MRSA contains the nucleotide dUTP instead of dTTP. The heat-labile Uracil-DNA Glycosylase (UNG) is suitable to prevent carry-over contamination between PCRs. This technique relies on the incorporation of deoxyuridine triphosphate (dUTP) during all amplified reactions and the pretreatment of all successive PCR mixtures with the heat-labile UNG. UNG is not provided with Onar<sup>®</sup>MRSA.

## 3. Test Protocol

### 3.1 Preparation of Sample Material

Clinical samples should be from all kinds of infection from the patients. Primarily swabs from skin infections, furuncles, ichors and blood as well as medical equipment such as catheter are suitable immediate testing. Swabs should be sterile and without medium (e.g. Darcan<sup>®</sup>). It is necessary to prepare templates for PCR analysis by DNA extraction using a suitable DNA extraction kit (e.g. MB DNA Extraction Kit, Cat. No. 56-1100). 10 µl of the DNA extract can be used directly as PCR template. The extracts can be stored at temperatures of below -18 °C for a period of one year. Repeated freezing and thawing or storage in the refrigerator for more than 12 hours should be avoided. The sample should not contain more than 100 µg/ml DNA. To avoid false positive results, we recommend the use of aerosol-preventive filter tips and gloves.

### 3.2 Rehydration of the Reagents

1. centrifuge tubes with lyophilized components (5 sec at maximum speed)
2. add 365 µl (2 ml for HT version) of Rehydration Buffer to *Primer/Probe/Nucleotide Mix*
3. incubate for 5 minutes at room temperature
4. vortex and centrifuge again

**Keep reagents on ice and store below -18°C after rehydration.**

### 3.3 Experimental Protocols and Result Reading

#### 3.3.1 LightCycler®2.0

##### Program 1: Pre-incubation

Cycles	1
Analysis Mode	None

##### **Temperature Targets [°C] Segment 1**

Target Temperature [°C]	95
Incubation time [min]	2:00
Temperature Transition Rate [°C/s]	20.0
Secondary Target Temperature [°C]	0
Step Size [°C]	0.0
Step Delay [Cycles]	0
Acquisition Mode	None

##### **LightCycler® 2.0:**

**Please check the correct settings for the „seek temperature“, this should be adjusted on at least 90°C.**

##### Program 2: Amplification

Cycles	45
Analysis Mode	Quantification

##### **Temperature Targets [°C] Segment 1 Segment 2 Segment 3 Segment 4**

Target Temperature [°C]	95	55	60	72
Incubation time [s]	0	5	7	5
Temperature Transition Rate [°C/s]	20.0	20.0	20.0	20.0
Secondary Target Temperature [°C]	0	0	0	0
Step Size [°C]	0.0	0.0	0.0	0.0
Step Delay [Cycles]	0	0	0	0
Acquisition Mode	None	None	Single	None

**Segment 3 is essential for the performance of the kit. Do not skip or modify this step.**

##### Program 3: Cooling

Cycles	1
Analysis Mode	None

##### **Temperature Targets [°C] Segment 1**

Target Temperature [°C]	40
Incubation time [s]	30
Temperature Transition Rate [°C/s]	20.0
Secondary Target Temperature [°C]	0
Step Size [°C]	0.0
Step Delay [Cycles]	0
Acquisition Mode	None

**The incubation time depends on the polymerase used. Hot start enzymes need to be activated at 95°C. Please see polymerase data sheet for duration**

### Result Reading:

- Select the fluorescence channel 1, 2, and 3.
- Click on *Quantification* to generate the amplification plots and the specific *ct*-values
- The threshold will be generated automatically.
- Samples showing no significant increase in the amplification plot can be considered as negative.

target	<i>mecA</i>	<i>femA</i>	Internal control
channel for LC 2.0	channel 1 (530)	channel 3 (610)	channel 2 (560)

### 3.3.2 *Rotorgene 6000 (5-plex)*

#### Program Step 1: Pre-incubation

Setting	Hold
Hold Temperature	95 °C
Hold Time	3 min 0 sec

**Please check the correct settings for the filter combination:**

**green filter (470-510):** MRSA (*mecA*) target

**orange filter (585-610):** *S. aureus* (*femA*)

**yellow filter (530-555):** internal control

#### Program Step 2: Amplification

Setting	Cycling
Cycles	45
Denaturation	95°C for 5 sec
Annealing	55°C for 10 sec —> <b>acquiring to Cycling A (green, orange and yellow)</b>
Elongation	72°C for 7 sec
Gain setting	automatic (auto Gain)
Slope Correct	activated
Ignore First	deactivated

### Result Reading:

- Open the menu *Analysis*
- Select *Quantitation*
- Check the required filter set (green, yellow or orange) according to the following table and start data analysis by double click.
- The following windows will appear:
  - Quantitation Analysis - Cycling A* (green/orange or yellow)
  - Quant. Results - Cycling A* (green/orange or yellow)
  - Standard Curve - Cycling A* (green/orange or yellow)
- In window *Quantitation Analysis*, select first *linear scale* and then *slope correct*
- Threshold setup: (not applicable if a standard curve was carried with the samples and auto threshold was selected)
  - In window *CT Calculation* set the threshold value to 0-1
  - Pull the threshold line into the graph. Adapt the threshold line to the initial linear section of the positive control reaction.
- The *ct*-values can be taken from the window *Quant. Results*.
- Samples showing no *ct*-value can be considered as negative.

target	<i>mecA</i>	<i>femA</i>	Internal control
channel	green	orange	yellow

### 3.3.3 ABI 7500

#### Detector Settings:

MRSA ( <i>mecA</i> ) Target Probe:	Reporter - FAM	Quencher - none
<i>S. aureus</i> ( <i>femA</i> ) Target Probe:	Reporter - ROX	Quencher - none
Internal Control Probe:	Reporter - VIC	Quencher - none

**The ROX Reference needs to be disabled. Activate both detectors for each well.  
Measurement of fluorescence during extension.**

#### Program Step 1: Pre-incubation

Setting	Hold
Temperature	95 °C
Incubation time	3:00 min

**The incubation time depends on the polymerase used. Hot start enzymes need to be activated at 95 °C. Please see polymerase data sheet for duration.**

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#### Result Reading:

- Enter the following basic setting at the right task bar:

Data:	Delta RN vs. Cycle
Detector:	FAM, ROX, VIC etc.
Line Color:	Well Color

target	<i>mecA</i>	<i>femA</i>	Internal control
channel	FAM	ROX	VIC

- Open a new window with for the *Graph settings* by clicking the right mouse button  
Select the following setting and confirm with *ok*:

Real Time Settings:	Linear
Y-Axis Post Run Settings:	Linear and Auto Scale
X-Axis Post Run Settings:	Auto Scale
Display Options:	2
- Initiate the calculation of the ct-values and the graph generation by clicking on *Analyze* within the report window.
- Pull the threshold line into the graph. Adapt the threshold line to the initial linear section of the positive control reaction.
- Samples showing no ct-value can be considered as negative.

### 3.4 PCR Master Mix Setup

Total volume per reaction is 25  $\mu$ l. When setting up reactions, calculations should also include positive (PC) and negative controls (NC). Prepare the master mix in a 1.5 ml reaction tube and mix gently. Do not vortex!

	for 1 reaction	for 25 reactions-
Primer/Probe/Nucleotide Mix	14,0 $\mu$ l	350,0 $\mu$ l
Internal Control DNA	1,0 $\mu$ l	25,0 $\mu$ l
polymerase (5 U/ $\mu$ l)	0,2 $\mu$ l	5,0 $\mu$ l
	15,2 $\mu$ l	380,0 $\mu$ l

**The total duration from master mix preparation to PCR cycling must not exceed 45 minutes to avoid a decrease in the fluorescent signal.**

Aliquot 15  $\mu$ l of master mix into each PCR reaction tube. Pipett 10  $\mu$ l of sample elution buffer or DNA-free water per tube as negative control. Seal the tubes. Add 10  $\mu$ l of sample per PCR reaction tube. Seal these tubes as well before proceeding with the positive controls or standards. Add 10  $\mu$ l of the positive control to the remaining tubes as the last step in order to avoid cross contamination.

master mix	15,0 $\mu$ l
template DNA, NC or PC	10,0 $\mu$ l
total test volume	25,0 $\mu$ l

### 3.5 Result Interpretation

The presence of MRSA in the sample is indicated by an increasing fluorescence signal in the *mecA* channel during PCR. The concentration of the contaminant can be calculated by a software comparing the crossing cycle number of the sample with a standard curve created in the same run. An increase of the fluorescence signal in *femA* channel shows the presence of *S. aureus* and in combination with the amplification of the *mecA* specific gene fragment confirms the presence of MRSA in the sample.

A successfully performed PCR without inhibition is indicated by an increasing fluorescence signal in the internal control channel, provided the Internal Control was added to the master mix. MRSA DNA and Internal Control DNA are competitors in PCR. Because of the very low concentration of Internal Control in the PCR mix, the signal strength in this channel is reduced with increasing MRSA DNA loads in the sample.

<i>mecA</i> PCR	<i>femA</i> PCR	Internal Control	Interpretation
positive	positive	irrelevant	MRSA or MRCoNS* with MSSA
positive	negative	irrelevant	MRCoNS
negative	positive	irrelevant	MSSA
negative	negative	negative	PCR inhibition
negative	negative	positive	MRSA negative

\*MRCoNS: methicillin resistant coagulase negative *Staphylococcus*

## 4. Trouble shooting

Before repeating a negative and a positive control run please check the cycler protocol and pipetting scheme. When using polymerases other than the MB TAQ DNA Polymerase, please note the comments made under chapter 1. The concentration can then be raised up to 2.5 U/reaction. Please note the complete change of the pipetting scheme.

Insufficient (inhibition) or no amplification of the internal control may be due to the following reasons:

- activity of *Taq* polymerase is insufficient or enzyme not compatible with the kit buffer
- programming mistake
- pipetting mistake
- sample elution buffer and sample used to complete the negative control are different

## 5. Product Specification

<i>Intended use</i>	Monitoring of methicillin resistant <i>Staphylococcus aureus</i> (MRSA) DNA in clinical samples
<i>Regulatory status</i>	CE labelled according to IVD Directive 98/79/EC Annex II for the complete process from sample preparation to result
<i>Analytical sensitivity</i>	Detection of less than 10 MRSA particles per test (10 µl).
<i>Diagnostic specificity</i>	98.6 % (1022 swab samples, compared with Chromagar plates)
<i>Diagnostic sensitivity</i>	86.7 % (1022 swab samples , compared with Chromagar plates)
<i>Dynamic range</i>	50 to 5x10 <sup>6</sup> GU/reaction
<i>Subtype inclusivity</i>	All known SCCmec variants
<i>Exclusivity</i>	No cross reactivity with <i>S. epidermidis</i> , <i>S. intermedius</i> and <i>S. pseudointermedius</i> and commensally present organisms in the oral flora or known community acquired pathogens, like <i>Legionella pneumophila</i> or <i>Mycoplasma pneumoniae</i>
<i>Quantitative range</i>	Using MB DNA extraction kit: 150 to 5x10 <sup>6</sup> GU/reaction

## Related Products

### MB Taq DNA Polymerase

53-0050/-0100/-0200/-0250 MB Taq DNA Polymerase 50/100/200/250 units

### Clinical Diagnostic Kits for qPCR

21-2025/-2100/-2250	Onar®Ls <i>Legionella species</i>	25/100/250 tests
21-3025/-3100/-3250	Onar®Lp <i>Legionella pneumophila</i>	25/100/250 tests
20-2025/-2100/-2250	Venor®Mp <i>Mycoplasma pneumoniae</i>	25/100/250 tests
24-2025/-2100/-2250	Onar®Ct <i>Chlamydia trachomatis</i>	25/100/250 tests
23-2025/-2100	Onar®Syphilis <i>Treponema pallidum</i>	25/100 tests
25-2025/-2100/-2250	Onar®Pertussis <i>Bordetella pertussis/parapertussis</i>	25/100/250 tests

### Quantification Standards, 100 µl each, 1x10<sup>8</sup> genomes/µl

52-0101	<i>Legionella pneumophila</i> DNA Standard
52-0119	<i>Mycoplasma pneumoniae</i> DNA Standard

### Genomic DNA Extracts, 100 µl each, +/- 10 ng / 100 µl

51-0566	<i>Streptococcus pneumoniae</i> , DSMZ 20566
51-5571	<i>Bordetella pertussis</i> , DSMZ 5571
51-3415	<i>Bordetella parapertussis</i> , DSMZ 13415
51-0177	<i>Ureaplasma urealyticum</i> , NC 10177
51-0111	<i>Mycoplasma hominis</i> , NC 010111
51-0112	<i>Mycoplasma orale</i> , NC 010112
51-0117	<i>Mycoplasma fermentans</i> PG19, NC 010117
51-0119	<i>Mycoplasma pneumoniae</i> , NC 010119
51-0129	<i>Mycoplasma arginini</i> , NC 010129
51-0162	<i>Mycoplasma arthritidis</i> , NC 010162
51-0195	<i>Mycoplasma genitalium</i> , NC 010195
51-1746	<i>Mycoplasma penetrans</i> , NC 11746

### DNA Remove™

15-2025	DNA Decontamination Reagent, spray bottle	250 ml
15-2200	DNA Decontamination Reagent, refill bottles	4 x 500 ml