Art. Nr. / Order No. 34-2025/-2100/-2250

AquaScreen® *Legionella pneumophila*

PCR Detection Kit for qPCR

version 2.0

Gebrauchsinformation / Instructions for Use

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Reagents for 25/100/250 reactions

Hersteller / Manufacturer:

Minerva Biolabs GmbH, Koepenicker Strasse 325, 12555 Berlin, Germany

FOR IN VITRO USE ONLY

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INDICATION

AquaScreen® *Legionella pneumophila* PCR Detection Kit for the quantitative detection of *Legionella pneumophila* in water samples, prepared with the AquaScreen® FastExtract kit.

EXPLANATION OF THE TEST

Legionella pneumophila is ubiquitous in surface water and moist soil and is known as an opportunistic pathogen of Legionnaires' disease and Pontiac fever. The incubation period is 20-10 days. The airborne transmission is unlikely. Most infections occur through inhalation of aerosols bearing amoeba, the natural habitat of legionella.

AquaScreen® *Legionella pneumophila* is an in vitro assay for the quantitative detection of *Legionella pneumophila* in clinical specimens. Compared to culture methods that can detect legionella first after a few days, the PCR-based detection of this pathogen is completed within 6 hours. PCR also allows for a higher specificity without annoying background growth simplifying result reading drastically. The quantitative result of PCR is linear up to 10⁶ particles/sample with no need for preparing dilutions of the sample material.

The design of the kit complies with the requirements of AFNOR T90-471.

TEST PRINCIPLE

The amplified *mip* gene is highly conserved within the species *Legionella*. It allows detection of all *L. pneumophila* serogroups 1-14 and 16. Cross activities to other water born microorganisms of are not known. The detection limit is 10 genome equivalents per PCR. Because of the wide linear range reliable results without a time-consuming conditioning of the sample material are possible.

The kit contains the nucleotide dUTP instead of dTTP. The heat-labile Uracil-DNA Glycosylase (UNG) is suitable to prevent carry-over contamination between PCRs. This technique relies on the incorporation of deoxyuridine triphosphate (dUTP) during all amplified reactions and the pretreatment of all successive PCR mixtures with the heat-labile UNG. UNG is not included.

REAGENTS

Each kit contains reagents for 25, 100 or 250 reactions. The expiry date of the unopened package is marked on the package label. The kit components are stored until use at +2 to +8 °C. The lot specific Certificate of Analysis can be downloaded from our website (www.minerva-biolabs.com).

Kit Component Label Information	Quantity			Cap Color
	25 reactions Order No. 34-2025	100 reactions Order No. 34-2100	250 reactions Order No. 34-2250	
Lp Mix	1 x lyophilized	4 x lyophilized	10 x lyophilized	red
Rehydration Buffer	1 x 1.8 ml	1 x 1.8 ml	2 x 1.8 ml	blue
Positive Control DNA	1 x lyophilized	1 x lyophilized	1 x lyophilized	green
Internal Control DNA	1 x lyophilized	1 x lyophilized	2 x lyophilized	yellow
PCR grade Water	1 x 2.0 ml	1 x 2.0 ml	1 x 2.0 ml	white

NEEDED, BUT NOT INCLUDED IN THE KIT

The kit contains all reagents for analysis. General industrial supplies and reagents, usually available PCR laboratories are not included:

- qPCR device with filter sets for the detection of the fluorescence dyes FAM™ and ROX™
- PCR reaction tubes for the specific qPCR device
- 1.5 ml reaction tubes, DNA- and RNA-free
- Micro centrifuge for 1.5 ml PCR reaction tubes
- Pipettes with corresponding filter tips to prepare and dispense the reaction mix (10, 100 und 1000 µl)
- Optional: For calibration we recommend our *Legionella pneumophila* DNA Quantification Standard (Cat.-No. 52-0101).

SPECIMEN

The preparation of the sample material is explained in the AquaScreen® FastExtract manual.

The extracts can be stored at a temperature of at least -18 °C for a period of one year. Repeated freezing and defrosting, or storage in the refrigerator for more than 12 hours should be avoided. The sample should not contain more than 100 µg/ml DNA.

PRECAUTIONS

For in vitro diagnostic use. This kit should be used only by trained persons. All samples should be considered potentially infectious and handled at the local or national regulations. This kit does not contain hazardous substances and may be disposed of according to local regulations.

TEST PROCEDURE

The test should be carried out with negative and positive controls and patient samples in duplicate. For quantification, a dilution series of an appropriate standard (Cat. No. 52-0101) is kept. All reagents and samples must be equilibrated to 2 to 8 °C prior use.

1. Rehydration of the Reagents

After reconstitution, the reagents should be stored at 2 to 8 °C for up to 4 days or below -18 °C till the expiration date stated on the box label. Repeated freezing and thawing should be avoided and reconstituted controls (internal amplification control and the positive control) stored in aliquots.

1.	<i>Lp</i> Mix Internal Control Positive Control	red cap yellow cap green cap	Spin all lyophilized components for 5 sec at maximum speed of the mini centrifuge
2.	<i>Lp</i> Mix	red cap	Add 365 µl rehydration buffer (blue cap)
3.	Internal Control	yellow cap	Add 300 µl PCR grade Water (white cap)
4.	Positive Control	green cap	Add 300 µl PCR grade Water (white cap)
5.	<i>Lp</i> Mix Internal Control Positive Control	red cap yellow cap green cap	Incubate 5 min at room temperature
6.	<i>Lp</i> Mix Internal Control Positive Control	red cap yellow cap green cap	Vortex DNA briefly and spin for 5 sec

2. Preparation of the Reaction Mix

Prepare the required amount of master mix at room temperature in a 1.5 ml reaction tube for all control and test reactions.

1.		for 1 reaction	for 25 reactions
	<i>Lp</i> Mix	15.0 µl	350.0 µl
	Internal Control DNA	1.0 µl	25.0 µl
2.	Homogenize the reaction mix by carefully snapping the tube		
3.	Add 15 µl to each PCR tube, discard remaining material.		

3. Loading the test tubes

- Negative Control:
1. Add 10 µl elution buffer from DNA extraction (ref. Chapter sample material) or PCR Grade Water (white cap).
 2. Sample: Add 10 µl of each sample.
 3. Positive Control: Add 10 µl Positive Control (green cap).
 4. Spin all PCR tubes briefly, load the qPCR cyclers and start the program.

Preparation of the master mix and sample loading should not take more than 45 minutes to avoid a reduction in the fluorescent signal. The pipetting sequence should be respected and the tubes closed after each sample load.

4. Starting the reaction

1. Load the cycler, check each PCR tube and the cycler lid for tight fit.
2. Program the qPCR cycler or check stored temperature profiles. See Appendix for temperature profiles of selected qPCR cyclers. Programs for additional cyclers might be available on request.
3. Start the program and data reading.

5. Result reading

1. Save the data at the end of the run.
2. Read the channels for the fluorescence dyes FAM™ and ROX™ and show the 2nd deviation of the data.
3. Read the calculation of the Ct values for the negative controls, the positive controls and the samples.
4. Evaluate in accordance with the laboratory and instrument-specific reference ranges.

Reference range for selected qPCR devices:

qPCR device	Manufacturer	Evaluating the threshold cycle		
		positive	borderline	negative
LightCycler® 1.2	Roche Diagnostics	< 40	≥ 40	no Ct
ABI Prism® 7500	Applied Biosystems	< 40	≥ 40	no Ct
RotorGene® 6000	Corbett Research	< 40	≥ 40	no Ct
Mx3005P®	Agilent Technologies	< 40	≥ 40	no Ct

NOTES ON THE TEST PROCEDURE

1. This leaflet must be widely understood for a successful use of the kit. The reagents supplied should not be mixed with reagents from different lots and used as an integral unit. The reagents of the kit should not be used beyond its shelf life.
2. Any deviation from the test method can affect the results.
3. Inhibition may be caused by the sample matrix, but also by sample elution buffer of DNA extraction kits not recommended or validated. Negative controls should always be completed with the use sample elution buffer.
4. For each test setup, at least one negative control may be included also including sample preparation. Positive controls facilitate the evaluation of the test. Typical Ct values for the internal control and positive control are shown on the Certificate of Analysis and can be used for orienting quality control.
5. The use of control samples is advised to secure the day-to-day validity of results. The controls should be carried out in the same manner as the samples. It is recommended to run laboratory specific control samples with a high, medial and low (e.g. 3x LOD₉₅) level or establish commercial controls, for example, Minerva Biolabs Quantification Standard *Legionella pneumophila* (Order No. 52-0101). Participation in external quality control programs is recommended.

INTERPRETATION OF RESULTS

The presence of *Legionella pneumophila* in the sample is indicated by an increasing fluorescence signal in the *Legionella pneumophila* channel during PCR. The concentration of the contaminant can be calculated by a software comparing the crossing cycle number of the sample with a standard curve created in the same run.

A successfully performed PCR without inhibition is indicated by an increasing fluorescence signal in the internal control channel, provided the Internal Control was added to the master mix. *Legionella pneumophila* DNA and Internal Control DNA are competitors in PCR. Because of the very low concentration of Internal Control in the PCR mix, the signal strength in this channel is reduced with increasing *Legionella pneumophila* DNA loads in the sample.

If a 10 µl test volume is used, then the total quantity of genomic equivalents present in the water sample can be calculated by multiplying the qPCR result by 10 (elution volume of extraction 100 µl).

Example: A 500 ml water sample is used for the analysis. On the basis of the real-time PCR, a result of 60 DNA copies is determined using a test volume of 10 µl. 60 multiplied by 10 totals 600. Result: in 500 ml sample water, there were 600 intact legionella particles (viable culturable, viable but not culturable (VBNC), and dead yet still intact legionella). That results in 120 legionella particles per 100 ml of water.

Detection of <i>Legionella pneumophila</i> FAM™ channel	Internal control ROX™ channel	Interpretation
positive	irrelevant	<i>Legionella pneumophila</i> positive
negative	negative	PCR inhibition
negative	positive	<i>Legionella pneumophila</i> negative


ANALYTICAL CHARACTERISTICS OF THE TEST

The precision of the kit was determined using genomic DNA preparations (*Legionella pneumophila*) at different concentrations generated on a RotorGene® 6000 instrument.

The intra-assay variation shown is based on the scattering of the Ct values of samples of the same concentrations within a test batch, each with four replicates per concentration. Given is the mean coefficient of variation (CV) from five repetitions. The inter-assay variation shown is based on the scattering of the Ct values from samples of equal concentrations of two independent assays, each with four replicates per concentration. Given is the average of the retail price of 3 reps.

Concentration [GU/PCR]	Intraassay Variation VC [%]	Interassay Variation VC [%]	Lot Variation VC [%]
10	2,51	5,27	4,08
100	0,90	1,97	2,97
1000	0,37	1,73	2,80
10000	0,46	1,86	2,93
100000	0,45	2,12	3,50
1000000	0,29	2,14	4,87
10000000	0,33	1,21	2,07

A typical standard curve is shown below. The linear regression was determined from four replicates per dilution level (1 to 10⁶ GU/PCR) with the RotorGene® 6000.

Analytical Sensitivity: 2 GU/μl

Analytical Specificity: 100 % for *Legionella pneumophila*, sero group 1-14 and 16

No cross reactivity with other *Legionella* species, relevant commensals or DNA of human origin

Linearity 50 to 10⁶ GU/PCR

Pipetting time < 60 minutes without influence

Storage of components after rehydration at 2 to 8 °C up to 4 days without effect

Appendix: Programming and data recording devices of different qPCR

1. LightCycler® 1.2, 1.5, 2.0 and 480

Program 1: Pre-incubation

Cycles	1
Analysis Mode	None
Temperature Targets	Segment 1
Target Temperature [°C]	95
Incubation time [min]	2:00
Temperature Transition Rate [°C/s]	20.0
Secondary Target Temperature [°C]	0
Step Size [°C]	0.0
Step Delay [Cycles]	0
Acquisition Mode	None

LightCycler® 2.0:

Please check the correct settings for the „seek temperature“, this should be adjusted on at least 90°C.

Program 2: Amplification

Cycles	45			
Analysis Mode	Quantification			
Temperature Targets	Segment 1	Segment 2	Segment 3	Segment 4
Target Temperature [°C]	95	55	60	72
Incubation time [s]	0	5	7	5
Temperature Transition Rate [°C/s]	20.0	20.0	20.0	20.0
Secondary Target Temperature [°C]	0	0	0	0
Step Size [°C]	0.0	0.0	0.0	0.0
Step Delay [Cycles]	0	0	0	0
Acquisition Mode	None	None	Single	None

Program 3: Cooling

Cycles	1
Analysis Mode	None
Temperature Targets	Segment 1
Target Temperature [°C]	40
Incubation time [s]	30
Temperature Transition Rate [°C/s]	20.0
Secondary Target Temperature [°C]	0
Step Size [°C]	0.0
Step Delay [Cycles]	0
Acquisition Mode	None

Result Reading:

- Select the fluorescence channels F1 and F2 or Channel 1 and 3 at the LC 2.0
- Click on *Quantification* to generate the amplification plots and the specific ct-values
- The threshold will be generated automatically.
- Samples showing no significant increase in the amplification plot can be considered as negative.

target	<i>L. pneumophila</i>	Internal control
channel for LC 1.2, 1.5	F1 (530)	F2 (640)
channel for LC 2.0, 480	channel 1 (530)	channel 3 (610)

2. RotorGene® 6000 (5-plex)

Program Step 1: Pre-incubation

Setting	Hold
Hold Temperature	95 °C
Hold Time	3 min 0 sec

Program Step 2: Amplification

Setting	Cycling
Cycles	45
Denaturation	95°C for 5 sec
Annealing	55°C for 10 sec → acquiring to Cycling A (green and orange)
Elongation	72°C for 7 sec
Gain setting	automatic (auto Gain)
Slope Correct	activated
Ignore First	deactivated

Please check the correct settings for the filter combination:

Target	<i>L. pneumophila</i>	Internal Control
Channel	Green	Orange
Wave length	470-510 nm	585-610 nm

Result Reading:

- Open the menu *Analysis*
- Select *Quantitation*
- Check the required filter set (green and orange) according to the following table and start data analysis by double click.
- The following windows will appear:
 - Quantitation Analysis - Cycling A (green or orange)*
 - Quant. Results - Cycling A (green or orange)*
 - Standard Curve - Cycling A (green or orange)*
- In window *Quantitation Analysis*, select first *linear scale* and than *slope correct*
Threshold setup (not applicable if a standard curve was carried with the samples and auto threshold was selected)
 - In window *CT Calculation* set the threshold value to 0-1
 - Pull the threshold line into the graph. Adapt the threshold line to the initial linear section of the positive control reaction.
- The ct-values can be taken from the window *Quant. Results*.
- Samples showing no ct-value can be considered as negative.

3. ABI Prism® 7500

Detector Settings:

Legionella Target Probe: Reporter - FAM™, Quencher - none
Internal Control Probe: Reporter - ROX™, Quencher - none

The ROX™ Reference needs to be disabled. Activate both detectors for each well.
Measurement of fluorescence during extension.

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Result Reading:

- Enter the following basic setting at the right task bar:

Data:	Delta RN vs. Cycle	Target	<i>L. pneumophila</i>	Internal Control
Detector:	FAM™ and ROX™	Channel	FAM™	ROX™
Line Color:	Well Color			

- Open a new window with for the *Graph settings* by clicking the right mouse button
Select the following setting and confirm with *ok*:
 - Real Time Settings: Linear
 - Y-Axis Post Run Settings: Linear and Auto Scale
 - X-Axis Post Run Settings: Auto Scale
 - Display Options: 2
- Initiate the calculation of the ct-values and the graph generation by clicking on *Analyze* within the report window.
- Pull the threshold line into the graph. Adapt the threshold line to the initial linear section of the positive control reaction.
- Samples showing no ct-value can be considered as negative.

4. Mx3005P®

- Mark the wells in use in the setup menu, tab “Plate Setup”
- Select FAM™ und ROX™ for “Collect Fluorescence Data”
- Check that the Reference Dye function is set to “none” according to the ground settings of the device
- Specify the type of sample (negative control, positive control, standards, sample) by using the tab “well type”
- Program the temperature profile by using the tab “Thermal Profile Design“:
 - Segment 1 (Pre-Melt): 2 min, 95 °C
 - Segment 2:
 - 30 sec, 95 °C
 - 30 sec, 55 °C - data collection end
 - 45 sec, 60 °C - data collection end
 - 45 cycles
- Start the program in the menu “Run” by clicking on “Start” in the window “Run Status”

Interpretation of the raw data:

- Select the “Analysis” menu, click the tab “Analysis Selection/Setup” and mark all positions for analysis
- Activate the following options in window “algorithm enhancement“:
 - Amplification-based threshold
 - Adaptive baseline
 - Moving average
- Set the threshold to automatic in tab “Results” and “Amplification Plots”
- Read the *Ct values* in “Text Report”

APPENDIX

Limited Product Warranty

This warranty limits our liability for replacement of this product. No warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided. Minerva Biolabs shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

Limited License

The use of this product for clinical diagnostics of pathogens is covered by a Scorpions license from DxS Ltd. Further information on Scorpions licenses can be obtained from DxS Ltd, 48 Grafton St, Manchester, UK M13 9XX.

Trademarks

LightCycler is a registered trademark of a member of the Roche Group. ABI Prism is a registered trademark of Applied Biosystems Corporation or its subsidiaries in the US and certain other countries. FAM™ and ROX™ are trademarks of Applied Biosystems Corporation or its subsidiaries in the US and certain other countries. Mx3005P is a trademark of Agilent Technologies. RotorGene is a registered trademark of Qiagen GmbH. AquaScreen is a registered trademark of Minerva Biolabs GmbH.

Related Products

Quantification Standards, 100 µl each, 1x10⁶ genomes/µl

52-0101 *Legionella pneumophila* DNA Standard

DNA Remover™

15-2025 DNA Decontamination Reagent, spray bottle 250 ml
15-2200 DNA Decontamination Reagent, refill bottles 4 x 500 ml

Aqua Screen Fast Extract Kit

32-1010/1050/1200 Extraction Kit 10/50/250 extractions

Diagnostic Kits for qPCR

33-2025/-2100/-2250 AquaScreen® *Legionella specc.* 25/100/250 reactions
34-6025/-6100/-6250 AquaScreen® *TotalBacteria* 25/100/250 reactions

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Made in Germany

Minerva Biolabs GmbH develops and manufactures products in accordance with EN ISO 9001:2008 and EN ISO 13485:2003 quality system requirement. Reg.No. SY 60026567 0001 & SX 60025009 0001

...as precise as diagnostics should be™