***Legionella species* Diagnostic Kit for qPCR**

Version 2.0

-Instructions for Use of Type 1 and Type 2-

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Contents

1. Reagents and Materials	9
1.1 Test Kit Components	9
1.2 Stability and Storage	9
1.3 Supplemental Requirements	9
2. Application and Test Principle	10
3. Test Protocol	10
3.1 Preparation of Sample Material	10
3.2 Rehydration of the Reagents	11
3.3 Experimental Protocol and Result Reading	11
3.3.1 LightCycler® 1.2, 1.5 and 2.0 (for type 1)	11
3.3.2 Rotorgene 6000 (5-plex) (for type 1)	12
3.3.3 ABI 7500 (for type 2)	13
3.4 PCR Master Mix Setup	14
3.5 Result Interpretation	15
4. Trouble Shooting	15
5. Instrument Compatibility	15

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v = validiert

1. Reagents and Materials

Please do not use the product after the date of expiry. Discard remaining reagents.

1.1 Test Kit Components

Primer/Probe/Nucleotide Mix (PPN Mix), type 1 only

Primer set, probes and deoxynucleotide triphosphates dATP, dCTP, dGTP and dUTP, aliquoted for 25 reactions, lyophilized

red cap

REAC	B
------	---

Primer/Probe/Nucleotide Mix (PPN Mix), type 2 only

Primer set, probes and deoxynucleotide triphosphates dATP, dCTP, dGTP and dUTP, aliquoted for 25 reactions, lyophilized

orange cap

REAC	B
------	---

Rehydration Buffer

1.8 ml

blue cap

REAC	A
------	---

Positive Control DNA

Fragment of *Legionella pneumophila* DNA, non-infectious, lyophilized

green cap

CONTROL	+
---------	---

Internal Control DNA

Plasmid DNA, non-infectious, lyophilized

yellow cap

CONTROL	i
---------	---

PCR Grade Water

2 ml

white cap

REAC	H
------	---

1.2 Stability and Storage

Kit components are stable during shipping. Upon receipt, store at +2 °C to +8 °C. After rehydration of the *Primer/Probe/Nucleotide Mix*, the *Internal Control Probe*, the *Positive Control* and the *Internal Control*, store below -18 °C and avoid repeated freezing and thawing. Protect the *Primer/Probe/Nucleotide Mix* and *Internal Control Probe* from light. For repeated testing of low sample numbers, *Primer/Probe/Nucleotide Mix*, *Internal Control Probe* and controls should be aliquoted after rehydration. By following these recommendations, the kit is stable until the expiration date stated on the Guarantee Certificate (please find the Certificate on our web-site).

1.3 Supplemental Requirements

- real-time PCR thermal cycler (check list on page 15 for compatibility)
- vials required for the particular thermal cycler used
- micro centrifuge, micropipettes and filtered tips
- polymerase

The test provides excellent results with MB Taq DNA-Polymerase (Cat # 53-0050; 53-0100; 53-0200; 53-0250). We can neither guarantee a high level of sensitivity nor compatibility with other polymerases. If you intend to test our MB Taq DNA-Polymerase in parallel with your in-house polymerase, please feel free to contact us and get a gratis MB Taq sample (10 units). If you want to use your own polymerase, it may be necessary to use the specific buffer provided with this polymerase.

2. Application and Test Principle

Legionella spp. ubiquitously colonize surface water and wet soil and is the cause of the legionnaires' disease. Transmission is most likely aerogenic. Typical symptoms are high fever, cough, thorax ache, diarrhea, discomposure, sometimes with heavy progression.

The Onar[®]Ls test system is an in vitro test for the qualitative and the quantitative diagnosis of *Legionella* spp. in clinical samples. The test is based on the polymerase chain reaction. The supplied primer set is specific for a segment of the 16S rRNA region of the *Legionella* genome. The selected template is highly preserved within the *Legionella* species. Clinically relevant species being detected are e.g. *L. pneumophila* SG 1-14 and 16, *L. longbeachae*, *L. dumoffii*, *L. bozemanii* and *L. gormanii*. Cross activity with pharyngeal commensals or other pathogens in the respiratory tract, like *Bordetella*, *Chlamydia*, *Mycoplasma*, *Mycobacteria*, *Pseudomonas* and *Streptococci* is not known. The detection limit is 10 genome equivalents per 10 µl sample volume.

The target probe emits fluorescent light at 520 nm. The kit also provides internal control DNA, which can be added to the reaction. When running the PCR with the internal control DNA, a successful performed reaction with a negative result is indicated by a distinct fluorescent signal. The internal control is detected by another probe, which emits light at ~610 nm (type 1) and ~560 nm (type 2).

The Onar[®]Ls PCR kit contains the nucleotide dUTP instead of dTTP. The heat-labile Uracil-DNA Glycosylase (UNG) is suitable to prevent carry-over contamination between PCRs. This technique relies on the incorporation of deoxyuridine triphosphate (dUTP) during all amplified reactions and the pretreatment of all successive PCR mixtures with the heat-labile UNG. UNG is not provided with Onar[®]Ls.

3. Test Protocol

3.1 Preparation of Sample Material

Nasopharyngeal swabs, nasal and pharyngeal secretions, sputum, provoked sputum, bronchial lavage, tissue and infected cells or cultures can be used as test materials. The quality of sample-taking has a major influence on the reliability of the test results. Material from the lower respiratory tract is particularly suitable for detecting legionella. A fine catheter and a vacuum pump with secretion trap can be used for nasopharyngeal secretions. The oral cavity can be rinsed beforehand, and saliva and rinsing water discarded. Taking a sample from the oral cavity should not be preceded by eating or gargling. The swab is turned upwards behind the arch of the palate and wiped. When taking a sample through the nasal cavity, the thin swab is advanced into the nasopharynx and wiped repeatedly. Ensure that adequate material is obtained. The swab should be washed out in 1 ml transfer buffer (e.g. 50 mM Tris-HCl, pH 7.5, 0,1% sodium azide, 0,5% sodium dodecylsulfate). The sample material is stable for a few days and should be transported cooled.

DNA extraction with a suitable DNA extraction kit (e.g. Minerva Biolabs, MB DNA Extraction Kit) is always advisable when preparing the samples in order to remove inhibitors of the PCR safely and to concentrate the legionella DNA at greater sample volumes. The obtained DNA extract can be used directly for the Onar[®]Ls test. The extracts can be stored at a temperature of at least -18 °C for a period of one year. Repeated freezing and defrosting, or storage in the refrigerator for more than 12 hours should be avoided. The sample should not contain more than 100 µg/ml DNA.

3.2 Rehydration of the Reagents

1. centrifuge tubes with lyophilized components (5 sec at maximum speed)
2. add Rehydration Buffer

<i>Primer/Probe/Nucleotide Mix</i>	365 μ l
------------------------------------	-------------
4. add appropriate amount of deionized, DNA-free water

<i>Internal Control DNA</i>	300 μ l
<i>Positive Control DNA</i>	300 μ l
3. incubate for 5 minutes at room temperature
4. vortex and centrifuge again
5. Store reagents below -18 °C after rehydration.

3.3 Experimental Protocols and Result Reading

Protocols for other qPCR machines are available from our web page.

This detection kit is designed as a dual color detection experiment. Due to the spectral overlap of the emission spectra of the two dyes, with some qPCR machines (e.g. LightCycler from Roche Diagnostics) a residual cross talk of the target channel into the internal control channel can be observed. To correct for this cross talk, color compensation files have to be generated in a separate run using a Color Compensation Kit offered by the company providing the real-time PCR instrument. In this case the test results are still interpretable as an amplification of the target DNA took necessarily place.

3.3.1 LightCycler® 1.2, 1.5 and 2.0 (type 1)

Program 1: Pre-incubation

Cycles	1
Analysis Mode	None
Temperature Targets [°C]	Segment 1
Target Temperature [°C]	95
Incubation time [min]	2:00
Temperature Transition Rate [°C/s]	20.0
Secondary Target Temperature [°C]	0
Step Size [°C]	0.0
Step Delay [Cycles]	0
Acquisition Mode	None

LightCycler® 2.0:

Please check the correct settings for the „seek temperature“, this should be adjusted on at least 90 °C.

Program 2: Amplification

Cycles 45
Analysis Mode Quantification

Temperature Targets [°C]	Segment 1	Segment 2	Segment 3	Segment 4
Target Temperature [°C]	95	55	60	72
Incubation time [s]	0	5	7	5
Temperature Transition Rate [°C/s]	20.0	20.0	20.0	20.0
Secondary Target Temperature [°C]	0	0	0	0
Step Size [°C]	0.0	0.0	0.0	0.0
Step Delay [Cycles]	0	0	0	0
Acquisition Mode	None	None	Single	None

Segment 3 is essential for the performance of the kit. Do not skip or modify this step.

Program 3: Cooling

Cycles 1
Analysis Mode None

Temperature Targets [°C]	Segment 1
Target Temperature [°C]	40
Incubation time [s]	30
Temperature Transition Rate [°C/s]	20.0
Secondary Target Temperature [°C]	0
Step Size [°C]	0.0
Step Delay [Cycles]	0
Acquisition Mode	None

Result Reading:

- Select the fluorescence channels F1 and F2 or Channel 1 and Channel 3 at the LC 2.0
- Click on *Quantification* to generate the amplification plots and the specific ct-values
- The threshold will be generated automatically.
- Samples showing no significant increase in the amplification plot can be considered as negative.

target	L. species	Internal control
channel for LC 1.2, 1.5	F1 (530)	F2 (640)
channel for LC 2.0	channel 1 (530)	channel 3 (610)

3.3.2 Rotorgene 6000 (5-plex) (type 1)

Program Step 1: Pre-incubation

Setting Hold
Hold Temperature 95 °C
Hold Time 3 min 0 sec

Please check the correct settings for the filter combination:

green filter (470-510): L. species
filter orange (585-610): internal control

Program Step 2: Amplification

Setting	Cycling
Cycles	45
Denaturation	95°C for 5 sec
Annealing	55°C for 10 sec —> acquiring to Cycling A (green and orange)
Elongation	72°C for 7 sec
Gain setting	automatic (auto gain)
Slope Correct	activated
Ignore First	deactivated

Result Reading:

- Open the menu *Analysis*
- Select *Quantitation*
- Check the required filter set (green or orange) according to the following table and start data analysis by double click.
- The following windows will appear:
 - Quantitation Analysis - Cycling A* (green or orange)
 - Quant. Results - Cycling A* (green or orange)
 - Standard Curve - Cycling A* (green or orange)
- In window *Quantitation Analysis*, select first *linear scale* and than *slope correct*.
 - Threshold setup (not applicable if a standard curve was carried with the samples and auto threshold was selected)
 - In window *CT Calculation* set the threshold value to 0-1
 - Pull the threshold line into the graph. Adapt the threshold line to the initial linear section of the positive control reaction.
- The ct-values can be taken from the window *Quant. Results*.
- Samples showing no ct-value can be considered as negative.

target	L. species	Internal control
channel	green	orange

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Please check the correct settings of the detectors!

target probe: FAM
internal control: VIC
quencher: none

The ROX Reference needs to be disabled. Activate both detectors for each well. Measurement of fluorescence during extension.

Result Reading:

- Enter the following basic setting at the right task bar:

Data: Delta RN vs. Cycle
Detector: FAM and VIC
Line Color: Well Color

target	L. species	Internal control
channel	FAM	VIC

- Open a new window with for the *Graph settings* by clicking the right mouse button
Select the following setting and confirm with *ok*:
Real Time Settings: Linear
Y-Axis Post Run Settings: Linear and Auto Scale
X-Axis Post Run Settings: Auto Scale
Display Options: 2
- Initiate the calculation of the ct-values and the graph generation by clicking on *Analyze* within the report window.
- Pull the threshold line into the graph. Adapt the threshold line to the initial linear section of the positive control reaction.
- Samples showing no ct-value can be considered as negative.

3.4 PCR Master Mix Setup

Total volume per reaction is 25 μl . When setting up reactions, calculations should also include positive (PC) and negative controls (NC). Pipette master mix on ice into a 1.5 ml reaction tube and mix gently.

Pipetting scheme:

	for 1 reaction	for 25 reactions*
<i>Primer/Probe/Nucleotide Mix</i>	14.0 μl	350.0 μl
<i>Internal Control DNA</i>	1.0 μl	25.0 μl
polymerase (5 U/ μl)	0.2 μl	5.0 μl

+ template DNA, NC or PC 10.0 μl

* the equivalent of the content of one red-capped vial

The total duration from master mix preparation to PCR cycling must not exceed 45 minutes to avoid a decrease in the fluorescent signal.

Aliquot 15 μl of master mix into each PCR reaction tube. The excess of master mix (1 μl) is rejected. Pipette 10 μl of sample elution buffer or DNA-free water per tube as negative control. Seal the tubes. Add 10 μl of sample per PCR reaction tube. Seal these tubes as well before proceeding with the positive controls or standards. Add 10 μl of the positive control to the remaining tubes as the last step in order to avoid cross contamination.

3.5 Result Interpretation

The presence of *Legionella* DNA in the sample is indicated by an increasing fluorescence signal in the target channel during PCR. The concentration of the contaminant can be calculated by a software comparing the crossing cycle number of the sample with a standard curve created in the same run.

A successfully performed PCR without inhibition is indicated by an increasing fluorescence signal in the internal control channel, provided the Internal Control was added to the master mix. *Legionella* DNA and Internal Control DNA are competitors in PCR. Because of the very low concentration of Internal Control in the PCR mix, the signal strength in this channel is reduced with increasing *Legionella* DNA loads in the sample.

L. species	Internal Control	Interpretation
positive	irrelevant	<i>Legionella species</i> positive
negative	negative	PCR inhibition
negative	positive	<i>Legionella species</i> negative

4. Trouble shooting

Before repeating a negative and a positive control run please check the cycler protocol and pipetting scheme. Insufficient (inhibition) or no amplification of the internal control may be due to the following reasons:

- programming mistake
- pipetting mistake
- sample elution buffer and sample used to complete the negative control are different

5. Instrument Compatibility

Instrument	type 1	type 2
LightCycler® 1.2	+v	–
LightCycler® 1.5	+	–
LightCycler® 2.0	+	o
LightCycler® 480	+	o
Rotorgene 3000	+	o
Rotorgene 6000 (5-plex)	+v	o
ABI Prism® 7000	–	+
ABI Prism® 7300	–	+
ABI Prism® 7500	–	+v
ABI Prism® 7700	–	+
ABI Prism® 7900	–	+
iCycler iQ®	+v	o
iQ™5	o	o
Opticon 2	–	o
Chromo 4	o	o
MX3000P	o	o
MX4000	o	o

+ = recommended Kit version

– = the internal control is not detectable

o = untested but presumed to be compatible

v = validated

Appendix

Limited Product Warranty

This warranty limits our liability for replacement of this product. No warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided. Minerva Biolabs shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

§Notice to Purchaser

This product is optimized for use in the Polymerase Chain Reaction („PCR”) covered by patents owned by Hoffmann-La Roche, Inc., and F. Hoffmann-La Roche Ltd. („Roche”). No license under these patents to use the PCR process is conveyed expressly or by implication to the purchaser of this product. Minerva Biolabs does not encourage or support the unauthorized or unlicensed use of the PCR process. Use of this product is restricted to persons that either have a license to perform PCR or are not required to obtain a license.

Limited License

The use of this product for the detection of mycoplasma is covered by a Scorpions license from DxS Ltd. Further information on Scorpions licenses can be obtained from DxS Ltd, 48 Grafton St, Manchester, UK M13 9XX.

Trademarks

LightCycler is a registered trademark of a member of the Roche Group. ABI Prism is a registered trademark of Applied Biosystems Corporation or its subsidiaries in the US and certain other countries. Venor and Onar are registered trademarks of Minerva Biolabs GmbH.

Related Products

MB Taq DNA Polymerase

53-0050/-0100/-0200/-0250 MB Taq DNA Polymerase 50/100/200/250 units

Clinical Diagnostic Kits for qPCR

20-2025/-2100/-2250	Venor [®] Mp-QP <i>Mycoplasma pneumoniae</i>	25/100/250 tests
21-2025/-2100/-2250	Onar [®] Ls-QP <i>Legionella species</i>	25/100/250 tests
21-3025/-3100/-3250	Onar [®] Lp-QP <i>Legionella pneumophila</i>	25/100/250 tests
22-2025/-2100/-2250	Onar [®] MRSA Methicillin-resistant <i>S. aureus</i>	25/100/250 tests
23-2025/-2100	Onar [®] Tp <i>Treponema pallidum</i>	25/100 tests
24-2025/-2100/-2250	Onar [®] Ct <i>Chlamydia trachomatis</i>	25/100/250 tests
25-2025/-2100/-2250	Onar [®] Pertussis <i>Bordetella pertussis/parapertussis</i>	25/100/250 tests

Quantification Standards, 100 µl each, 1x10⁶ genomes/µl

52-0101	<i>Legionella pneumophila</i> DNA Standard
52-0119	<i>Mycoplasma pneumoniae</i> DNA Standard

Genomic DNA Extracts, 100 µl each, +/- 10 ng / 100 µl

51-0566	<i>Streptococcus pneumoniae</i> , DSMZ 20566
51-5571	<i>Bordetella pertussis</i> , DSMZ 5571
51-3415	<i>Bordetella parapertussis</i> , DSMZ 13415
51-0177	<i>Ureaplasma urealyticum</i> , NC 10177
51-0111	<i>Mycoplasma hominis</i> , NC 010111
51-0112	<i>Mycoplasma orale</i> , NC 010112
51-0117	<i>Mycoplasma fermentans</i> PG19, NC 010117
51-0119	<i>Mycoplasma pneumoniae</i> , NC 010119
51-0129	<i>Mycoplasma arginini</i> , NC 010129
51-0162	<i>Mycoplasma arthritidis</i> , NC 010162
51-0195	<i>Mycoplasma genitalium</i> , NC 010195
51-1746	<i>Mycoplasma penetrans</i> , NC 11746

DNA Remover™

15-2025	DNA Decontamination Reagent, spray bottle	250 ml
15-2200	DNA Decontamination Reagent, refill bottles	4 x 500 ml

Extraktion Kit

56-1100	MB DNA Extraction Kit	100 extractions
---------	-----------------------	-----------------