## Experimental Protocol and Evaluation for Cepheid Smart Cycler® with Minerva Biolabs real-time PCR kits

### 1. Experimental Protocol

#### Phase 1:

Stopp

temp	sec	optic
95	180	off

#### Phase 2:

Repeat 40


4-temperaturecycle

temp	sec	optic
95	30	off
55	30	off
60	30	on
72	30	off

#### Phase 3:

Stopp

temp	sec	optic
40	30	off

The screenshot shows the 'Smart Cycler Dx' software interface. The main window displays 'Analyse: Experimentelle Untersuchung' with the following table:

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Please check the correct settings for the filter combination of Cepheid Smart Cycler®. For the target you must choose the filter FAM (483-533) and for the Internal Control-Target the filter TxR (558-610).

## 2. Evaluation and interpretation of the results

The analysis of the obtained data is divided into two parts:

- Quantitative analysis of DNA in fluorescence channel FAM-Signals (483-533)
- Quantitative analysis of Internal Control DNA in fluorescence channel TxR-Signals (558-610).

The amplification curves were obtained by performing the described procedure with a dilution series of a Quantification Standard and the Cepheid Smart Cycler<sup>®</sup> instrument. The fluorescence values versus cycle number are displayed. In the same run the amplification of the Internal Control DNA was shown in channel TxR (558-610).