

ExtractNow™ RNA Mini Kit

Protocol 1: RNA extraction from tissue/biopsy samples



reconstitute with absolute ethanol



disrupt tissue



+ ≤ 20 mg tissue
+ 450 μl DL



Proceed with the next step or store the sample in Lysis Buffer D at -20 °C.

⊗ max. speed / 1 min proceed with complete supernatant

Protocol 2: RNA extraction from eukaryotic cells



reconstitute with absolute ethanol ≤ 5 x 10⁶ cells



+ 400 μl DL
2 min
⊗ resuspend 20 sec
⊗ 3 min



Note: ensure that all cells are completely lysed (no visible cell clumps) before you proceed to the next step.

Protocol 3: RNA extraction from bacteria



reconstitute with absolute ethanol ≤ 5 x 10⁹ bacteria

prepare stock solution of lysozyme

pellet cells for 2 - 5 min at 5000 × g
discard supernatant

+ resuspend in 100 μl TE
+ 2 μl (20 mg/ml) or
6 μl (50 mg/ml) lysozym stock solution




mix
⊗ incubate until the solution becomes clear or viscous

+ 450 μl DL
⊗ 3 min



+ supernatant
⊗ 10,000 × g for 2 min

Important:
The flow-through contains the RNA.
Thus do not discard the flow-through!



+ equal volume 70 % ethanol to the flow-through

mix

+ add
⊗ vortex
⊗ incubate
⊗ centrifuge

⊗ 6000 × g for 1 min

store the extracted RNA at -20 °C or at -80 °C for long time storage.

⊗ 10,000 × g for 3 min
change tube

+ 30 to 80 μl RNase-free water
⊗ 1 min

+ 500 μl CW
⊗ 10,000 × g for 1 min
change tube

+ 700 μl DW
⊗ 10,000 × g for 1 min
change tube



+ load flow-through/ ethanol mix
⊗ 10,000 × g for 1 min
change tube