**Art. Nr. / Order No. 11-8025/-8050/-8100/-8250**

# Quantification Standard

*Pseudomonas aeruginosa Titrated Genomic DNA*

## **Gebrauchsinformation / Instructions for Use**

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### INDICATION

Genomic DNA of *Pseudomonas aeruginosa* can be used as amplification and sensitivity control of end point (gel-based evaluation) or quantitative, real-time PCR.

### EXPLANATION OF THE TEST

This product provides isolated genomic DNA from *Pseudomonas aeruginosa* (DSMZ 50071). To obtain the microorganisms nutrient broth was inoculated with *Pseudomonas aeruginosa*. The microorganisms were harvested in the late log phase of growth by repeated washing and centrifugation. The DNA was extracted with a classic phenol-chloroform protocol and further purified with adsorption chromatography. The concentration of the DNA was quantified photometrically (calibrated to weight standards) and by qPCR (compared against exactly quantified calibrator plasmids). The DNA concentration was adjusted with regular TE80 buffer.

### TEST PRINCIPLE

Pure and titrated genomic DNA can be used as amplification and sensitivity control of end point PCR (gel-based evaluation). For quantitative PCR the Quantification Standard can serve for creating standard curves by using dilutions of the material as sample in the PCR. The software of various devices will be able to calculate from qPCR data corresponding concentrations and will generate a standard curve, which can be used to determine unknown DNA concentrations.

### REAGENTS

Each kit contains reagents for 10 dilution series. The expiry date of the unopened package is marked on the package label. The kit components are stored below  $-18\text{ }^{\circ}\text{C}$  until use. Repeated freezing and thawing should be avoided. The opened kit is stable for 6 months at  $+2$  to  $+8\text{ }^{\circ}\text{C}$  if used under sterile conditions. The lot specific quality control certificate (*Certificate of Analysis*) can be downloaded from our website ([www.minerva-biolabs.com](http://www.minerva-biolabs.com)).

Kit Component	Quantity	Cap Color
Quantification Standard $1 \times 10^6$ genomes/ $\mu\text{l}$	100 $\mu\text{l}$	green
Tris Buffer 10 mM Tris-HCl, pH 8.5	3 x 1 ml	white

## NEEDED, BUT NOT INCLUDED IN THE KIT

The kit contains the reagents for the preparation of the dilution series. General industrial supplies and reagents, usually available in molecular biology laboratories are not included:

- 1.5 ml reaction tubes, DNA- and RNA-free
- Micro centrifuge for 1.5 ml reaction tubes
- Pipettes with corresponding filter tips to prepare and dispense the mix (10, 100 und 1000  $\mu$ l)
- vortex

## PRECAUTIONS

The kit is intended for research use only. Not for clinical diagnostics or testing of human samples without extensive validation.

This kit should be used only by trained persons.


This kit does not contain hazardous or infectious substances and may be disposed of according to local regulations.

Cross contamination may lead to false-positive results. The test should be performed according to good laboratory practice.

## PROCEDURE

### 1. Preparation of the Dilution Series

1. Thaw *Quantification Standard* and provided *Tris Buffer*.
2. Label six 1.5 ml reaction tubes consecutively and fill each with 90  $\mu$ l of *Tris Buffer*.
3. Vortex *Quantification Standard* briefly (1 to 2 seconds) at medium speed.
4. Add 10  $\mu$ l of the *Quantification Standard* to reaction tube no. 1, close the tube and vortex briefly at medium speed.
5. Add 10  $\mu$ l of the content of reaction tube no. 1 to reaction tube no. 2, close the tube and vortex briefly at medium speed.
6. Add 10  $\mu$ l of the content of reaction tube no. 2 to reaction tube no. 3, close the tube and vortex briefly at medium speed.
7. Proceed with the following reaction tubes of the dilution series in the same way.

### 2. PCR

Performing DNA amplification with the use of validated PCR assays is highly recommended (eg. Minerva Biolabs' Onar<sup>®</sup>qEUB Bacteria Detection Kit for qPCR, order no. 12-2xxx). Please follow the PCR kit manual while treating the dilutions as samples. For qPCR these samples have to be edit in the qPCR cyclor program as standards with their dedicated concentrations.

The amount of solution used as sample material defines the number of genome copies per PCR reaction:

Reaction tube no.	2 $\mu$ l sample volume	5 $\mu$ l sample volume	10 $\mu$ l sample volume
<b>1</b>	2x10 <sup>5</sup> genome copies	5x10 <sup>5</sup> genome copies	1x10 <sup>6</sup> genome copies
<b>2</b>	2x10 <sup>4</sup> genome copies	5x10 <sup>4</sup> genome copies	1x10 <sup>5</sup> genome copies
<b>3</b>	2x10 <sup>3</sup> genome copies	5x10 <sup>3</sup> genome copies	1x10 <sup>4</sup> genome copies
<b>4</b>	200 genome copies	500 genome copies	1000 genome copies
<b>5</b>	20 genome copies	50 genome copies	100 genome copies
<b>6</b>	2 genome copies	5 genome copies	10 genome copies

### 3. Evaluation

In qPCR the Ct-values should decrease linearly with ascending DNA concentration per reaction, when a suitable PCR assay is used. The software of the qPCR device calculates a standard curve and slope using the DNA concentrations stated by the user and the appendant Ct-values. Also the Ct-values of samples with unknown DNA concentrations are automatically compared to the standard curve and concentrations are assigned.

The following figures were generated using Minerva Biolabs Onar<sup>®</sup>qEUB Bacteria Detection Kit for qPCR and the RotorGene<sup>®</sup> 6000 instrument.


Fig. 1: Amplified dilution series of  $10^4$  down to 10 genome copies of *Pseudomonas aeruginosa* DNA.

Fig. 2: Standard curve generated with the RotorGene<sup>®</sup> 6000 instrument using second derivative maximum method and the data from Fig. 1

### Limited Product Warranty

This warranty limits our liability for replacement of this product. No warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided. Minerva Biolabs shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

### Trademarks

Onar is a registered trademark of Minerva Biolabs GmbH, Germany.

RotorGene is a registered trademark of Corbett Research, Australia.

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## Related Products

### MB Taq DNA Polymerase

53-0050/-0100/-0200/-0250	MB Taq DNA Polymerase (5 U/ $\mu$ l)	50/100/200/250 units
53-1050/-1100/-1200/-1250	MB Taq DNA Polymerase (1 U/ $\mu$ l)	50/100/200/250 units

### Diagnostic Kits for conventional PCR

11-1025/-1050/-1100/-1250	Venor <sup>®</sup> GeM Mycoplasma Detection Kit	25/50/100/250 tests
11-2025/-1050/-1100/-1250	Venor <sup>®</sup> GeM OneStep Mycoplasma Detection Kit	25/50/100/250 tests
11-7024/-7048/-7096/-7240	Venor <sup>®</sup> GeM Advance Mycoplasma Detection Kit	24/48/96/240 tests
12-1025/-1050/-1100/-1250	Onar <sup>®</sup> EUB Bacteria Detection Kit	25/50/100/250 tests

### Diagnostic Kits for qPCR

11-4025/100/250	Venor <sup>®</sup> GeM-qEP Mycoplasma Detection Kit	25/100/250 tests
12-1025/-1100/-1250	Onar <sup>®</sup> qEUB Bacteria Detection Kit	25/100/250 tests

### Mycoplasma Elimination

10-0200/0500/1000	Mynox <sup>®</sup> Mycoplasma Elimination Reagent	2/5/10 treatments
10-0201/0501/1001	Mynox <sup>®</sup> Gold Mycoplasma Elimination Reagent	2/5/10 treatments

### Quantification Standards, 100 $\mu$ l each, 1x10<sup>8</sup> genomes/ $\mu$ l

52-0116	<i>Acholeplasma laidlawii</i>
52-0129	<i>Mycoplasma arginini</i>
52-0117	<i>Mycoplasma fermentans</i>
52-0130	<i>Mycoplasma hyorhinis</i>
52-0112	<i>Mycoplasma orale</i>
52-0119	<i>Mycoplasma pneumoniae</i>
52-0124	<i>Mycoplasma synoviae</i>
52-0164	<i>Spiroplasma citri</i>
52-0115	<i>Mycoplasma gallisepticum</i>

### Genomic DNA Extracts, 100 $\mu$ l each, +/- 10 ng / 100 $\mu$ l (selection)

51-0177	<i>Ureaplasma urealyticum</i> , NC 10177
51-0111	<i>Mycoplasma hominis</i> , NC 010111
51-0112	<i>Mycoplasma orale</i> , NC 010112
51-0117	<i>Mycoplasma fermentans</i> PG19, NC 010117
51-0119	<i>Mycoplasma pneumoniae</i> , NC 010119
51-0129	<i>Mycoplasma arginini</i> , NC 010129
51-0162	<i>Mycoplasma arthritidis</i> , NC 010162
51-0195	<i>Mycoplasma genitalium</i> , NC 010195
51-1746	<i>Mycoplasma penetrans</i> , NC 11746

### DNA Remover<sup>™</sup>

15-2025	DNA Decontamination Reagent, spray bottle	250 ml
15-2200	DNA Decontamination Reagent, refill bottles	4 x 500 ml

### Extraktion Kit

56-1100	MB DNA Extraction Kit	100 extractions
---------	-----------------------	-----------------

### Mycoplasma Off<sup>®</sup>

15-1000	Surface Disinfectant Spray, spray bottle	1000 ml
15-5000	Surface Disinfectant Spray, refill bottles	5 x 1000 ml

### Mycoplasma Off<sup>®</sup>-Wipes

15-1001	Surface disinfectant Wipes in dispenser box	120 wipes
15-5001	Surface Disinfectant Wipes, refill pack	5 x 120 wipes

### ZellShield<sup>™</sup>

13-0050/-0150	Contamination Prevention Reagent 100x concentrate	1000 ml/ 5 x 1000 ml
---------------	---	----------------------

### WaterShield<sup>™</sup>

13-0050/-0150	Contamination Prevention Reagent 100x concentrate	1000 ml/ 5 x 1000 ml
---------------	---	----------------------

### PCR Ultra Pure Water

55-0015/-0075	PCR Grade Water for all PCRs and RT-PCRs, cDNA synthesis sequencing and cycle sequencing	1.5 ml/ 5 x 1.5 ml
---------------	--	--------------------

## Manufacturer

Minerva Biolabs GmbH  
Koeppenicker Str. 325  
D-12555 Berlin  
Germany

## Ordering

Tel. +49 (0)30 2000 437-0  
Fax +49 (0)30 2000 437-9  
order@minerva-biolabs.com

## Product Information

www.minerva-biolabs.com  
info@minerva-biolabs.com

## Technical Service

Tel. +49 (0)30 2000 437-0  
support@minerva-biolabs.com

## Made in Germany

Minerva Biolabs GmbH develops and manufactures products in accordance with EN ISO 9001:2008 and EN ISO 13485:2003 quality system requirement. Reg.No. SY 60026567 0001 & SX 60025009 0001

...as precise as diagnostics should be<sup>™</sup>