Quantification Standard

Mycoplasma fermentans

genomic DNA isolate, 10⁶ genomes/ μ l, 100 μ l

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1.2 Stability and Storage

The *Quantification Standard* is stable during shipping on cool packs and should be stored below –18 °C. The *Quantification Standard* should be aliquoted in order to avoid repeated freezing and thawing. By following these recommendations, the product is stable until the expiration date stated on the box label. A *Certificate of Analysis* is available from our web page.

1.3 Supplemental Requirements

microcentrifuge
micro pipettes and filter tips
vortex

2. Application

This product provides isolated genomic DNA from *Mycoplasma fermentans*. To obtain the microorganisms a special liquid culture medium was inoculated with an early passage strain of *Mycoplasma fermentans* (NCTC 10117) and harvested at the end of the logarithmic growth phase by repeated washing and centrifugation. The DNA was extracted with a classic phenol-chloroform protocol and further purified with adsorption chromatography. The concentration of the DNA was quantified photometrically (calibrated to weight standards) and by qPCR (compared against exactly quantified calibrator plasmids). The DNA concentration was adjusted with regular TE80 buffer.


Genomic DNA of *Mycoplasma fermentans* can be used as amplification and sensitivity control of end point PCR (gel-based evaluation). For quantitative PCR the *Quantification Standard* can serve for creating standard curves by using dilutions of the material as sample in the PCR. The software of various devices will be able to calculate from qPCR data corresponding concentrations and will generate a standard curve, which can be used to determine unknown DNA concentrations.

3. Protocol

3.1 Preparation of the Dilution Series

1. Thaw *Quantification Standard* and provided *Tris Buffer*.
2. Label six 1.5 ml reaction tubes consecutively and fill each with 90 μl of *Tris Buffer*.
3. Vortex *Quantification Standard* briefly (1 to 2 seconds) at medium speed.
4. Add 10 μl of the *Quantification Standard* to reaction tube no. 1, close the tube and vortex briefly at medium speed.

5. Add 10 μl of the content of reaction tube no. 1 to reaction tube no. 2.
6. Close the tube and vortex briefly at medium speed.
7. Add 10 μl of the content of reaction tube no. 2 to reaction tube no. 3.
8. Close the tube and vortex briefly at medium speed.
9. Proceed with the following reaction tubes of the dilution series in the same way.

3.2 PCR

Performing DNA amplification with the use of validated PCR assays is highly recommended (eg. Minerva Biolabs' Venor[®]GeM-qEP *Mycoplasma* Detection Kit for q PCR, order no. 11-4025). Please follow the PCR kit manual.

The amount of solution used as sample material defines the number of genome copies per PCR reaction:

| Reaction tube no. | 2 μl sample volume | 5 μl sample volume | 10 μl sample volume |
|-------------------|-------------------------------|-------------------------------|--------------------------------|
| 1 | 2×10^5 genome copies | 5×10^5 genome copies | 1×10^6 genome copies |
| 2 | 2×10^4 genome copies | 5×10^4 genome copies | 1×10^5 genome copies |
| 3 | 2×10^3 genome copies | 5×10^3 genome copies | 1×10^4 genome copies |
| 4 | 200 genome copies | 500 genome copies | 1000 genome copies |
| 5 | 20 genome copies | 50 genome copies | 100 genome copies |
| 6 | 2 genome copies | 5 genome copies | 10 genome copies |

3.3 Evaluation

In qPCR the Ct-values should decrease linearly with ascending DNA concentration per reaction, when a suitable PCR assay is used. The software of the qPCR device calculates a standard curve and slope using the DNA concentrations stated by the user and the appendant Ct-values. Also the Ct-values of samples with unknown DNA concentrations are automatically compared to the standard curve and concentrations are assigned.

The following figures were generated using Minerva Biolabs Venor®GeM-qEP Mycoplasma Detection Kit for qPCR and the LightCycler® instrument.


Fig. 1: Amplified dilution series of 10^4 down to 10^0 genome copies of *Mycoplasma fermentans* as starting template. The fluorescent channel was set to F1.


Fig. 2: Standard curve generated with the LightCycler® instrument using second derivative maximum method and the data from Fig. 1

Limited Product Warranty

This warranty limits our liability to replacement of this product. No warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided. Minerva Biolabs shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

Notice to Purchaser

This product is optimized for use in the Polymerase Chain Reaction („PCR”) covered by patents owned by Hoffmann-La Roche, Inc., and F. Hoffmann-La Roche Ltd. („Roche”). No license under these patents to use the PCR process is conveyed expressly or by implication to the purchaser of this product. Minerva Biolabs does not encourage or support the unauthorized or unlicensed use of the PCR process. Use of this product is restricted to persons that either have a license to perform PCR or are not required to obtain a license.

LightCycler is a registered trademark of a member of the Roche Group.

Venor is a registered trademark of Minerva Biolabs.

Related Products

Polymerases

| | | | |
|------------------------|--------------------------|----------------|-----------|
| 53-0050/0100/0200/0250 | MB TAQ DNA Polymerase | 50/100/200/250 | units |
| 54-0100/0500 | EUB Polymerase, DNA-free | | 100 units |

Diagnostic Kits for conventional PCR

| | | | |
|---------------------|---|---------------|-------|
| 11-1025/050/100/250 | Venor [®] GeM Mycoplasma Detection Kit | 25/50/100/250 | tests |
| 11-7025/050/100/250 | Venor [®] GeM Advance Mycoplasma Detection Kit | 24/48/96/240 | tests |
| 12-1025/050/100/250 | Onar [®] EUB Eubacteria Detection Kit | 25/50/100/250 | tests |

Diagnostic Kits for qPCR

| | | | |
|-----------------|--|------------|-------|
| 11-4025/100/250 | Venor [®] GeM-qEP Mycoplasma Detection Kit Type 1 | 25/100/250 | tests |
| 11-5025/100/250 | Venor [®] GeM-qEP Mycoplasma Detection Kit Type 2 | 25/100/250 | tests |

Mycoplasma Elimination

| | | | |
|-------------------|--|--------|------------|
| 10-0200/0500/1000 | Mynox [®] Mycoplasma Elimination Reagent | 2/5/10 | treatments |
| 10-0201/0501/1001 | Mynox [®] Gold Mycoplasma Elimination Reagent | 2/5/10 | treatments |

Genomic DNA Extracts, +/- 10 ng/100 µl, 100 µl each

| | |
|---------|--|
| 51-1723 | <i>Lactobacillus acidophilus</i> , NC 001723 |
| 51-0566 | <i>Streptococcus pneumoniae</i> , DSMZ 20566 |
| 51-0792 | <i>Clostridium acetobutylicum</i> , DSMZ 792 |
| 51-5571 | <i>Bordetella pertussis</i> , DSMZ 5571 |
| 51-3415 | <i>Bordetella parapertussis</i> , DSMZ 13415 |
| 51-0177 | <i>Ureaplasma urealyticum</i> , NC 10177 |
| 51-0111 | <i>Mycoplasma hominis</i> , NC 010111 |
| 51-0112 | <i>Mycoplasma orale</i> , NC 010112 |
| 51-0113 | <i>Mycoplasma salivarium</i> , NC 010113 |
| 51-0115 | <i>Mycoplasma gallisepticum</i> , NC 010115 |
| 51-0116 | <i>Acholeplasma laidlawii</i> , NC 010116 |
| 51-0117 | <i>Mycoplasma fermentans</i> PG19, NC 010117 |
| 51-0119 | <i>Mycoplasma pneumoniae</i> , NC 010119 |
| 51-0124 | <i>Mycoplasma synoviae</i> , NC 0100124 |
| 51-0129 | <i>Mycoplasma arginini</i> , NC 010129 |
| 51-0130 | <i>Mycoplasma hyorhinis</i> , NC 010130 |
| 51-0162 | <i>Mycoplasma arthritidis</i> , NC 010162 |
| 51-0195 | <i>Mycoplasma genitalium</i> , NC 010195 |

Quantification Standard, 1x10⁶ genomes/µl, 100 µl each

| | |
|---------|-------------------------------|
| 52-0112 | <i>Mycoplasma orale</i> |
| 52-0116 | <i>Acholeplasma laidlawii</i> |
| 52-0119 | <i>Mycoplasma pneumoniae</i> |
| 52-0124 | <i>Mycoplasma synoviae</i> |

Mycoplasma Off[®]

| | | | |
|---------|------------------------------------|---------|----|
| 15-1000 | Surface Disinfectant Spray | 1000 | ml |
| 15-5000 | Surface Disinfectant Spray, refill | 5x 1000 | ml |

ZellShield™

| | | | |
|---------|---|----|----|
| 13-0050 | Microbial Contamination Preventive Reagent for Cell Cultures, 100x solution | 50 | ml |
|---------|---|----|----|

DNA Remover™

| | | | |
|---------|---|--------|----|
| 15-2025 | DNA Decontamination Reagent, spray bottle | 250 | ml |
| 15-2200 | DNA Decontamination Reagent, refill bottles | 4x 500 | ml |