Art. Nr. / Order No. 52-0440

Quantification Standard

Chlamydia trachomatis Titrated Genomic DNA

Gebrauchsinformation / Instructions for Use

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INDICATION

Genomic DNA of *Chlamydia trachomatis* can be used as amplification and sensitivity control of end point (gel-based evaluation) or quantitative, real-time PCR.

EXPLANATION OF THE TEST

This product provides isolated genomic DNA from *Chlamydia trachomatis*. To obtain the microorganisms BGM cells were inoculated with *Chlamydia trachomatis* (serovar E, strain DK20). The inclusion bodies were harvested by repeated washing and centrifugation. The DNA was extracted with a classic phenol-chloroform protocol and further purified with adsorption chromatography. The concentration of the DNA was quantified photometrically (calibrated to weight standards) and by qPCR (compared against exactly quantified calibrator plasmids). The DNA concentration was adjusted with regular TE80 buffer.

TEST PRINCIPLE

Pure and titrated genomic DNA can be used as amplification and sensitivity control of end point PCR (gel-based evaluation). For quantitative PCR the Quantification Standard can serve for creating standard curves by using dilutions of the material as sample in the PCR. The software of various devices will be able to calculate from qPCR data corresponding concentrations and will generate a standard curve, which can be used to determine unknown DNA concentrations.

REAGENTS

Each kit contains reagents for 10 dilution series. The expiry date of the unopened package is marked on the package label. The kit components are stored below -18 °C until use. Repeated freezing and thawing should be avoided. The opened kit is stable for 6 months at +2 to +8 °C if used under sterile conditions. The lot specific quality control certificate (*Certificate of Analysis*) can be downloaded from our website (www.minervabiolabs.com).

Kit Component	Quantity	Cap Color
Quantification Standard 1x10 ⁶ genomes/μl	100 μl	green
<i>Tris Buffer</i> 10 mM Tris-HCl, pH 8.5	3 x 1 ml	white

NEEDED, BUT NOT INCLUDED IN THE KIT

The kit contains the reagents for the preparation of the dilution series. General industrial supplies and reagents, usually available in molecular biology laboratories are not included:

- 1.5 ml reaction tubes, DNA- and RNA-free
- Micro centrifuge for 1.5 ml reaction tubes
- Pipettes with corresponding filter tips to prepare and dispense the mix (10, 100 und 1000 μ l)
- vortex

PRECAUTIONS

The kit is intended for research use only. Not for clinical diagnostics or testing of human samples without extensive validation.

This kit should be used only by trained persons.


This kit does not contain hazardous or infectious substances and may be disposed of according to local regulations.

Cross contamination may lead to false-positive results. The test should be performed according to good laboratory practice.

PROCEDURE

1. Preparation of the Dilution Series

1. Thaw *Quantification Standard* and provided *Tris Buffer*.
2. Label six 1.5 ml reaction tubes consecutively and fill each with 90 μl of *Tris Buffer*.
3. Vortex *Quantification Standard* briefly (1 to 2 seconds) at medium speed.
4. Add 10 μl of the *Quantification Standard* to reaction tube no. 1, close the tube and vortex briefly at medium speed.
5. Add 10 μl of the content of reaction tube no. 1 to reaction tube no. 2, close the tube and vortex briefly at medium speed.
6. Add 10 μl of the content of reaction tube no. 2 to reaction tube no. 3, close the tube and vortex briefly at medium speed.
7. Proceed with the following reaction tubes of the dilution series in the same way.

2. PCR

Performing DNA amplification with the use of validated PCR assays is highly recommended (eg. Minerva Bi-olabs' Onar[®]Ct *Chlamydia trachomatis* Diagnostic Kit for qPCR, order no. 24-2xxx). Please follow the PCR kit manual while treating the dilutions as samples. For qPCR these samples have to be edit in the qPCR cyclor program as standards with their dedicated concentrations.

The amount of solution used as sample material defines the number of genome copies per PCR reaction:

Reaction tube no.	2 μl sample volume	5 μl sample volume	10 μl sample volume
1	2×10^5 genome copies	5×10^5 genome copies	1×10^6 genome copies
2	2×10^4 genome copies	5×10^4 genome copies	1×10^5 genome copies
3	2×10^3 genome copies	5×10^3 genome copies	1×10^4 genome copies
4	200 genome copies	500 genome copies	1000 genome copies
5	20 genome copies	50 genome copies	100 genome copies
6	2 genome copies	5 genome copies	10 genome copies

3. Evaluation

In qPCR the Ct-values should decrease linearly with ascending DNA concentration per reaction, when a suitable PCR assay is used. The software of the qPCR device calculates a standard curve and slope using the DNA concentrations stated by the user and the appendant Ct-values. Also the Ct-values of samples with unknown DNA concentrations are automatically compared to the standard curve and concentrations are assigned.

The following figures were generated using Minerva Biolabs Onar[®]Ct *Chlamydia trachomatis* Diagnostic Kit for qPCR and the LightCycler[®] 1.5 (Roche Diagnostics GmbH, Germany) instrument.


Fig. 1: Amplified dilution series of 10^4 down to 10 genome copies of *Chlamydia trachomatis* as starting template.


Fig. 2: Standard curve generated with the LightCycler 1.5 instrument using second derivative maximum method and the data from Fig. 1

Limited Product Warranty

This warranty limits our liability for replacement of this product. No warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided. Minerva Biolabs shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

Trademarks

Venor is a registered trademark of Minerva Biolabs GmbH, Germany.

LightCycler is a registered trademark of Roche Diagnostics GmbH, Germany.

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Related Products

MB Taq DNA Polymerase

53-0050/-0100/-0200/-0250	MB Taq DNA Polymerase (5 U/ μ l)	50/100/200/250 units
53-1050/-1100/-1200/-1250	MB Taq DNA Polymerase (1 U/ μ l)	50/100/200/250 units

Diagnostic Kits for conventional PCR

11-1025/-1050/-1100/-1250	Venor [®] GeM Mycoplasma Detection Kit	25/50/100/250 tests
11-7024/-7048/-7096/-7240	Venor [®] GeM Advance Mycoplasma Detection Kit	24/48/96/240 tests
12-1025/-1050/-1100/-1250	Onar [®] EUB Eubacteria Detection Kit	25/50/100/250 tests

Diagnostic Kits for qPCR

11-4025/100/250	Venor [®] GeM-qEP Mycoplasma ssp. Detection Kit	25/100/250 tests
-----------------	--	------------------

Mycoplasma Elimination

10-0200/0500/1000	Mynox [®] Mycoplasma Elimination Reagent	2/5/10 treatments
10-0201/0501/1001	Mynox [®] Gold Mycoplasma Elimination Reagent	2/5/10 treatments

Quantification Standards, 100 μ l each, 1x10⁶ genomes/ μ l

52-0116	<i>Acholeplasma laidlawii</i>
52-0129	<i>Mycoplasma arginini</i>
52-0117	<i>Mycoplasma fermentans</i>
52-0130	<i>Mycoplasma hyorhinis</i>
52-0112	<i>Mycoplasma orale</i>
52-0119	<i>Mycoplasma pneumoniae</i>
52-0124	<i>Mycoplasma synoviae</i>
52-0164	<i>Spiroplasma citri</i>
52-0115	<i>Mycoplasma gallisepticum</i>

Genomic DNA Extracts, 100 μ l each, +/- 10 ng / 100 μ l (selection)

51-0177	<i>Ureaplasma urealyticum</i> , NC 10177
51-0111	<i>Mycoplasma hominis</i> , NC 010111
51-0112	<i>Mycoplasma orale</i> , NC 010112
51-0117	<i>Mycoplasma fermentans</i> PG19, NC 010117
51-0119	<i>Mycoplasma pneumoniae</i> , NC 010119
51-0129	<i>Mycoplasma arginini</i> , NC 010129
51-0162	<i>Mycoplasma arthritidis</i> , NC 010162
51-0195	<i>Mycoplasma genitalium</i> , NC 010195
51-1746	<i>Mycoplasma penetrans</i> , NC 11746

DNA Remover™

15-2025	DNA Decontamination Reagent, spray bottle	250 ml
15-2200	DNA Decontamination Reagent, refill bottles	4 x 500 ml

Extraktion Kit

56-1100	MB DNA Extraction Kit	100 extractions
---------	-----------------------	-----------------

Mycoplasma Off[®]

15-1000	Surface Disinfectant Spray, spray bottle	1000 ml
15-5000	Surface Disinfectant Spray, refill bottles	5 x 1000 ml

Mycoplasma Off[®]-Wipes

15-1001	Surface disinfectant Wipes in dispenser box	120 wipes
15-5001	Surface Disinfectant Wipes, refill pack	5 x 120 wipes

ZellShield™

13-0050/-0150	Contamination Prevention Reagent 100x concentrate	1000 ml/ 5 x 1000 ml
---------------	---	----------------------

PCR Ultra Pure Water

55-0015/-0075	PCR Grade Water for all PCRs and RT-PCRs, cDNA synthesis sequencing and cycle sequencing	1.5 ml/ 5 x 1.5 ml
---------------	--	--------------------

Manufacturer

Minerva Biolabs GmbH
Koenpinner Str. 325
D-12555 Berlin
Germany

Ordering

Tel. +49 (0)30 2000 437-0
Fax +49 (0)30 2000 437-9
order@minerva-biolabs.com

Product Information

www.minerva-biolabs.com
info@minerva-biolabs.com

Technical Service

Tel. +49 (0)30 2000 437-0
support@minerva-biolabs.com

Made in Germany

Minerva Biolabs GmbH develops and manufactures products in accordance with EN ISO 9001:2008 and EN ISO 13485:2003 quality system requirement.
Reg.No. SY 60026567 0001 & SX 60025009 0001

...as precise as diagnostics should be™