Art. Nr. / Order No. 11-7024/-7048/-7096/-7240

Venor[®] GeM Advance

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1. Reagents and Materials

1.1 Kit Components

Test Reaction Tubes tubes pre-coated with lyophilized primers, dNTPs, internal control DNA, and gel loading buffer/dye	3, 6, 12, or 30 strips of 8 transparent tubes each, depending on package size
Positive Control Reaction Tubes Test Reaction Tubes containing DNA-fragments of <i>Mycoplasma orale</i> genome	1, 2, 3, or 5 strips of 8 red tubes each, depending on package size
Caps for PCR Tubes	4, 6, 15, or 35 strips of 8 tubes each, depending on package size
Rehydration Buffer	1.6 ml

1.2 Stability and Storage

Kit components are stable during shipping. Upon receipt, store at +2 to +8 °C. The kit is stable until the expiration date stated on the Guarantee Certificate (online: <http://www.minerva-biolabs.com/>) or the product box.

1.3 Supplemental Requirements

PCR thermal cycler
agarose gel electrophoresis apparatus
micro centrifuge, micropipettes and filtered tips
polymerase

The test provides excellent results with MB *Taq* DNA Polymerase (Cat # 53-0050/0100/0200/0250). We can neither guarantee a high level of sensitivity nor compatibility with other polymerases. If you intend to test our MB *Taq* DNA Polymerase in parallel with your in-house polymerase, please feel free to contact us and get a gratis MB *Taq* sample (10 units). However, if you want to use your own polymerase, it may be necessary to use the specific buffer provided with this polymerase.

2. Application and Test Principle

Venor®GeM Advance utilizes the polymerase chain reaction (PCR), which was established as the method of choice for highest sensitivity in the detection of *Mycoplasma* and *Acholeplasma* contamination in cell cultures and other cell culture derived biologicals. Detection requires as little as 1 to 5 fg of mycoplasma DNA corresponding to 2-5 mycoplasma per sample volume. The primer set is specific to the highly conserved rRNA operon, or more specifically, the 16S rRNA coding region in the mycoplasma genome. This allows for detection of *M. orale*, *M. hyorhinis*, *M. arginini*, *M. fermentans*, *M. salivarium*, *M. hominis*, usually encountered as contaminants in cell cultures, but also *M. pneumoniae*, *Acholeplasma laidlawii*, *M. synoviae* and *Ureaplasma* species. Cross-detection of bacteria with close phylogenetic relation to *Mycoplasma* is not monitored. The „European Pharmacopoeia“ recommends to check for unspecific detection of *Clostridium*, *Lactobacillus*, and *Streptococcus*. None of the following species is detected with Venor®GeM Advance: *Clostridium acetobutylicum*, *Lactobacillus acidophilus* and *Streptococcus pneumoniae*. Furthermore there is no positive signal with human DNA as template.

Because the reaction tubes included with the kit are pre-coated with appropriate dNTPs and primers, the total assay time is greatly reduced compared to general protocols that require individual loading of reaction tubes. After thermal cycling the PCR can be loaded directly on the agarose gel.

Venor®GeM Advance also contains internal control DNA in the PCR mix. When running the PCR a successfully performed reaction is indicated by a 191 bp band on the agarose gel.

Venor®GeM Advance is intended for research use only. Not for clinical diagnostics or testing of human samples.

3. Test Protocol

3.1 Preparation of Sample Material

Samples should be derived from cultures which are at 90-100 % confluence. PCR inhibiting substances may accumulate in the medium of older cultures. For these sample materials a DNA extraction (e.g. MB DNA Extraction Kit Cat # 56-1100) is strictly recommended prior testing. Penicillin and streptomycin in the culture media do not inhibit mycoplasma or affect test sensitivity. Only cell culture supernatant should be applied to test for mycoplasma. Cell pellets should only be tested after suitable DNA extraction, since debris will interfere with the PCR reaction. With average titers at 10^6 and a maximum titer at 10^8 you will find sufficient mycoplasma in the supernatant to guarantee a sensitive PCR. However, other materials that can be tested are Fetal Calf Serum, vaccines, and paraffin-embedded samples following DNA extraction. If necessary, templates for PCR analysis are prepared by DNA extraction using commercially available extraction kits (MB DNA Extraction Kit Cat. No. 56-1100). Please be sure to remove any alcohol containing wash buffer from the preparation to avoid coelution of alcohol and sample material. Any remaining alcohol may inhibit the PCR. 2 µl of the extract can be used directly as PCR template. To avoid false positive results, we recommend the use of the PCR grade water delivered with the kit, aerosol-preventive filter tips and gloves.

The preparation of sample material could be performed by one of the following methods:

Heat-inactivation of the sample material

The templates for the PCR analysis are prepared by direct heating of the cell culture supernatant or the biological sample material:

- 100 µl liquid supernatant of the sample material is transferred into a sterile reaction tube;
- the supernatant is incubated at 95 °C for 10 minutes;
- the supernatant is centrifuged briefly (5 seconds, 1000 x g) to remove cellular debris.
- The supernatant is used in the PCR. Alternatively, the DNA can be purified with a commercial extraction kit (e.g. MB DNA Extraction Kit Cat # 56-1100).

Enrichment of mycoplasma by centrifugation

- 1 ml liquid supernatant of the sample material is transferred into a sterile reaction tube;
- the supernatant is centrifuged (15 minutes, 10.000 x g) to sediment mycoplasma particles. Alternatively: centrifuge the supernatant 6 min with 13.000 x g.
- The supernatant is rejected and the pellet is suspended into 50 µl buffer (10 mM Tris, pH 8.4).
- The sample should be vortexed and finally heated up to 95 °C for 10 min.

The extracts can be stored at a temperature of at least -18 °C for a period of one year. Repeated freezing and defrosting, or storage in the refrigerator for longer than 12 hours should be avoided. The sample should not contain more than 100 µg/ml DNA.

3.2 Thermal Profile

The programming process of your cycler is explained in the manual of the instrument.

Program

1 cycle	94 °C for 2 min
39 cycles	94 °C for 30 sec
	55 °C for 30 sec
	72 °C for 30 sec

cool down to 4 °C to 8 °C

The incubation time depends on the polymerase used. Some hot start enzymes need to be activated at 94 °C for more than 2 minutes. Please see polymerase data sheet for duration.

3.3 PCR Setup

1. Mix Rehydration Buffer and Polymerase

For optimal comparison of all samples tested the polymerase should not be added separately to the PCR tubes but premixed with the Rehydration Buffer.

Determine the total number of reactions including negative and positive controls. Calculate the required amount of Rehydration Buffer and Polymerase (1 unit/reaction) including an additional reaction (23 µl) to compensate for pipetting losses:

	Test Reaction	Negative Control	Positive Control
Rehydration Buffer	22.8 μ l	22.8 μ l	24.8 μ l
Polymerase (5 U/ μ l)	0.2 μ l	0.2 μ l	0.2 μ l

For other polymerase concentrations the amount of water needs to be adjusted.

Vortex the Rehydration Buffer vigorously (at least 30 sec), transfer the calculated amount into a fresh minicentrifuge tube and add the required amount of polymerase. Mix the Rehydration Buffer/Polymerase Mix by flicking the tube.

Do not vortex diluted polymerase.

2. Rehydration of PCR reagents

Remove and cut off the required amount of Test Reaction Tubes (transparent) and Positive Control Tubes (red) from the bag. Replace remaining tubes in bag and seal properly. Label tubes as appropriate. Centrifuge or tip the tubes on the table to collect the lyophilized material at the bottom of the tube. Peel off protective film from tubes. Aliquot 23 μ l of the prepared Rehydration Buffer/Polymerase Mix into each PCR reaction tube. Rehydrate Positive Control Tubes with 25 μ l of the mix.

Avoid cross contamination by preparing the Positive Controls after finishing the Negative Controls and Test Reactions.

The appearance of the reagent pellets can vary on technical reasons (yellowish pinhead-sized to small, dark red spot), but has no effect on the quality of the test.

3. Sample addition

The total volume per reaction is 25 μ l. Add the corresponding sample material directly to the PCR tubes:

Negative Control: 2 μ l of deionized water or elution buffer if extraction kit used

Test Reaction: 2 μ l of sample

Positive Control: no sample required

Close tubes with Caps for PCR Tubes included in kit. Mix contents of each tube thoroughly by flicking tubes. Do not vortex! Collect liquid contents at the bottom of the PCR tube by brief centrifugation. Incubate at room temperature for 5 minutes. Proceed immediately to thermal cycling.

3.4 Agarose Gel Run

Prepare a 1.5% standard agarose gel, approx. 5 mm thick, with a 5 mm-comb.

Load 5 μ l of each PCR reaction. No loading buffer and running dye are required.

Stop electrophoresis after 2 cm run distance (depending on the electrophoresis chamber used e.g. run for 20 minutes at 100 V).

3.5 Gel Evaluation

A distinct 191 bp band should appear in every lane indicating a successfully performed PCR. This band may fade out with increased amounts of amplicons formed, caused by mycoplasma DNA loads of $> 5 \times 10^6$ copies/ml. The initial concentration of positive control DNA exceeds 5×10^6 copies/ml in order to account for DNA loss resulting from repeated freezing and thawing.

Relevant amplicon sizes:

Internal control	191 bp
<i>Mycoplasma spp.</i> (see table in the appendix)	265-278 bp

Controls:

negative control band at 191 bp

positive control band at 267 bp, possibly an additional band at 191 bp


Interpretation of Test Reaction band patterns:

Band pattern	Interpretation
band at 191 bp	negative sample
band at 270 bp and at 191 bp	mycoplasma-positive sample with weak contamination
strong band at 270 bp	mycoplasma-positive sample, highly contaminated
no band	PCR inhibition or activity of polymerase is insufficient

With Venor[®]GeM designed for high sensitivity and therefore prone to nonspecific annealing, bands of various length that are less intensive can be produced, but do not indicate positive results. Possible primer self-annealing produces another band of 80-90 bp in length, but also does not affect the precision or results of the test.

If the PCR of a sample is inhibited, PCR inhibitors can easily be removed from the sample by performing a DNA extraction with a commercially available kit. A list of recommended DNA extraction kits is provided in the appendix.

Figure 1: Amplicon pattern

4. Trouble Shooting

No amplification of control DNA may be due to the following reasons:

- activity of polymerase is insufficient
- control DNA tubes have not been spun down before rehydration
- programming mistake
- pipetting mistake

Before rerun of a negative and a positive control please check thermocycler protocol and pipetting scheme. When using polymerases other than the MB TAQ DNA Polymerase, please note the comments made under chapter 1.3. The enzyme concentration can then be raised up to 2.5 U/reaction. Please note the complete change of the pipetting scheme.

Appendix

Detection Range and Sizes of Amplicons

No.	species	amplicon size (bp)
1	<i>Mycoplasma orale</i> ^{1,2}	266
2	<i>Mycoplasma pneumoniae</i> ²	273
3	<i>Mycoplasma penetrans</i>	274
4	<i>Mycoplasma pirum</i>	274
5	<i>Acholeplasma laidlawii</i> ²	273
6	<i>Mycoplasma fermentans</i>	267
7	<i>Ureaplasma urealyticum</i>	273
8	<i>Mycoplasma hyorhini</i> ²	268
9	<i>Mycoplasma pulmonis</i>	268
10	<i>Mycoplasma falconis</i>	268
11	<i>Mycoplasma arthritidis</i>	267
12	<i>Mycoplasma arginini</i>	267
13	<i>Mycoplasma spermatophilum</i>	267
14	<i>Mycoplasma opalescens</i>	266
15	<i>Mycoplasma primatum</i>	267
16	<i>Mycoplasma maculosum</i>	267
17	<i>Mycoplasma bovis</i>	267
18	<i>Mycoplasma cloacale</i>	266
19	<i>Mycoplasma hyosynoviae</i>	265
20	<i>Mycoplasma synoviae</i> ²	266
21	<i>Mycoplasma salivarium</i>	266
22	<i>Mycoplasma faucium</i>	265
23	<i>Mycoplasma hominis</i>	266
24	<i>Mycoplasma genitalium</i>	273
25	<i>Mycoplasma bovigenitalium</i>	267
26	<i>Mycoplasma sp. ovine/caprine</i>	267
27	<i>Mycoplasma agalactica</i>	267
28	<i>Mycoplasma timone</i>	266

¹ provided as positive control DNA

² Test strain according to European Pharmacopoeia 2.6.7. Mycoplasmas

Limited Product Warranty

This warranty limits our liability for replacement of this product. No warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided. Minerva Biolabs shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

Trademarks

Venor[®] is a registered trademark of Minerva Biolabs GmbH, Germany.

Last technical revision: September 2009

Related Products

MB Taq DNA Polymerase

53-0050/-0100/-0200/-0250	MB Taq DNA Polymerase (5 U/ μ l)	50/100/200/250 units
53-1050/-1100/-1200/-1250	MB Taq DNA Polymerase (1 U/ μ l)	50/100/200/250 units

Detection Kits for endpoint PCR

11-1025/-1050/-1100/-1250	Venor [®] GeM Mycoplasma Detection Kit	25/50/100/250 tests
11-7024/-7048/-7096/-7240	Venor [®] GeM Advance Mycoplasma Detection Kit	24/48/96/240 tests
12-1025/-1050/-1100/-1250	Onar [®] EUB Bacteria Detection Kit	25/50/100/250 tests

Diagnostic Kits for qPCR

11-4025/100/250	Venor [®] GeM-qEP Mycoplasma Detection Kit	25/100/250 tests
12-2025/100/250	Onar [®] -qEUB Bacteria Detection Kit	25/100/250 tests

Mycoplasma Elimination

10-0200/0500/1000	Mynox [®] Mycoplasma Elimination Reagent	2/5/10 treatments
10-0201/0501/1001	Mynox [®] Gold Mycoplasma Elimination Reagent	2/5/10 treatments

Quantification Standards, 100 μ l each, 1x10⁶ genomes/ μ l

52-0116	<i>Acholeplasma laidlawii</i>
52-0129	<i>Mycoplasma arginini</i>
52-0117	<i>Mycoplasma fermentans</i>
52-0130	<i>Mycoplasma hyorhinis</i>
52-0112	<i>Mycoplasma orale</i>
52-0119	<i>Mycoplasma pneumoniae</i>
52-0124	<i>Mycoplasma synoviae</i>
52-0164	<i>Spiroplasma citri</i>

Genomic DNA Extracts, 100 μ l each, +/- 10 ng / 100 μ l (selection)

51-0177	<i>Ureaplasma urealyticum</i> , NC 10177
51-0111	<i>Mycoplasma hominis</i> , NC 010111
51-0112	<i>Mycoplasma orale</i> , NC 010112
51-0117	<i>Mycoplasma fermentans</i> PG19, NC 010117
51-0119	<i>Mycoplasma pneumoniae</i> , NC 010119
51-0129	<i>Mycoplasma arginini</i> , NC 010129
51-0162	<i>Mycoplasma arthritidis</i> , NC 010162
51-0195	<i>Mycoplasma genitalium</i> , NC 010195
51-1746	<i>Mycoplasma penetrans</i> , NC 11746

DNA Remover™

15-2025	DNA Decontamination Reagent, spray bottle	250 ml
15-2200	DNA Decontamination Reagent, refill bottles	4 x 500 ml

Extraction Kit

56-1100	MB DNA Extraction Kit	100 extractions
---------	-----------------------	-----------------

Mycoplasma Off[®]

15-1000	Surface Disinfectant Spray, spray bottle	1000 ml
15-5000	Surface Disinfectant Spray, refill bottles	5 x 1000 ml

Mycoplasma Off[®] Wipes

15-1001	Surface Disinfectant Wipes	120 wipes
15-5001	Surface Disinfectant Wipes, refill sachets	5 x 120 wipes

ZellShield[®]

13-0050/-0150	Contamination Prevention Reagent 100x concentrate	50 ml/ 3 x 50 ml
---------------	---	------------------

Manufacturer

Minerva Biolabs GmbH
Koeppenicker Str. 325
D-12555 Berlin
Germany

Ordering

Tel. +49 (0)30 2000 437-0
Fax +49 (0)30 2000 437-9
order@minerva-biolabs.com

Product Information

www.minerva-biolabs.com
info@minerva-biolabs.com

Technical Service

Tel. +49 (0)30 2000 437-0
support@minerva-biolabs.com

Made in Germany

Minerva Biolabs GmbH develops and manufactures products in accordance with EN ISO 9001:2008 and EN ISO 13485:2003 quality system requirement. Reg.No. SY 60026567 0001 & SX 60025009 0001

...as precise as diagnostics should be™