

# Genomic & Microbial Reference Materials



- Genomic DNA extracts of defined microorganisms for specificity testing in conventional and qPCR
- PCR quantification standards of mycoplasma and bacterial genomic DNA for applications in conventional and qPCR
- 10CFU™ Sensitivity Standards for validating robustness and detection limit of molecular mycoplasma test methods as required by the *European Pharmacopoeia 2.6.7/Japanese Pharmacopoeia, 17. edition, chapter G3*

# Genomic DNA Extracts



## Application

- Conventional and qPCR
- Specificity testing

These preparations contain genomic DNA extracted from low passage and defined microorganisms, extracted by subsequent column absorption methods. The DNA extract was partially sequenced to confirm identity. Titration was done by optical density measurement against a weight calf thymus DNA standard.

## Content

1 vial with genomic DNA extract with 10 ng +/- 2 ng, freeze-dried  
1 vial with 200 µl 10 mM Tris-HCl buffer, pH 8.5, for dissolving the DNA

## Ordering Information

51-0116	<i>Acholeplasma laidlawii</i>	51-3361	<i>Methicilin-resistant Staphylococcus aureus</i> (MRSA)
2127-30007	<i>Acinetobacter bamanii</i>	2134-30164	<i>Morganella morganii</i>
2128-30187	<i>Aeromonas hydrophila</i>	51-0030	<i>Micrococcus luteus</i>
2101-00819	<i>Aspergillus fumigatus</i>	51-0129	<i>Mycoplasma arginini</i>
51-0031	<i>Bacillus cereus</i>	51-0162	<i>Mycoplasma arthritis</i>
51-0010	<i>Bacillus subtilis</i>	51-0117	<i>Mycoplasma fermentans</i>
2129-02046	<i>Bacillus thuringiensis</i>	51-0115	<i>Mycoplasma gallisepticum</i>
51-5571	<i>Bordetella pertussis</i>	51-0195	<i>Mycoplasma genitalium</i>
2130-07288	<i>Burkholderia capacia</i>	51-0111	<i>Mycoplasma hominis</i>
2102-04688	<i>Campylobacter jejuni</i>	51-0130	<i>Mycoplasma hyorhinis</i>
51-1386	<i>Candida albicans</i>	51-0112	<i>Mycoplasma orale</i>
2103-11226	<i>Candida glabrata</i>	51-1746	<i>Mycoplasma penetrans</i>
2104-11947	<i>Candida guilliermondii</i>	51-0119	<i>Mycoplasma pneumoniae</i>
2105-70624	<i>Candida haemulonii</i>	51-0124	<i>Mycoplasma synoviae</i>
2106-30039	<i>Citrobacter freundii</i>	2135-10036	<i>Neisseria meningitidis</i>
2125-90874	<i>Candida tropicalis</i>	2136-13387	<i>Proteus vulgaris</i>
2107-04595	<i>Citrobacter koseri</i>	51-0071	<i>Pseudomonas aeruginosa</i>
51-0792	<i>Clostridium acetobutylicum</i>	2116-04479	<i>Proteus mirabilis</i>
2108-00756	<i>Clostridium perfringens</i>	51-7058	<i>Salmonella enterica</i>
51-0053	<i>Enterobacter aerogenes</i>	2117-30121	<i>Serratia marcescens</i>
2110-30054	<i>Enterobacter cloacae</i>	2137-04782	<i>Shigella flexneri</i>
2111-20680	<i>Enterococcus casseliflavus</i>	2138-05570	<i>Shigella sonnei</i>
2112-30633	<i>Enterococcus durans</i>	51-0164	<i>Spiroplasma citri</i>
51-0478	<i>Enterococcus faecalis</i>	51-0231	<i>Staphylococcus aureus</i>
2113-20477	<i>Enterococcus faecium</i>	51-0044	<i>Staphylococcus epidermidis</i>
2114-20160	<i>Enterococcus hirse</i>	2118-20328	<i>Staphylococcus hominis</i>
51-0083	<i>Escherichia coli</i>	2119-20263	<i>Staphylococcus haemolyticus</i>
2115-08579	<i>Escherichia coli</i> O157:H7	2122-20480	<i>Streptococcus bovis</i>
51-1368	<i>Fluoribacter bozemaniae</i> (syn. <i>Legionella bozemaniae</i> )	2123-06176	<i>Streptococcus dysgalactiae</i>
2131-05934	<i>Geobacillus stearothermophilus</i>	2126-20523	<i>Streptococcus mutans</i>
2132-30104	<i>Klebsiella pneumoniae</i>	51-0566	<i>Streptococcus pneumoniae</i>
2133-70603	<i>Klebsiella pneumoniae</i> , ESBL+	2139-20068	<i>Streptococcus sanguinis</i>
51-1723	<i>Lactobacillus acidophilus</i>	2120-20229	<i>Staphylococcus saprophyticus</i>
51-1370	<i>Legionella dumofii</i>	2121-50170	<i>Stenotrophomonas maltophilia</i>
51-1533	<i>Legionella jordanis</i>	51-0177	<i>Ureaplasma urealyticum</i>
51-0101	<i>Legionella pneumophila</i>	2140-04780	<i>Yersinia enterocolitica</i>
51-1514	<i>Legionella pneumophila subsp. fraseri</i>	2141-08992	<i>Yersinia pseudotuberculosis</i>
51-1515	<i>Legionella pneumophila subsp. pascuelleri</i>		

# PCR Quantification Standards



## Application

- Conventional and qPCR
- Standard curves, evaluation of assay performance
- Preparation of dilution series for quantification
- Low titer (e.g.  $3 \times \text{LOD}_{95}$ ) controls

Standardization and quantification of nucleic acid detection is a difficult task as there are no reliable standards for low nucleic acid copy numbers. Detection methods such as PCR and NASBA are prone to inhibition caused by many substances commonly found in plant and animal tissues, food matrices and extraction solutions. When co-purified with DNA or RNA, inhibitors reduce amplification efficiency, causing an underestimation of the quantity of target nucleic acid or even false-negative results.

Minerva Biolabs calibration reagents contain genomic DNA extract from low passage and defined microorganisms. The DNA is extracted by subsequent column absorption methods. The DNA extract was partially sequenced and the sequence aligned to confirm identity. Titration was done after photometric quantification of the preparation standard and dsDNA fluorometric quantification against a synthetic standard.

qPCR user are able to include a precise low count of DNA copies in their assays, providing them with a true value for the estimation of detection limits and the comparison of different detection methods.

## Content


1 vial with DNA,  $1 \times 10^8$  genomes, freeze-dried  
3 vials with 2 ml of Tris-HCl buffer, 10 mM, pH 8.5., for dissolving the DNA and preparing dilutions

## Ordering Information

52-0116	<i>Acholeplasma laidlawii</i>
52-5571	<i>Bordetella pertussis</i>
52-0083	<i>Escherichia coli</i>
52-0101	<i>Legionella pneumophila</i>
52-0129	<i>Mycoplasma arginini</i>
52-0117	<i>Mycoplasma fermentans</i>
52-0115	<i>Mycoplasma gallisepticum</i>

52-0130	<i>Mycoplasma hyorhinis</i>
52-0112	<i>Mycoplasma orale</i>
52-0119	<i>Mycoplasma pneumoniae</i>
52-0124	<i>Mycoplasma synoviae</i>
52-0164	<i>Spiroplasma citri</i>
52-0071	<i>Pseudomonas aeruginosa</i>
52-0103	<i>Mycoplasma salivarium</i>

Standard Curve



Real-time Amplification Plot




Fig. Quantification of *Mycoplasma pneumoniae* DNA. Logarithmic plot of fluorescence vs cycle number (Venor®GeM qEP, platform: ABI Prism® 7500). Template DNA ranging from  $2 \times 10^5$  - 2 genome equivalents.

# 10CFU™ Sensitivity Standards

For validating robustness and detection limit of molecular mycoplasma test methods in presence of the sample matrix.



## Application

*European Pharmacopoeia 2.6.7/Japanese Pharmacopoeia, 17. edition, chapter G3 "Mycoplasma"* requires a sensitivity of 10 CFU/ml sample volume for NAT-based methods like PCR to replace the traditional culture method. This feature of the test method must be shown by the performing lab as part of the robustness testing in presence of the sample matrix. As most cell culture labs and production facilities cannot accept vital mycoplasma in their facility or do not have access to a microbiology lab able to cultivate mycoplasma, these preparations allow safe and reliable validation of the procedure.

The mycoplasma have been cultivated in culture broth described in EP 2.6.7/JP, 17. edition, chapter G3, titrated immediately in culture broth and plated for quantification in colony forming units (CFU/ml). Each dilution series has been performed in multiple by different operators for highest precision. The mycoplasma broth was harvested in the early logarithmic phase of the growth to avoid a high ratio of dead mycoplasma particles and correspondingly a high GU\*:CFU ratio. All strains have been obtained from official culture collections and cultivated in low passages.

Each vial contains 10 CFU of inactivated mycoplasma. By adding the sample matrix of interest a sample according to EP 2.6.7/JP, 17. edition, chapter G3, is prepared which has to be tested positive by the method applied. Obviously, the inactivated sample material is not suitable for the culture method anymore. As a result of proficiency tests on DNA amplification methods for mycoplasma detection it became obvious that in means of highest sensitivity DNA extraction is indispensable. The extract can directly be used for PCR.

\* Please note: This standard material was not titrated for genome copies (GU) as EP 2.6.7/JP, 17. edition, chapter G3 does not provide sensitivity limits on DNA level. No guarantee for a particular GU:CFU ratio is provided with this product and the ratio may vary from lot to lot.

## Package Content

3 vials with 10 CFU of the corresponding mycoplasma species  
2 negative control vials

For the mycoplasma set:  
2 vials with 9 CFU of each mycoplasma species listed in the EP 2.6.7 and 2 negative controls

## Ordering Information

102-1003 <i>Mycoplasma arginini</i>	102-7003 <i>Mycoplasma hyorhinis</i>
102-2003 <i>Mycoplasma orale</i>	102-8003 <i>Acholeplasma laidlawii</i>
102-3003 <i>Mycoplasma gallisepticum</i>	102-9003 <i>Spiroplasma citri</i>
102-4003 <i>Mycoplasma pneumoniae</i>	102-1103 <i>Mycoplasma salivarium</i>
102-5003 <i>Mycoplasma synoviae</i>	102-0002 Mycoplasma Set
102-6003 <i>Mycoplasma fermentans</i>	

## How to order

Tel.: +49-30-2000437-0  
E-mail: [order@minerva-biolabs.com](mailto:order@minerva-biolabs.com)  
Internet: [www.minerva-biolabs.com](http://www.minerva-biolabs.com)